

activeMEDICAL

SUPPLIES > SOLUTIONS > SUPPORT

School Infrastructure Portfolio

Issue 1



Contents

Capability Statement	4
Ceiling Hoist System	6
Equipment Offering	7
Fixed Manual Traverse	9
Hoyer Pool Hoist	10
Alto Patient Hoist + Pivot Frame	11
Pivot Frame Sling - Allegro	12
General Purpose Sling with Head Support - Allegro	12
Change Table - Electric Hi-Lo	14
Wall Mounted Change Table	14
Examination Couch - Electric	15
Staebel Patient Transfer	15
Staebel Shower Trolley	16
Juvo Attendant Propelled Shower Commode with Sliding Footrest	16
Juvo Self Propelled Shower Commode with Swing Away Footrest	17
Tilting In Space Shower Commode	17
Regency Care Chair	18
Tom 2 Cot	19
Nino Adolescent Bed	20

CAPABILITY STATEMENT

Active Medical Supplies is a pioneer in medical device and equipment manufacturing.

The experienced team at Active Medical understand both professionally and personally, the importance of maintaining independence within the comfort of your own home or the community including educational settings. A lot of the people we help are still independent, however, just need assistive solutions as result of an injury, Hospital visit, Healthcare professional's advice or to help loved ones with their daily life.

With over 30 years of experience helping with solutions in community care, we understand sometimes that daily tasks may become difficult. Often things that we have taken for granted like getting out of bed, lifting yourself from a chair or just even moving around the house, inherit more effort than they have before. For those who come into our showroom, we can replicate the scenarios you would face at home with our daily living set ups.

Our experienced and courteous staff gladly provide all the information and assistance you require.

We work extensively with Physiotherapists, Occupational Therapists and Aged Care Clinicians so that we can find the right form of solution for you. Over the course of 30 years, we have provided solutions to community care with diligence and professionalism. Our customers' assurance comes from the fact that we help you find the most suitable solution.

We take great pride in the assembling and manufacturing of all of our products. We culminate to the highest standards in care, comfort for the patient and project cost efficiency.

Below are pictures of jobs we have completed throughout Queensland schools. These include, Merrimac State High School, Darling Point Special School, Loganlea State High School, Pine Rivers Special School, Currumbin Special school and many more.

Procurement & Distribution

Active Medical has developed a network of manufacturers in Australia and around the world – from Germany and Sweden to the United States of America and consults directly with industry professionals to ensure the products we provide responds to unique patient specifications.

The success of our business is supported by our ability to collaborate with manufacturers that not only provide the highest quality of products and services, but are aligned with our mission and commitment to our social, corporate and

environmental responsibilities.

Our efficient distribution systems are underpinned by our strong relationships with manufacturers and two centralised warehouses with a capacity to hold up to 50,000 items. This ensures the fast turn around and on time delivery demanded by the health care industry and the needs of people in care. We ensure that all our products are ethically sourced and supplied and meet environmental standards.



DID YOU KNOW...



We sanitize **all products** that leave our warehouses



Our staff are **double vaxxed** with boosters



We offer full product **demonstrations and video trials**



Our showroom is a **controlled environment** that we can also offer trials



We offer **contact free delivery**

CAPABILITY STATEMENT

OUR MISSION

To care for those who need care and those who provide care.

OUR VISION

To supply and serve the hospital, aged and home care community with world class solutions and unequalled professionalism. We will be defined by the quality of our products, the depth of our knowledge, the sincerity of our service, the practice of good ethics and our respect for patients, customers and suppliers.

OUR VALUES


Care and Kindness – for each other
Hard Work – to achieve each other's goals
Respect – for each other's needs and contributions
Doing the Right Thing – by each other
Innovation and Inspiration – to lift each other's horizon

COME IN & HAVE A CHAT!

Book a consultancy today!

Active Medical facilitates a showroom on site to replicate the home care environment. Our showrooms are used to stimulate a residential care bedroom setting and allows us to offer unique demonstrations and training tailored to your needs.

We are located in Cannon Hill and Toowoomba with wheelchair accessibility. We showcase our extensive range of electric beds, wheelchairs, mobility, patient handling, safety equipment, bathroom solutions and accessories.



Cannon Hill

35 Southgate Ave, Cannon Hill QLD 4170



Toowoomba

11 Grange St Toowoomba City QLD 4350



Showroom Opening Hours

Monday to Friday | 8:00am to 4:00pm



Showroom Opening Hours

Monday to Friday | 8:00am to 4:00pm

*Products and information listed in this catalogue are subject to change. Website and digital information supercedes all print. For the most up to date information, please contact us directly via the details listed above.


*Disclaimer

Active Medical Supplies does not take responsibility for incorrect application of products. This publication is a guide only to assist in selecting products suitable to your clients. Please use in conjunction with professional advice.

CEILING HOIST SYSTEM


Sugarglider® Basic Rail Systems

A basic rail system consists of one or more straight rails, with or without curves for increased room coverage, flexibility, and utilisation.





Fully Integrated, Sugarglider® Multi-Room System

Fully integrated systems are the ultimate in ceiling hoist solutions. By using a comprehensive range of components, we can provide opportunity to cover multiple rooms with a single system including traversing doorways, hallways, and corridors.



Sugarglider® Full Room Coverage X-Y Rail Systems

A rail system with full room coverage consists of a pair of fixed (primary) rails and one or two (secondary) cross rails/traverse rails. Each cross-rail is mounted on a pair of traverse trolleys and can travel sideways in the fixed rails along the total width of the lift system. This rail system allows the ceiling lift unit, which is mounted on a cross-rail, to travel in all directions and be used anywhere in the room. Two adjacent full room coverage rail systems can be connected by means of transition gate.



*Disclaimer

Active Medical Supplies does not take responsibility for incorrect application of products. This publication is a guide only to assist in selecting products suitable to your clients. Please use in conjunction with professional advice.

EQUIPMENT OFFERING



Bath, Shower & Toilet

A full range of bathroom necessities, including bathing aids, transfer benches, shower trolleys, commodes and accessories, toileting aids and change tables.



Standard Yoke



4 Point Yoke



Patient Handling

Patient handling, weighing and transfer equipment, such as floor, ceiling and stand-up lifters, slings and belts, slide, top and turn sheets and positioning and lifting aids.



Pressure Care

For treatment and prevention of pressure injuries, we provide low, medium and high risk cushions and mattresses, heel and elbow protectors and Action gel products.



Wheelchairs

Included in our range are general purpose and lightweight adult and paediatric wheelchairs, Bariatric wheelchairs, ward chairs, accessories and positioning aids.

EQUIPMENT OFFERING



Seating

An extensive range of chairs including flotation, care, Bariatric, mobile and electric lifter chairs and accessories, as well as orthopaedic and specialist chairs.



Mobility Solutions

Axillary, elbow and forearm crutches, walking sticks and accessories, walking frames, mobility accessories, seat walkers and forearm walkers are all on offer.



Bedroom Solutions

High quality beds and bedroom solutions for Hospital, Paediatric and Bariatric needs, plus safety and comfort accessories and commodes.



Patient Safety

Ensure patient monitoring with Invisabveam and Cura1 alarms, systems and accessories and increase safety with hip protectors, footwear and impact mats.

PATIENT HANDLING



SugarGlider Ceiling Hoist

The SugarGlider Ceiling Hoist delivers innovation, affordability and quality. Ceiling hoists dramatically reduce injuries resulting from the handling of patients and residents in both institutional and home care settings. The SugarGlider Ceiling Hoist has been designed with a full set of standard features and can be safely operated by a single caregiver.

Built for quiet operation, the hoist provides a smooth & effortless transfer. The SugarGlider Ceiling Hoist will lift up to a maximum of 205kg, allowing the greatest flexibility in the type of patients that can be transferred. For additional flexibility, powered traversing options can be included.

- Smallest Hoist of its class which allows for 180mm of additional lifting height
- Quick lifting speed
- Emergency pull cord
- Modular design for simplified maintenance
- Quick charge batteries for reduced service time
- Audible low battery indicator & visual battery/charge level display Digital readout
- Auto shut off when not in use, to save battery life
- 2.2 strap vertical movement allows you to easily lift from the floor

- Pneumatic Hand Controls – 100% waterproof
- Manual option – for emergency raising or lowering
- Available with return to charge option
- Lift Motor 24V DC
- Traverse Motor 24V DC
- Charger Input 240V AC, 1.0 Amps
- Charger Output 24V DC, 1.5 Amps
- Battery 24V (2 x 12V DC), 5.0AH
- Hoist Case Flame Retardant ABS
- Hand Control Pneumatic 100% waterproof

SPECIFICATIONS:

Lift Speed	46 mm
Lifting Height	2130 mm
Unit Weight Under	10 kg
SWL	205kg

Patient Handling



Hoyer Pool Hoist

- This hydraulic pool lifter makes recreational and therapeutic swimming safer and more enjoyable for both the person being lifted and the attendant. Designed for the person that requires assistance in transfer.
- Easily operated hydraulic power operation.
- Full 158mm lifting range and 360 degree swivel to bring user gently from chair to below water.
- Heavy duty construction.
- Stainless steel painted finish to resist the humid and corrosive swimming pool environment.
- Flush mounted recessed socket (included)
- No under water installation.
- Maximum Weight Capacity 182kg
- Lift Range 1570mm


SPECIFICATIONS:

Product Code	2835
Lift Range (Length of Lift Strap)	1570mm
SWL	182kg

PATIENT HANDLING

**SPACE LINK™ CONCEPT**

This breakthrough technology is for optimized patient space and enhances floor lift capacity.



Standard Yoke



4 Point Yoke



Tri Hook Yoke



Pivot Frame



Powered Pivot Frame



Weigh Scale

Alto Patient Hoist + Powered Pivot Frame


Alto Lift 200: A medium sized heavy duty mobile patient lifting hoist with a safe working load of 200kg and a compact ergonomic design. The Alto Lift 200, whilst being compact, provides superior patient lift space through the use of SPACE LINK™ technology.

The Alto Lift 200 is the lifter of choice for institutional use or home care where a high capacity lift is required. Designed and tested to comply with: AS ISO 10535-2002 Hoists for transfer of disabled persons - Requirements and test methods.

SPECIFICATIONS:

Product Code	2892	Minimum Storage Height	1350 mm
Base Length	1150 mm	Under Base Height	50 mm
External Base Width	645 mm	Maximum Hoist Reach	1060 mm
Maximum Internal Length	860 mm	Minimum Hoist Reach	610 mm
Internal Base Width Closed	578 mm	Max Hoist Reach at	41 kg
Internal Base Width Open	1130 mm	Max Height	4400 mm
Base Width Adjustment	Electric	Min Hoist Reach at Min Height	430 mm
Under Bed Clearance	125 mm	Handle Height	1000 mm
Minimum Height of	420 mm	Turning Radius	1210 mm
Spreader Bar	300 mm	Mass of Heaviest Part (Base Assembly)	17.5 kg
Maximum Height of	870 mm	Lifter Weight	42 kg
Spreader Bar	1920 mm	SWL	200 kg
Hoisting Range*	1620 mm		

PATIENT HANDLING



Pivot Frame Sling - Allegro

The Pivot sling is designed for use with the pivot/tilt frame. It attaches with the keyhole tags engineered to the highest standards in strength and size tolerances and feature auto clip on action (no safety straps required). It also features sling centre marking and positioning handles to assist correct positioning and woven label that will not wash out or fade.

Available in sizes from x-small through to XX-large and in fabric or mesh styles.

General Purpose Sling with Head Support - Allegro


The General Purpose sling is a standard sling providing good support for those with some upper body control. Suitable for general lifting for patient lifters with a yoke style attachment, and available in sizes from X-small to XX-large. Also with the option of fabric or mesh construction, this is a versatile sling that is simple and safe to use.

Also available:
Allegro General Purpose Hygiene Slings With Head Support

Allegro Mesh		Allegro Fabric				Allegro Mesh		Allegro Fabric					
XS	27P02XS	XS	27P01XS	XS	27GPHS02XS	XS	27GPHS01XS	XS	27GPHS01XS				
S	27P02S	S	27P01S	S	27GPHS02S	S	27GPHS01S	S	27GPHS01S				
M	27P02M	M	27P01M	M	27GPHS02M	M	27GPHS01M	M	27GPHS01M				
L	27P02L	L	27P01L	L	27GPHS02L	L	27GPHS01L	L	27GPHS01L				
XL	27P02XL	XL	27P01XL	XL	27GPHS02XL	XL	27GPHS01XL	XL	27GPHS01XL				
XXL	27P02XXL	XXL	27P01XXL	XXL	27GPHS02XXL	XXL	27GPHS01XXL	XXL	27GPHS01XXL				
	XS	S	M	L	XL	XXL							
A	720	870	1000	1000	1170	1320	A	820	790	920	1110	1270	1360
B	190	230	270	260	270	290	B	210	230	240	260	270	290
C	580	830	920	1010	1090	1180	C	750	780	920	1000	1020	1050
D	870	1220	1310	1410	1540	1710	D	1310	1280	1410	1630	1680	1780

Allegro Mesh		Allegro Fabric				Allegro Mesh		Allegro Fabric					
XS	27P02XS	XS	27P01XS	XS	27GPHS02XS	XS	27GPHS01XS	XS	27GPHS01XS				
S	27P02S	S	27P01S	S	27GPHS02S	S	27GPHS01S	S	27GPHS01S				
M	27P02M	M	27P01M	M	27GPHS02M	M	27GPHS01M	M	27GPHS01M				
L	27P02L	L	27P01L	L	27GPHS02L	L	27GPHS01L	L	27GPHS01L				
XL	27P02XL	XL	27P01XL	XL	27GPHS02XL	XL	27GPHS01XL	XL	27GPHS01XL				
XXL	27P02XXL	XXL	27P01XXL	XXL	27GPHS02XXL	XXL	27GPHS01XXL	XXL	27GPHS01XXL				
	XS	S	M	L	XL	XXL							
A	720	870	1000	1000	1170	1320	A	820	790	920	1110	1270	1360
B	190	230	270	260	270	290	B	210	230	240	260	270	290
C	580	830	920	1010	1090	1180	C	750	780	920	1000	1020	1050
D	870	1220	1310	1410	1540	1710	D	1310	1280	1410	1630	1680	1780

NEED A CUSTOM SLING? LET'S HAVE A CHAT!



At Active Medical we specialise in finding the perfect fit for any individual!

We understand that each person is unique and have 30 years' experience specialising in custom slings to suit their specific need. A custom sling can benefit both carer & individual reach their full potential in everyday life.

REACH OUT TODAY

TO ENQUIRE ABOUT CUSTOMISABLE SLINGS.



1800 267 267



homecare@activemedicalsupplies.com.au



35 Southgate Ave
Cannon Hill QLD 4170



www.activemedicalsupplies.com.au

PATIENT HANDLING



Change Table - Electric Hi-Lo

This new and improved Change Table includes the following features making the ideal option for multiple environments. With a 450-900mm height adjustment range and a 250kg weight limit, it suits a larger range of clients.

The vinyl padded top is easy to clean and the base has been designed for easy hoist access. Fold downside rails are an optional extra.

- Vinyl upholstered top 700mm wide x 1880mm long
- Electric hi-lo (450mm – 900mm)
- 125mm braking castors
- Designed for hoist access
- Handset operated
- 250kg capacity
- Quality Australian made product

SPECIFICATIONS:

Product Code	7805
Low Height	450 mm
High Height	900 mm
Lying surface Width	700 mm
Lying Surface Length	1880 mm
SWL	250 kg



Wall Mounted Change Table

- Electric Hi Low
- Padded Top
- Battery Backup
- Anti-Bacterial Powder Coating
- Infa Red Handset (Cordless)
- Safety Side Rails
- Fasteners for Fixing to Wall



SPECIFICATIONS:

Product Code	7808
Low Height	350 mm
High Height	950 mm
Lying surface Width	720 mm
Lying Surface Length	1420 mm
SWL	170 kg

PATIENT HANDLING



Examination Couch - Electric

The heavy duty examination table has been designed for stability for heavier patients, and the look is centered on well- designed compact simplicity. It is manufactured with an electric height adjustment with twin towers, which are straight forward to operate, enhancing ease for care activities with its range of heights.

- Electric hi-lo 960-560mm
- Electric backrest raise/lower
- 250kg capacity
- Vinyl upholstered top 700mm x 1880mm
- Handset operated
- 125mm braking castors
- Quality Australian made product

SPECIFICATIONS:

Product Code	7840
Low Height	560 mm
High Height	960 mm
Lying surface Width	700 mm
Lying Surface Length	1880 mm
SWL	250 kg



Staebel Patient Transfer

The Staebel Patient Transfer Device is the perfect solution for safe transfer of a patient or resident within the home or facility. This is your all-in-one solution without the need for additional lifting straps and belts, and promoting independence. The width adjustable base allows good access to most seating, whilst also designed for good access beneath beds etc.



SPECIFICATIONS:

Product Code	212145
Seat Width	450 mm
Overall Width	640 mm
Overall Length	930 mm
SWL	150 kg

BATH, SHOWER & TOILET



Stabel Shower Trolley

- Shower Trolley With Two-Section Top
- Trolley Top Length: 190 Cm, Trolley Top Width: 65 Cm
- Adjustable 5° Trolley Top Inclination
- Electrically Lowering/Raising Of Trolley Top Via Console From 54 Cm To 97 Cm
- Electrically Adjustable Head Section Inclination 45° (Head Section Length: 70 Cm)
- Castors With Diameter Of 150 Mm
- Central Braking System
- Easy Steering System (For Ease Of Manoeuvrings Trolley's Front Left Castor Can Be Locked In Direction Of Steering)
- Side Rails – Foldable, Height: 18,5 Cm With Safety Release Mechanism (Material – Powder Coated Stainless Steel)
- Metal Base Cover
- Water-Resistant Shower Mattress (Height: 27 Cm)
- Mattress Integrated Drain With 645 – 2000 Mm Long Flexible Hose
- Water-Resistant Head Pillow – Semi-Cylinder
- Trolley Weight: 110 Kg
- Safe Working Load: 170 Kg
- Comes In Multiple Colour

SPECIFICATIONS:

Product Code	1209ST
Shower trolley top Width	650 mm
Shower trolley top Length	1900 mm
Shower trolley top Height	540 - 970 mm
Height with Casters	600 - 1030 mm
Trolley top inclination	0° - 5°
Length	2070 mm
Width	810 mm
SWL	170 kg



Juvo Attendant Propelled Shower Commode with Sliding Footrest

The ultimate in strength, quality, and stability the Juvo mobile shower chair is ideal for the Petite client! The wheelbase design prevents forward tipping and has proven to be an effective solution in all areas from home to aged care to hospital.

- Attendant propelled
- Quality stainless-steel frame for longevity in a demanding environment
- Flip up armrests for ease of side transfers
- Armrests can be flipped so that they are positioned in front of occupant to reduce risk of slipping forward when being wheeled
- Smooth, slatted backrest design to facilitate easy cleaning
- Generous size sliding stainless steel footplate, weight bearing to utilise as a commode you must source the slide-in pan rack

* Please note that the seat for this product is sold separately

Seat colours can be varied.



SPECIFICATIONS:

Product Code	1904J
Seat Width	410 mm
Width Between Arms	440 mm
Overall Height	1020 mm
Backrest Height	340 mm
Seat to Heel	400 mm
Anterior Seat Height	595 mm
Toilet Clearance – To Pan	315 mm
Toilet Clearance – Under Pan	480 mm
Toilet Clearance – Height	500 mm
Toilet Clearance – Width	360 mm
SWL	200 kg

MORE SIZES AVAILABLE

BATH, SHOWER & TOILET



Juvo Self Propelled Shower Commode with Swing Away Footrest

The ultimate in strength, quality, and stability the Juvo mobile shower commode is built for heavy use. The rear self-propel wheels are ideal for the independent user.

- Assists transportation of the patient for showering and toileting
- Can be positioned over the toilet
- Durable total stainless-steel frame
- Selection of seamless PU padded seats & pressure care seats (open or closed fronts)
- Removable seat
- Reinforced seat/back frame
- Swing away, removable & height adjustable leg rests with flip-up footplates.
- Patented TEARSAFE design, no protrusions
- Now with braking stainless steel 125mm castors
- Swing back independent armrests to assist inside transfer and patient access
- Two position height adjustable armrests
- Safety arms provide user security & safety
- Stainless steel pan rack guides
- This product is tested and certified to Australian Standards
- Please note that the seat, pan rack and pan for this product are sold separately

Seat colours can be varied.



SPECIFICATIONS:

Product Code	1925J
Seat Width	460 mm
Anterior Seat Height	590 mm
Posterior Seat Height	565 mm
Width Between Arms	465 mm
SWL	200 kg

MORE SIZES AVAILABLE



Tilting In Space Shower Commode

- Assists transportation of the patient for showering and toileting
- Gas strut tilt mechanism for smooth and easy tilting
- Multi-position 35 deg. tilt
- Can be positioned over the toilet
- Durable total stainless-steel frame
- Removable seat for ease of cleaning, no tools required
- Reinforced seat/back frame
- Swing away, removable and height adjustable leg rests with flip up footplates
- Patented TEARSAFE design, no protrusions
- Four Stainless steel 125mm braking castors
- *Swing back independent armrests to assist inside transfer and patient access*
- Adjustable seat height
- Two position height adjustable armrests
- Safety arms provide user security and safety for transit
- Stainless steel pan rack guides
- * Please note that the seat for this product is sold separately Selection of seamless PU padded seats & pressure care seats (open or closed front)

Seat colours can be varied.



SPECIFICATIONS:

Product Code	1951J
Seat Width	460 mm
Seat to Heel	410 mm
Overall Width	570 mm
Overall Depth Upright	900 mm
Overall Height	1490 mm
Backrest Height from Seat	470 mm
Width Between Arms	470 mm
Armrest Height	200 - 240 mm
SWL	200 kg

MORE SIZES AVAILABLE

SEATING



Regency Care Chair

The Regency Care Chair is designed to give the patient comfort by combining the original tilt system with an independent multi-positional back rest. Its many features take into account the needs of both carer and user of the chair. The Regency R2900 derives its status from a proven and successful track record across a broad spectrum of the industry.

- Reclining backrest – full recline to supine position
- Tilting seat
- Pressure management seating – gel cushion in seat, self-air adjusting back, profiled head cushion
- 2-position armrests

- Forward seat tilt for ease to stand, legrest angles under seat
- Elevating contoured legrest complete with fold out, extendible footplate
- Low skin shear 2-way stretch fabric on seating surfaces
- Drop down armrests and swingaway headwings – side transfer friendly
- Covers replaceable
- Patient lifter friendly castor base
- 5 year limited warranty
- 175mm castors, braking on rear, directional lock on front
- Quality Australian made construction

Seat colours can be varied. Enquire today for more colour options!



SPECIFICATIONS:

Product Code	5050G	Legrest Length	480 - 650 mm
Seat Width	490 mm	Overall Width	660 mm
Seat Depth	520 mm	Overall Depth	1040 mm
Seat Height (Tilt/Normal)	570 - 620 mm	Overall Height	1335 mm
Backrest Height	700 mm	SWL	180 kg
Armrest Height	175 mm		

TOM 2



Tom 2 Cot

A state-of-the-art children's cot must be safe and be able to provide excellent access for nursing staff. The electric positioning of TOM 2 children's cot responds to these and further demands in its original concept of telescopic siderails, column construction mechanism and above standard safety parameters.

The Tom 2 bed fully complies with the new standard for hospital children's beds. DIN 32623:2009-11 and its construction guarantees perfect protection against children fall. Also complies with standards: EN 60601-1, EMC 60601-1-2, EMC 60601-1-6

- 3 height positions for side rails, 0cm, 40cm, 80cm
- Backrest retracts when it raises, meaning zero gap
- Integrated collapsible lockable side rails
- 4 x Tente Motion 125mm castors

- Central operated castors (one castor with steer system)
- 4 x holders for IV in corners of bed
- Undercarriage cover
- Removable head/foot ends
- Accessory holders: Pair of accessory rails (located on both sides under side rails)
- Battery Backup
- Horizontal Wall/corner bumpers
- Pastal Blue Trim with Teddy Decor
- Quality European product

SPECIFICATIONS:

Product Code	7216
Mattress Platform	1400 x 700 mm
Electric hi-lo	580 – 830 mm
Electric Tilt/Reverse Tilt	15°
Backrest Adjustment	40°
SWL	75 kg




NINO



Nino Adolescent Bed

The Nino Series allows children to perceive the world around them while at the same time providing protection and security. The transparent Perspex merges with the room allowing occupants to have an unrestricted view of their immediate surroundings yet still protects against falls and injuries. Also available in different powdercoat colours, or with higher door elements or doors on both sides.


- Battery backup
- Underbed lighting
- Backlit handset
- White laminate



Standard features also include:

- Nino Series – Single Width
- Electric Hi-low
- Electric Backrest
- Electric Kneebreak
- Standard Locking Castors
- Gates one side with Perspex infills
- Fixed Panel one side with Perspex Infills
- 750mm High Sides (from Mattress Base)
- Head & Foot Infills in Natural Oak Matt Finish
- Battery Backup
- Doors one side (750mm high perspex sides)
- 125mm braking and directional castors

Powdercoat Colours Available:



NINO



Configuration Possibilities



SPECIFICATIONS:

Overall dimensions 2200 x 1014 mm

Sleeping surface dimensions 2000 x 870 mm

NET Bed weight 140 kg


Underbed clearance 100 mm (at lowest point)

Low Position 390 mm

High Position 830 mm

SWL 200 kg





SEE SOMETHING YOU LIKE? LET'S HAVE A CHAT!

Book your consultancy today with one of our friendly Team Members. See below information for ways you can get in contact...



1800 267 267



homecare@activemedicalsupplies.com.au



35 Southgate Ave
Cannon Hill QLD 4170



www.activemedicalsupplies.com.au

The background is a solid teal color. In the top left, there is a stylized sun with a spiral center and radiating lines. In the top right, there is a white grid pattern. In the center, there is a simple line drawing of an umbrella. On the right side, there are several overlapping circles of varying shades of teal, some with smaller circles inside them, resembling a pattern of bubbles or cells. In the bottom left, there is a stylized flower with five petals.

activeMEDICAL

S U P P L I E S > S O L U T I O N S > S U P P O R T

35 Southgate Ave, Cannon Hill QLD 4170
P. 1800 267 267 | F. 07 3899 8825
E. hello@activemedicalsupplies.com.au

www.activemedicalsupplies.com.au