

**activeMEDICAL**

SUPPLIES > SOLUTIONS > SUPPORT

# Cura 1

Falls Prevention &  
Monitoring Devices

*We believe in partnering with healthcare professionals within allied health, to provide better patient outcomes through quality of our equipment solutions.*

## Contents

Capability Statement	4
Monitors	6
ActiveCare Solutions	10
iBeam	12
Cordless Pads, Mats & Transmitters	13
HomeCare Alarms	14

## CAPABILITY STATEMENT

**Active Medical Supplies is a pioneer in medical device and equipment manufacturing.** The experienced team at Active Medical understand both professionally and personally, the importance of maintaining independence within the comfort of your own home or the community including educational settings. A lot of the people we help are still independent, however, just need assistive solutions as result of an injury, Hospital visit, Healthcare professional's advice or to help loved ones with their daily life.

With over 30 years of experience helping with solutions in community care, we understand sometimes that daily tasks may become difficult. Often things that we have taken for granted like getting out of bed, lifting yourself from a chair or just even moving around the house, inherit more effort than they have before.

For those who come into our showroom, we can replicate the scenarios you would face at home with our daily living set ups. Our experienced and courteous staff gladly provide all the information and assistance you require.

We work extensively with Physiotherapists, Occupational Therapists and Aged Care Clinicians so that we can find the right form of solution for you. Over the course of 30 years, we have provided solutions to community care with diligence and professionalism. Our customers' assurance comes from the fact that we help you find the most suitable solution.

We take great pride in the assembling and manufacturing of all of our products. We culminate to the highest standards in care, comfort for the patient and project cost efficiency.

## CAPABILITY STATEMENT

### Did you know...



We sanitize **all products** that leave our warehouses



We offer full product **demonstrations and video trials**



We offer **contact free delivery**



Our staff are **double vaxxed** with boosters



Our showroom is a **controlled environment** that we can also offer trials

## Procurement & Distribution

Active Medical has developed a network of manufacturers in Australia and around the world – from Germany and Sweden to the United States of America and consults directly with industry professionals to ensure the products we provide responds to unique patient specifications.

The success of our business is supported by our ability to collaborate with manufacturers that not only provide the highest quality of products and services, but are aligned with our mission and commitment to our social, corporate and environmental responsibilities.


Our efficient distribution systems are underpinned by our strong relationships with manufacturers and two centralised warehouses with a capacity to hold up to 50,000 items. This ensures the fast turn around and on time delivery demanded by the health care industry and the needs of people in care. We ensure that all our products are ethically sourced and supplied and meet environmental standards.



## Come in & have a chat! Book a consultancy today!

Active Medical facilitates a showroom on site to replicate the home care environment. Our showrooms are used to stimulate a residential care bedroom setting and allows us to offer unique demonstrations and training tailored to your needs.

We are located in Cannon Hill and Toowoomba with wheelchair accessibility. We showcase our extensive range of electric beds, wheelchairs, mobility, patient handling, safety equipment, bathroom solutions and accessories.



**Cannon Hill**  
35 Southgate Ave, Cannon Hill QLD 4170



**Toowoomba**  
11 Grange St Toowoomba City QLD 4350



**Showroom Opening Hours**  
Monday to Friday | 8:00am to 4:00pm



**Showroom Opening Hours**  
Monday to Friday | 8:00am to 4:00pm

**\*Disclaimer**

Active Medical Supplies does not take responsibility for incorrect application of products. This publication is a guide only to assist in selecting products suitable to your clients. Please use in conjunction with professional advice.

\*Products and information listed in this catalogue are subject to change. Website and digital information supercedes all print. For the most up to date information, please contact us directly via the details listed above.

# MONITORS

# MONITORS



## Six Channel Cordless Falls Monitor (Secure Reset)

## Under Mattress Bed Pads

No more cords between monitor and pads. Can be used with up to six cordless sensor pads or other triggering devices.

The ultimate in patient comfort and system reliability. This bed pad is discretely installed underneath the mattress. Please be certain that the bed frame is compatible (has a flat area for the pad to sit on - see diagram below). Recommended weight range for reliable operation: 40kg to 200kg.

**Features:**


- Monitor up to 6 wireless devices simultaneously.
- Works with a single Premium Pad or Mat.
- Alarm memory and auto re-arm.
- Amber LED visual alarm indicator.
- Battery operated or optional external plug pack.
- Up to 6 months battery life (nominal).
- Easy to access battery compartment.
- Can be connected to a nurse call system.
- 3 selectable alarm tones.
- Easy thumbwheel adjustable alarm volume.
- "Pause" feature allows caregiver to attend to the patient without setting the alarm off.

**Safety Features:**

- Alerts for lost nurse call connection, lost AC power supply & low battery
- All safety alerts sound even if alarm volume is muted, including low battery, AC lost, pad signal lost, nurse call lost & change pad.
- Cable strain relief & breakaway cord technology prevents damage & reduces risk of entanglement when used with Premium Pads & Mats.
- Non resettable "Pad Signal Lost" warning.
- Caregiver secure reset option.

**What's included:**

- 1 x 2700 Six Channel Cordless Falls Monitor
- 1 x 2742 Soft Monitor Cover
- 1 x 2642 Caregiver Reset Key
- 1 x Non-Absorbent Mounting Strap
- 3 x Alkaline AA Batteries



**SPECIFICATIONS:**

Product Code	Name
93156	ActiveCare Console
93515	ActiveCare Pager

**SPECIFICATIONS:**

Product Code	92700
--------------	-------



## Falls Monitor with SafeTime Support

Connect with any cura" Premium or Safetime Bed / Chair Pad or Floor Mat to alert a carer when a monitored patient gets out of bed. Compatible with Premium and SafeTime Bed & Chair Pads and Floor Mats.

Comes with mounting bracket, hanger strap, soft cover, reset key and batteries.

**SPECIFICATIONS:**

Product Code	92701
--------------	-------



## Cordless Multiport Expander Bracket Mount

For nurse call panel configurations that are too small or crowded to accept the 2544 Cordless Multiport Expander, there is the 2549 Cordless Multiport Expander Bracket Mount, which has its nurse call plug on a Curly Cord, allowing the unit to be installed o" the panel.

**SPECIFICATIONS:**

Product Code	92549
--------------	-------



## Economy Incontinence Monitor

Simple operation and maintenance. Works with all cura" Economy Incontinence Sensors. Red LED visual alarm indicator. 3 x AAA battery operated or optional 6V external plug pack. 3 selectable alarm tones. Can be connected to a nurse call system.

**SPECIFICATIONS:**

Product Code	94020
--------------	-------



## Wireless Call Button with Pull Cord

An IP 65 rated wireless call button for use with other cura" wireless equipment. Fully wireless - no cabling required.

**SPECIFICATIONS:**

Product Code	92558
--------------	-------



## Cordless Multiport Expander

Plugs into a nurse call system and can connect up to five cura" wireless devices (2 x pad & 3 x other transmitters eg PET), a 2638 Nurse Call Pendant and one other wired device (eg a Premium Floor Mat).

**SPECIFICATIONS:**

Product Code	92544
--------------	-------



## Personal Emergency Transmitter (PET)

A low-cost panic button with lanyard and a large splash resistant button for ease of use. Includes a lanyard or use with a hospital wristband.

**SPECIFICATIONS:**

Product Code	92588
--------------	-------



## MONITORS



### Wireless Handset

Wireless patient handset with call button and reassurance LED. Works with nurse call or direct to a pager. Includes wall mount bracket for when handset is not in use.

#### SPECIFICATIONS:

Product Code 927001



### Nurse Call Pendant

Simple operation and maintenance and includes strain relief and a metal clip on the cable for securing to a bed.

#### SPECIFICATIONS:

Product Code 92406




### Plug Pack

6Vdc 1A (regulated) suitable for many cura" products. Being white in colour makes it aesthetically suitable for most clinical environments.

#### SPECIFICATIONS:

Product Code 92638




### Bed Incontinence Sensor

Made from soft non-allergenic PVC. Reusable & easy to clean. For use with Economy Incontinence Monitor.

#### SPECIFICATIONS:

Product Code 94021



### Chair Incontinence Sensor

Made from soft non-allergenic PVC. Reusable & easy to clean. For use with Economy Incontinence Monitor.

#### SPECIFICATIONS:

Product Code 94022

## MONITORS



### Premium & Safetime Bed / Chair Pads & Floor Mats

For use with cura" Premium Falls Monitors. SafeTime Pads require a monitor with the SafeTime feature.

#### SPECIFICATIONS:

Product Code	Name
92623	Premium Floor Mat
92626	Premium Bed Pad
92628	Premium Chair Pad



### Cordless Bed / Chair Pads & Floor Mats

These Pads and Mats can be used with a 2719 Transmitter to connect to your cura" Cordless Monitor.

#### SPECIFICATIONS:

Product Code	Name
92624	Cordless Floor Mat
92659	Cordless Bed Pad
92630	Cordless Chair Pad

## ACTIVECARE SOLUTIONS



### Basic Falls Solution

When a person gets out of bed the Under Mattress Pad sends an alarm signal to the ActiveCare Pager.

Pressing the red button on the Personal Emergency Transmitter (PET) sends an alarm signal to the ActiveCare Pager. Multiple PETs can be used with a single Pager.

- Simple to use.
- Alarm sent to Pager.

### Advanced Falls Solution

The ActiveCare Pager provides the carer full mobility around the house with up to 9 devices connected.

Pressing the red button on the Handset sends an alarm signal to the ActiveCare Pager.

Wireless Waterproof Call and Reset devices can be installed in wet areas, with an emergency pull-cord designed to be accessible for someone that has fallen.

Wireless Magnetic Door Switches can be used to monitor exits such as doors as well as windows.

## ACTIVECARE SOLUTIONS



### Expanded Falls Solution

Add an ActiveCare Console for more expansion options.

For homes where more than 9 devices are used, ActiveCare Console can monitor up to 60 devices. (Also suitable for wards and small health facilities.)

Alarms displayed on the ActiveCare Console are also displayed on the ActiveCare Pager.

# THE EVOLUTION OF MODERN HEALTHCARE

Cura1 products are designed to augment well designed care programs & practices, Cura1 does not claim that any of its products are suitable as a substitute, or as a total solution for individuals that require care and constant supervision & are at risk. All our products are supplied with full warranty & are guaranteed to be free from defects in material & workmanship.



Cura1 has been supplying industry standard Fall Prevention and Anti-Wandering Solutions for almost two decades. Each of our systems must pass rigorous testing and exceeds industry standards for quality and reliability. Our world-class customer service team is on call to ensure your 100% satisfaction.



## IBeam



### iBeam

The narrow beam creates a precise detection zone.

This is especially useful for preventing unwanted alarms if other people frequent the same room, or where only a narrow area requires monitoring. A narrow beam allows a patient to move in their bed without triggering an alarm.

Whether across a door, stairs or beside a bed, there are countless ways that our iBeam can be used to protect those in your care and to make the job of being a caregiver easier.

When used with our sturdy, adjustable-height EzyStand, precise positioning of the iBeam could not be simpler.


### SPECIFICATIONS:

Product Code 92721



### USE WITH:

Product Code	Name
92700	Six Channel Cordless Falls Monitor
92544	Cordless Multiport Expander
93501	ActiveCare Pager
93500	ActiveCare Consol
92424	12 VDC Plug Pack



### Ezy Stand


The EzyStand / iBeam combination allows a caregiver to quickly deploy targeted motion detection alarms in a variety of positions such as next to a bed, in front of a chair, near a door or stairs.

Our innovative 2723 EzyStand magnetically attaches to the 2721 iBeam and significantly adds to its ease of use.

The stand is sturdy, attractive and perfect for use in both professional and home environments.

### SPECIFICATIONS:

Product Code 92723



## CORDLESS PADS, MATS & TRANSMITTER



### Universal Cordless Pads, Mats & Transmitter

- Guaranteed 18 month warranty
- Universal Bed Pads are foldable (single fold) Simple operation and maintenance
- When fitted with a Universal Transmitter (2720), Pads & Mats work with cura1 Cordless Receivers eg. 2700, 2544, 2549, 3500, 3501 & 3508
- When fitted with a USB to RJ10 Breakaway Cable (2693), Pads & Mats work with a Premium Monitor or receiver eg. 2701, 2544
- Easy to clean non-allergenic soft PVC
- Floor Mat has low profile edges to reduce tripping hazards & allow easy wheelchair access
- Foldable Floor Mats in beige colour to help some Alzheimer patients who have trouble with dark mats being perceived as a "dark hole in the floor"

### SPECIFICATIONS:

Floor Mat	
Product Code	92692
Width	610 mm
Length	1220 mm
Width (Single Fold)	600 mm
Length (Single Fold)	900 mm
Width (Double Fold)	610 mm
Length (Double Fold)	1220 mm

### SPECIFICATIONS:

Chair Pad	
Product Code	92671
Width	250 mm
Length	380 mm

### SPECIFICATIONS:

Product Code	Name
92691	Universal Chair Pad
92690	Universal Bed Pad
92693	Mini USB to RJ10 Breakaway Cable
92720	Universal Cordless Transmitter

### SPECIFICATIONS:

Bed Pad	
Product Code	92672
Width	250 mm
Length	760 mm



## HEMOCARE ALARMS



### Wireless Call Buttons With Caregiver Pager


When mobility is precious or restricted, a push button paging system preserves independence and privacy while staying secure.

- Super-simple setup and operation.
- Portable (no wires), versatile and long battery life.
- Pager range is suitable for typical residences.\*
- Audible chime / alarm pager options.
- Ideal for those recovering from illness, the elderly or people with disabilities.
- An affordable aid to assisted living for seniors.
- Call Buttons can be worn as a pendant, carried in a pocket, or mounted on a wall, bed or wheelchair.
- Lanyard straps and mounts included.
- Pager includes belt clip.

### Wireless Motion Sensor With LED Page

Ideal for carers monitoring the movements of frail, elderly or dementia sufferers, for when they approach dangerous areas, or get out of a bed or chair.

- Super-simple setup and operation, just position the sensor, then switch on and test.
- Exceptionally narrow detection zone reduces false alarms.
- Portable (no wires), versatile and long battery life.
- Pager range is suitable for typical residences.\*
- Audible chime / vibrate pager options.
- Motion Sensor can be silent when an alarm condition is active.
- Sensor detects motion up to 8m.
- Sensor can be powered by optional AC adaptor (not included).



Affordable home care solutions to increase quality of life for the elderly, frail or dementia sufferers and their carers.



## HEMOCARE ALARMS



### Optional Expansion Capabilities of LED Pager

- Can use multiple
- Pagers in the same system
- Long battery life
- Low battery warning
- Includes belt clip

### Seizure Alert Monitor

The cura1 SeizureAlert Alarm System provides peace of mind for carers and sufferers of convulsive seizures.

It detects & raises an alarm as a seizure occurs, providing the person being monitored a sense of privacy with the comfort and security that sleep time will be a safe and enjoyable experience.

- 12Vdc (standby current consumption < 100mA)
- Power fail LED indication
- 24 hours standby battery backup
- Fully user programmable
- Plug Pack included
- Adjustable sensitivity controls
- Durable ABS plastic construction
- Normally open relay alarm output

### SPECIFICATIONS:

Product Code	93010
Width	120 mm
Length	155 mm
Height	40 mm
Weight	350 g

### USE WITH:

Product Code	Name
92721	iBeam
92659	Cordless Bed Pad
92630	Cordless Chair Pad
92624	Cordless Floor Mat
92588	Personal Emergency Transmitter
2503	Wireless Magnetic Door Switch
93504	ActiveCare Bracket
92723	EzyStand

### USE WITH:

Product Code	Name
93508	LED Pager
92588	Personal Emergency Transmitter
	Cordless Pads / Mats





**activeMEDICAL**

S U P P L I E S > S O L U T I O N S > S U P P O R T

**P.** 1800 267 267 | **F.** 07 3899 8825

**E.** [hello@activemedicalsupplies.com.au](mailto:hello@activemedicalsupplies.com.au)

**[www.activemedicalsupplies.com.au](http://www.activemedicalsupplies.com.au)**