

activeMEDICAL
SUPPLIES > SOLUTIONS > SUPPORT

staebel 
Products

Issue 1

Contents

Capability Statement	4
Bed, Bath & Toilet	6
Patient Handling	12
Pressure Care	27
Seating & Mobility	30

CAPABILITY STATEMENT

Active Medical Supplies is a pioneer in medical device and equipment manufacturing. The experienced team at Active Medical understand both professionally and personally, the importance of maintaining independence within the comfort of your own home or the community including educational settings. A lot of the people we help are still independent, however, just need assistive solutions as result of an injury, Hospital visit, Healthcare professional's advice or to help loved ones with their daily life.

With over 30 years of experience helping with solutions in community care, we understand sometimes that daily tasks may become difficult. Often things that we have taken for granted like getting out of bed, lifting yourself from a chair or just even moving around the house, inherit more effort than they have before.

For those who come into our showroom, we can replicate the scenarios you would face at home with our daily living set ups. Our experienced and courteous staff gladly provide all the information and assistance you require.

We work extensively with Physiotherapists, Occupational Therapists and Aged Care Clinicians so that we can find the right form of solution for you. Over the course of 30 years, we have provided solutions to community care with diligence and professionalism. Our customers' assurance comes from the fact that we help you find the most suitable solution.


We take great pride in the assembling and manufacturing of all of our products. We culminate to the highest standards in care, comfort for the patient and project cost efficiency.

Procurement & Distribution

Active Medical has developed a network of manufacturers in Australia and around the world – from Germany and Sweden to the United States of America and consults directly with industry professionals to ensure the products we provide responds to unique patient specifications.

The success of our business is supported by our ability to collaborate with manufacturers that not only provide the highest quality of products and services, but are aligned with our mission and commitment to our social, corporate and environmental responsibilities.

Our efficient distribution systems are underpinned by our strong relationships with manufacturers and two centralised warehouses with a capacity to hold up to 50,000 items. This ensures the fast turn around and on time delivery demanded by the health care industry and the needs of people in care. We ensure that all our products are ethically sourced and supplied and meet environmental standards.



*Disclaimer

Active Medical Supplies does not take responsibility for incorrect application of products. This publication is a guide only to assist in selecting products suitable to your clients. Please use in conjunction with professional advice.

CAPABILITY STATEMENT

Did you know...



We sanitize **all products** that leave our warehouses



We offer full product **demonstrations and video trials**



We offer **contact free delivery**



Our staff are **double vaxxed** with boosters




Our showroom is a **controlled environment** that we can also offer trials

Come in & have a chat!

Book a consultancy today!

Active Medical facilitates a showroom on site to replicate the home care environment. Our showrooms are used to replicate a residential care bedroom setting and allows us to offer unique demonstrations and training tailored to your needs.

We are located in Cannon Hill and Toowoomba with wheelchair accessibility. We showcase our extensive range of electric beds, wheelchairs, mobility, patient handling, safety equipment, bathroom solutions and accessories.



Cannon Hill

35 Southgate Ave, Cannon Hill QLD 4170



Toowoomba

11 Grange St Toowoomba City QLD 4350



Showroom Opening Hours

Monday to Friday | 8:00am to 4:00pm



Showroom Opening Hours

Monday to Friday | 8:00am to 4:00pm

*Products and information listed in this catalogue are subject to change. Website and digital information supercedes all print. For the most up to date information, please contact us directly via the details listed above.

BED, BATH & TOILET



Staebel Aluminium Shower Chair with Arms

This Shower Chair provides a comfortable seat for those unable to stand whilst showering. The armrests provide security when entering or exiting the shower, and assist the user on and off the seat.

- 3 piece padded vinyl seat
- Hand rail for added support while transferring
- Reversible padded backrest to suit every situation
- Welded frame



Staebel ExtraCare Shower Chair with Arms

The Extracare Shower Chair provides a comfortable seat for those unable to stand whilst showering. The armrests provide security when entering or exiting the shower, and assist the user on and off the seat.

- Strengthened frame for users needing extra room and weight capacity
- Seat height and angle adjusts to suit the individuals requirements
- One piece moulded supportive seat and backrest
- Small drainage slots allow water to escape and remove risk of entrapment

DVA

SPECIFICATIONS:

Product Code	1315J
DVA	B4002 - A
Seat Width	440 mm
Seat Depth	430 mm
Seat Height	400 - 550 mm
Bath Rim Clear	460 mm
Seat-Back Height	350 mm
Overall Width	570 mm
Overall Depth	540 mm
SWL	200 kg

SPECIFICATIONS:

Product Code	1320J
Seat Width	440 mm
Seat Depth	430 mm
Seat Height	400 - 550 mm
Width Between Arms	500 mm
Backrest Height	350 mm
Overall Width	620 mm
Overall Depth	570 mm
SWL	200 kg

BED, BATH & TOILET



Staebel Aluminium Over Toilet Frame

This aluminium over toilet frame is adjustable in height and fits over the toilet. This frame has been designed to assist the user who has difficulty sitting down or standing up by raising the toilet seat height.

- Aluminium construction for long life
- No risk of rust
- Removable splash guard included
- Lift up seat that removes for cleaning
- Ideal for those with hip or knee injuries or who require an increased seat height
- Unique stacking design saves space when storing more than one frame
- Multi D leg adjustment ensures stability on uneven surfaces
- Push button height adjustment
- Non slip guides on base of legs

MASS

SPECIFICATIONS:

Product Code	1755J
Width Between Arms	450 mm
Seat Height	475 - 625 mm
Overall Width	560 mm
Armrest Height	700 - 850 mm
Footprint (W x D)	560 x 590 mm
SWL	150 kg



Shower Stool with Plastic Seat

The aluminium shower stool with arms, provides security when entering or exiting the shower and assist the user on and off the seat.

- Aluminium construction for long life
- Plastic Seat
- Compact size for the smaller shower recess
- Adjustable seat height for the individual's requirements

SPECIFICATIONS:

Product Code	1565J
Seat Width	370 mm
Seat Depth	290 mm
Seat Height	520 - 670 mm
Width Between Arms	460 mm
Armrest Height (From Floor)	640 - 840 mm
Overall Width	550 mm
Overall Depth	370 mm
SWL	175 kg

OPTIONS AVAILABLE:

Shower Stool with Padded Seat	
Product Code	1567

BED, BATH & TOILET



Staubel ExtraCare Over Toilet Frame

This ACare Extra Care over toilet frame, one of four available sizes, features an Aluminium frame with armrests and plastic seat and splash guard.

- 500mm wide seat width
- Ideal for people with hip or knee injuries or those who need an increased seat height
- Bowl available as an accessory for commode use.
- Splash guard included
- 150kg Safe Weight Limit

SPECIFICATIONS:

Product Code	1750J
Width Between Arms	510 mm
Seat Width	500 mm
Seat Height	500 mm
Armrest Height	700 - 850 mm
SWL	150 kg



Staubel Shower Trolley

- Shower trolley with one-section top
- Central braking system
- Easy steering system (for ease of maneuvering trolley's front left castor can be locked in direction of steering)
- Side rails – foldable, height: 185 mm with safety release mechanism (Material - powder coated stainless steel)
- Metal base cover
- Water-resistant shower mattress (height: 270 mm)
- Mattress integrated drain with 645 – 2000 mm long flexible hose
- Water-resistant head pillow – semi - cylinder

SPECIFICATIONS:

Product Code	1209ST
Shower trolley top Width	650 mm
Shower trolley top Length	1900 mm
Shower trolley top Height	540 - 970 mm
Height with Castors	600 - 1030 mm
Trolley top inclination	0° - 5°
Length	2070 mm
Width	1070 mm
SWL	170 kg

BED, BATH & TOILET



Staebel Junior Shower Trolley

- Shower trolley with one-section trolley top
- Adjustable 5° trolley top inclination
- Hydraulic lowering/raising of trolley top via a console from 54 cm to 97 cm
- castors with diameter of 150 mm
- Central braking system
- Easy steering system (for ease of maneuvering trolley's front left castor can be locked in direction of steering)
- Side rails – foldable, height: 18,5 cm with safety release mechanism (material - powder coated stainless steel)
- Metal base cover
- Water-resistant shower mattress (height: 27 cm)
- Mattress integrated drain from 645 – 2000 mm long flexible hose
- Water-resistant head pillow – semi-cylinder
- Trolley weight: 79 kg

SPECIFICATIONS:

Product Code	1214ST
Shower trolley top Width	600 mm
Shower trolley top Length	1500 mm
Trolley top inclination	0° - 5°
Length	1670 mm
Width	760 mm
SWL	170 kg



Staebel Bariatric Shower Trolley

- Hydraulic lowering/raising of the stretcher top via a console lift
- Lowering/raising range: 540 mm to 970 mm
- Castors with a diameter of 150 mm
- Central braking easy steering system (for ease of maneuvering the trolley's front left castor can be locked in the direction of steering directional wheel)
- Side rails – foldable, height: 185mm with safety release mechanism (Material - powder coated stainless steel)
- Metal base cover
- Water resistant shower mattress (height: 270mm)
- Mattress integrated drain with a 645 up to 2000mm long flexible hose
- Water resistant head pillow – semi – cylinder

SPECIFICATIONS:

Product Code	1216ST
Shower trolley top Width	900 mm
Shower trolley top Length	2000 mm
Shower trolley top Height	540 - 970 mm
Height with Castors	600 - 1030 mm
Trolley top inclination	0° - 5°
Length	2170 mm
Width	1070 mm
SWL	250 kg

BED, BATH & TOILET



Staebel Bariatric Electric Shower Trolley

- Electric lowering/raising of the stretcher top via a console lift
- Lowering/raising range: 540 mm to 970 mm
- Castors with a diameter of 150 mm
- Central braking easy steering system (for ease of maneuvering the trolley's front left castor can be locked in the direction of steering directional wheel)
- Side rails – foldable, height: 185 mm with safety release mechanism (Material - powder coated stainless steel)
- Metal base cover
- Water resistant shower mattress (Height: 270 mm)
- Mattress integrated drain with a 645 up to 2000mm long flexible hose
- Water resistant head pillow semi cylinder
- Trolley weight: 110 kg

SPECIFICATIONS:

Product Code	1217ST
Shower trolley top Width	900 mm
Shower trolley top Length	2000 mm
Shower trolley top Height	540 - 970 mm
Height with Castors	600 - 1030 mm
Trolley top inclination	0° - 5°
Length	2170 mm
Width	1060 mm
SWL	250 kg

Staebel Shower Trolley - Electric with Backrest

- Shower trolley with two-section top
- Adjustable 5° trolley top inclination
- Electrically lowering/raising of trolley top via console from 54 cm to 97 cm
- Electrically adjustable head section inclination 45° (head section length: 70 cm)
- Castors with diameter of 150 mm
- Central braking system
- Easy steering system (for ease of maneuvering trolley's front left castor can be locked in direction of steering)
- Side rails – foldable, height: 18,5 cm with safety release mechanism (material - powder coated stainless steel)
- Metal base cover
- Water-resistant shower mattress (height: 27 cm)
- Mattress integrated drain with 645 – 2000 mm long flexible hose
- Water-resistant head pillow – semi-cylinder
- Trolley weight: 110 kg

SPECIFICATIONS:

Product Code	1211ST
Shower trolley top Width	650 mm
Shower trolley top Length	1900 mm
Trolley top inclination	0° - 5°
Length	2070 mm
Width	810 mm
SWL	170 kg

PATIENT HANDLING



Stael Patient Transfer

The Stael Patient Transfer Device is the perfect solution for safe transfer of a patient or resident within the home or facility. This is your all-in-one solution without the need for additional lifting straps and belts, and promoting independence. The width adjustable base allows good access to most seating, whilst also designed for good access beneath beds etc.



SPECIFICATIONS:

Product Code	212145
Seat Width	450 mm
Overall Width	640 mm
Overall Length	930 mm
SWL	150 kg



Stael Transfer Platform

Stael Transfer Platform is a manually-operated device designed as a platform for use by healthcare professionals to assist in the physical moving (transfer) of a patient in an upright (standing) position over a short distance. The device consists of a wheeled platform, a surrounding framework, kneepad and handlebar. It can include safety belt, to hold the person being transferred in the standing position.



SPECIFICATIONS:

Product Code	212150
Width	620 mm
Height	1180 mm
SWL	150 kg

PATIENT HANDLING



Staebel Transfer Belt

Staebel Transfer Platform Belt is used with the Staebel Transfer Platform when more support and stability is needed, for instance when the user's weight-bearing ability is questionable or varied. Staebel Transfer Platform Belt can also be used stand alone for sit-to-stand transfers.

- Anti-slip inner surface
- Sizes available, XSmall, Small, Medium, Large, XLarge & XXLarge
- Each size is colour coded for easy identification




SPECIFICATIONS:

Product Code	212152M
Size	Medium
Width	210 mm
SWL	170 kg

SIZES AVAILABLE:

Product Code	Size	SWL
212152XS	Extra Small	170 kg
212152S	Small	170 kg
212152L	Large	170 kg
212152XL	Extra Large	170 kg
212152XXL	Extra Extra Large	170 kg



Staebel Bed Transfer Assist

The Bed Transfer Assistant has been designed to assist in supporting the lower limbs during functional transfers. Simple and brilliant, our Bed Transfer Assistant allows one carer to safely raise a person's legs and swivel them around into bed and vice versa! Cotton upper for comfort.

- Raising legs and swivelling in a sitting position by a single carer



SPECIFICATIONS:

Product Code	210315
Width	600 mm
Length	640 mm
Thickness	30 mm
SWL	100 kg

PATIENT HANDLING



Staebel Disposable Rolling Slide

The Disposable Procedure Slide is a simple slide sheet that allows health professionals to easily move patients into the lithotomy position, ready for surgery. It is ideal for head-down procedures where the patient is positioned at an incline of up to 25 degrees. Those who work in Gynaecology will find the Disposable Procedure Slide hugely beneficial as it allows medical practitioners to:

- Slide a patient into a precise position (fine adjustments achievable), usually the Trendelenburg position
- Lock patient position with the sheet's unique and integral locking panels
- Position patients at inclines up to 25 degrees, facing head-down
- Adjust patient position mid-surgery
- Slide patients back up the operating table post-surgery
- Disposable slides help teams uphold excellent infection control
- Slippery material enables easy patient repositioning
- Integral locking pads ensure absolute safety and security
- Can remain under the patient
- 25 per pack

SPECIFICATIONS:

Product Code	210134
Width	1000 mm
Length	2000 mm



Staebel Knee Cushion

The Knee Cushion will provide a comfortable cushion for those who spend a long time kneeling as part of their job and/or find it painful to kneel on the floor. Our Knee Cushions can be used on any floor surface, inside or out, and will provide a clean and comfortable surface to kneel on whilst protecting and insulating your knees.

They are completely sealed and have been welded to ensure that they are easy and safe to wipe down. There is no need to send these Knee Cushions to the laundry or worry about cross-infection anymore.



SPECIFICATIONS:

Product Code	210250
Width	450 mm
Length	40 mm

PATIENT HANDLING



Stael ComfortCare Bed Top Sheet

The ComfortCare is an in-bed sliding system for the care of the long term or 24 hour patient. The amazing ComfortCare offers patients 100% comfort and carers the ability to manoeuvre and reposition them in all directions with ease and then secure that position when needed.

- Moves a patient in all directions
- The locking panel will hold the patient steady
- Offers patients comfort and dignity
- Patients can remain on this indefinitely, either directly or under a cotton bed sheet
- Carers can re-position and manoeuvre patients without the need to disturb or handle directly
- Reduces risk of injury to carers
- Enhances patient care and clinical outcomes
- Suitable for use with hoists
- Easy to use
- Polycotton system, with low-friction inner linings
- Two panel system enables up/down manoeuvres, turning, swivelling, sitting up safely, transferring and turning into prone position
- Integral locking panel means safe position can be secured, even in sitting
- Suitable for single carers
- Washable at 40°C to 80°C and low temperature tumble dry
- Spare panels available

SPECIFICATIONS:

Product Code	210295
Width	1400 mm
Length	2000 mm
SWL	200 kg

Stael ComfortCare Bed Base Sheet

The ComfortCare is an in-bed sliding system for the care of the long term or 24 hour patient. The amazing ComfortCare offers patients 100% comfort and carers the ability to manoeuvre and reposition them in all directions with ease and then secure that position when needed.

- Moves a patient in all directions
- The locking panel will hold the patient steady
- Offers patients comfort and dignity
- Patients can remain on this indefinitely, either directly or under a cotton bed sheet
- Carers can re-position and manoeuvre patients without the need to disturb or handle directly
- Reduces risk of injury to carers
- Enhances patient care and clinical outcomes
- Suitable for use with hoists
- Easy to use
- Polycotton system, with low-friction inner linings
- Two panel system enables up/down manoeuvres, turning, swivelling, sitting up safely, transferring and turning into prone position
- Integral locking panel means safe position can be secured, even in sitting
- Suitable for single carers
- Washable at 40°C to 80°C and low temperature tumble dry
- Spare panels available

SPECIFICATIONS:

Product Code	210296
Width	900 mm
Length	2000 mm

PATIENT HANDLING



Stabel – One Way Slide & Lock

The One Way Slide & Lock are designed to allow carers to help a user slide to the back of a chair or wheelchair and not be able to slide down.

Handles are positioned equally along the length of the product and carers can use the handles directly or attach extension straps to the handles, to prevent them having to bend down.

- To help people to sit back properly in chairs or wheelchairs without sliding back down. People independently, or with assistance by one or two carers depending upon need
- Prevent people sliding down in wheelchairs, chairs and beds
- Improves people's seating posture and comfort
- Allows people to be positioned further back in a chair
- Increases people's independence
- Can be used unassisted or with carer assistance



SPECIFICATIONS:

Product Code	210190
Width	800 mm
Length	650 mm



Stabel Perfect Positioner

Do you hoist people into their wheelchairs and find it difficult to re-position them so that they are sitting properly?

Originally designed to help reposition a tetraplegic gentleman in a moulded wheelchair, this unique handling device has been redesigned and now enables carers to easily and comfortably reposition the user (once hoisted into the wheelchair) accurately at the back of the chair and ensure their pelvic tilt is correct.

- Especially useful in moulded and specialised wheelchairs
- Correct sitting position in wheelchairs
- Accurate pelvic tilt
- One carer or two



SPECIFICATIONS:

Product Code	210340
Width	900 mm
Length	2000 mm

PATIENT HANDLING



Staebel In Bed Turner 4-way with Handles

InBed Turner 4-Way with Handles allows safe and easy repositioning for people in bed – enabling them to turn from side to side and move up and down, with or without the aid of a hoist. They are designed for more dependent users who need assistance in bed, often from a single carer, and are designed to be used on top of one of our InBed Turner Base Sheets. They do not work on their own.

- Can remain under the patient at all times.
- Great for those who need multi positioning in bed – up / down and turning.
- Reduces risk of tissue and nerve damage.
- Helps to reduce incidence of pressure sores developing.
- The combination of both silky smooth satin panels working with each other mean easy, quiet and safe movements.
- Cotton and polyester mix for great comfort and repeated washing.
- Washable to 90C.
- Four-way movement in all directions for patient care.
- Allows more dependent people to be moved in comfort and safety.
- Cotton and satin fibres make for comfortable in-bed experience.
- Choice of torso or full length models.
- Available with or without handles for attachment to hoists.
- Main sizes are all 200cm wide, allowing the sheet to be tucked in after repositioning, increasing patient security.



SPECIFICATIONS:

Product Code	210322
Width	2000 mm
Length	1400 mm
SWL	300 kg

SIZES AVAILABLE:

Product Code	210324
Width	2000 mm
Length	2000 mm
Height	13 mm

USE WITH:

Product Code	210321
Product Code	210323
Product Code	210325

PATIENT HANDLING




Staebel In Bed Turner Base Sheet

The InBed Turner Base Sheet is part of our renowned satin range for in-bed repositioning and patient care. All InBed Turner Base Sheets allow safe and easy repositioning for people in bed enabling them to turn easily from side to side but not slip down. Extra safety is provided by the striped cotton edges on all of our Base Sheets. They can be used on their own or combined with one of our Satin Top Sheets for those who need more in bed repositioning care.

- Reduces risk of tissue and nerve damage
- Helps reduce pressure sores developing
- No training required
- Easy turning and repositioning with little effort
- Great for those having difficulty turning in bed
- Washable
- Full length or midi length satin
- Central satin panel is smooth for ease of movement
- Anti-slip cotton polyester on edges reduces risk of slipping
- Cotton polyester for enhanced comfort and repeat washing
- Single and king options are fitted sheets. The double is a flat sheet.

USE WITH:

Product Code	210322
Product Code	210324




SPECIFICATIONS:

Product Code	210321
Size	Single
Width	2000 mm
Length	900 mm
SWL	300 kg

SIZES AVAILABLE:

Product Code	210323	210325
Size	King Single	Double
Width	2000 mm	2000 mm
Length	1050 mm	2300 mm
SWL	300 kg	300 kg

PATIENT HANDLING




Staebel Slideezi with Handles

The Slideezi with Handles are manufactured using an ultra-slippery material developed by pioneering fabric specialists Slideezi®. While one side of the slide sheet is luxuriously soft to offer maximum comfort and quietness, technically advanced fibres give the other side a slippery finish, ideal for repositioning a patient in a bed or trolley.

The Slideezi with Handles allow health professionals to effortlessly manoeuvre patients/clients into a virtually endless list of positions so that they can remain comfortable at all times. Whether patients are using a bed or a trolley, health professionals can use the slide sheets to turn patients into comfortable positions and also transfer them from one surface to another.

- Superb slideability on both sides for ease of use
- Easy to insert and remove
- Superb comfort for patient during repositioning
- Quiet and soft for improved patient experience
- Fabric is manufactured using advanced technical fibres
- Effortless and easy repositioning
- Used in pairs for best performance



SPECIFICATIONS:

Product Code	210432
Width	900 mm
Length	2000 mm
SWL	300 kg

Staebel Slideezi Support Straps - Pair

These Support Straps are used to reduce the risk of over-stretching and reaching in frequent supine transfers. No need to launder reducing the chance of losses.

A new concept in support straps, these can be wiped down and/or washed in soapy water, rinsed and then left to dry. Used to reduce the risk of over-stretching and reaching in frequent supine transfers, they simply clip onto the handles of the flat Slideezi with handles and other patient handling equipment.

- For use in supine transfers
- Reducing the risk of overreaching and stretching



SPECIFICATIONS:

Product Code	210064
Width	50 mm
Length	750 mm
SWL	100 kg

Deluxe Slidesheet

The ideal solution when transferring or bed repositioning patients across, up or down a bed. Reduces risk injury by assisting sliding, rather than lifting. Made from strong, non flammable, nylon fabric our slide sheets will not use their slippery quality or colour even after numerous washes. Can be machine washed to 50 degrees Celcius, tumble dried to 60 degrees Celcius and even ironed on a cool setting.

Available colours:

- Blue 2434B
- Green 2434G
- Red 2434R

SPECIFICATIONS:

Product Code	2434
Width	1450 mm
Length	2000 mm

PATIENT HANDLING



Staelbel Handy Legs

An aid to help lift a leg/legs into and out of bed and/or cars. Can be used independently or by a single carer.

It can be used in most situations where you need to raise your legs into bed or onto the couch/sofa but you must have sufficient strength to be able to sit unaided on the edge of the bed.

- The Handy Legs can also be used to help one leg only



Staelbel Support Walk Belt

A supportive Belt to help with transferring or walking. The ergonomically designed handles enable the care giver to hold the patient close without the need to grab clothes.

Ranging from Small to Large (adult sizes) this wonderful range of comfortable and supportive belts means that getting close to your patient has never been easier.

Simply place the belt around the person's waist, adjust to fit and then buckle. The ergonomically designed handles are perfectly located for user comfort and functionality. The belt is also washable.

- Helping a person from sitting to standing and vice versa
- Helping whilst transferring on a sitting transfer board
- Comfortable to wear
- Give a feeling of added security
- Easily adjusted to fit



Staelbel Transfer Board Slide

The Transfer Board Slide is a simple low-cost device that allows for a sitting transfer across the board without risk of shearing/induced tissue damage.

This is a wonderful little product to use with any sitting transfer board and makes sure that the slide across is what it should be – effortless. Just pop it on top of your board, tuck it under your bottom a little and sliding from one surface to another will be plain sailing.

Simple and highly effective, this Transfer Board Slide is only as wide as the board itself and has a panel of antislip that grips the top of the board.

- Transferring from one sitting surface to another
- Helping people move, both independently or with assistance



SPECIFICATIONS:

Product Code	210050
Width	250 mm
Length	280 mm
Thickness	50 mm

SPECIFICATIONS:

Product Code	210082
Width	250 mm
Length	730 - 910 mm

SIZES AVAILABLE:

Product Code	210083	210084
Width	250 mm	250 mm
Length	880 - 1040 mm	1040 - 1170 mm



SPECIFICATIONS:

Product Code	210360
Width	300 mm
Length	700 mm

PATIENT HANDLING



Staelbel Anti Shear Mattress Cover

A simple answer to the problem of constantly putting slide sheets in and out – lightweight and easy to use. You will love the Anti Shear Mattress Cover’s simplicity as it makes it so much easier to adjust someone’s position without disturbing them and really reduces the risk of injury for staff.

Thinner than a bed sheet, this product fits over the mattress and does not impinge on any pressure reducing qualities the mattress may have. You can sit the client/patient up in bed without the risk of them slipping down, and it only becomes fully slippery when you introduce a flat slide sheet on top to perform the manoeuvre that needs to be performed.

Just place the Anti Shear Mattress Cover over the bed mattress and, if you want to, put the normal bed sheet on top, tucking it in. When you need to reposition someone, untuck the bed sheet and unroll a flat slide sheet in between the Anti Shear Mattress Cover and bed sheet (without disturbing the patient) and use this to achieve a new position. Then slip out the slide sheet easily, retuck the bed sheet and your patient will be comfortable.

- Specialist technical fabric is fully breathable, lightweight and quiet
- Fabric is thinner than a bed sheet and will not impinge on pressure reducing qualities of mattresses
- Light-weight and compact
- Washable at 40°C to 74°C
- When using the base sheet to lie on, independent people can turn more easily in bed
- Breathable fabric means that people will not sweat or get too hot – as can happen with other systems

SPECIFICATIONS:

Product Code	210290
Width	900 mm
Length	2000 mm



Staelbel Sure Foot Pad

The Sure Foot Pad is a really helpful, inexpensive answer to the problem of feet slipping, on the floor or in bed.

Suitable for use on the floor or in bed, this Sure Foot Pad is designed as a safety aid that helps to prevent feet slipping. It can be used to stand on when sitting or rising from a chair, or it can be used with a slide sheet to aid movement in bed or with a Slide & Lock to aid movement in a chair. A helpful aid for more independent living, it is hard wearing and easy to clean.

Simply place the foot pad on the surface you need to gain more grip on and your feet won’t slip. You will love its simplicity because it gives you so much independence and is also a really useful aid for rehabilitation.

Ideal for use with small slide sheets to slide further up the bed or on the floor, on its own or with products such as Slide & Locks. In fact there are so many places that you can use this product, some people even have two; those for the floor and those for use in bed!

- Gaining more floor grip with your feet
- Helping people be more independent
- A rehabilitation aid



SPECIFICATIONS:

Product Code	210001
Width	400 mm
Length	480 mm

PATIENT HANDLING



Staebel Easy Hands

These silicon coated nylon slippy gloves are designed for carers and can be used to make patients/service users much more comfortable in many ways. For example: skin smoothing, positional adjustments and inserting slings etc.

- Technically designed fabric means the gloves slide easily on the outside but will not slide off the users hands
- Easy to slide under patients for improved patient care
- Can be used for a multitude of patient care tasks:
 - Skin smoothing
 - Tucking in sheets
 - Adjusting patient clothing
 - Moving pillows
 - Inserting slings etc
- Supplied in pairs
- Slippy outer and grippy inner material
- Fits to elbow



SPECIFICATIONS:

Product Code	210339
Width	200 mm
Length	520 mm



Staebel Easy Dress Arms

A pair of long slippy sleeves that can slide over a person's arm and ease the putting on of upper body garments.

For best practice: The sleeves work especially well with a flat slide sheet underneath the bed bound person to facilitate the slide.

- Help carers put a shirt, blouse or jumper on a bed bound person/child or someone who requires a lot of assistance



SPECIFICATIONS:

Product Code	210159
Size	Medium
Width	1070 mm
Length	280 mm

SIZES AVAILABLE:

Product Code	210161
Size	Large
Width	270 mm
Length	1090 mm

PATIENT HANDLING



Staebel Easy Dress Legs

A pair of long slippy leg sleeves that can slide over a person's legs and ease the putting on of lower half garments.

For best practice: The sleeves work especially well with a flat slide sheet underneath the bed bound person to facilitate the slide.

- Help carers put trousers/pyjamas on a bed bound person/child or someone who requires a lot of assistance



SPECIFICATIONS:

Product Code	210160
Size	Medium
Width	340 mm
Length	1200 mm

SIZES AVAILABLE:

Product Code	210162
Size	Large
Width	360 mm
Length	1300 mm




Staebel In Bed Sit Support

This easy-to-attach bed ladder allows patients to maintain control as they sit up and lie down. Patients who would otherwise struggle to manoeuvre themselves gain a whole new level of independence with the bed rope ladder, allowing them to reposition themselves safely and easily.

You will love how simple this multicoloured functional webbing ladder is to attach. No knots to ties or worry about. Just loop it through itself and then use the different coloured rungs to help you as you sit up or lie down. Ideal if you want to maintain your independence and just need a little more help or if you are going through a process of rehabilitation.

- To help people sit up/lie down independently and in a controlled manner



SPECIFICATIONS:

Product Code	210440
Width	160 mm
Length	1880 mm
SWL	300 kg



Staebel Easy Dress Kit

Getting dressed in bed has never been easier – the Easy Dress Kit enables dignified dressing for patients who need carer assistance.

This Easy Dress Kit is a complete set of in-bed dressing aids that enables one or two carers to dress a patient/client fully in bed. Where dressing is difficult or uncomfortable, they enable easy dressing and promote patient dignity.

- 2 x Easy Dress Legs
- 2 x Easy Dress Arms
- 1 x Flat Sheet (folded in a handy storage bag)
- Easy dressing aid for carers to clothe a person in bed
- Suitable as an aid for an individual who can clothe themselves but require a little help
- Super Slippy fabric for extra ease
- Available in Washable and disposable material
- Available in Standard (one size fits all)
- Washable 40° to 74° and low tumble dry
- Replacement parts available if lost



SPECIFICATIONS:

Product Code	210153 (Complete Kit)
--------------	-----------------------

PATIENT HANDLING



Staebel Emergency Lifting Kit

An answer to the problem of handling a collapsed person in a number of situations – in a tight space, on the floor or in bed.

This specially designed pack contains all the equipment you will need, especially where there is no mechanical aid available. For example:

- Arrest in bed or in a chair
- Collapsed in an awkward corner or tight space.
- Collapsed on the floor and no mechanical hoist or lifting mattress available to lift off the floor.

Includes:

- 2 x Slideezi + Handles
- 1 x Rolling Slide
- 1 x Lifting Sheet
- 2 x Support Straps

Staebel Lifting Sheet

The Lifting Sheet offers the ultimate falls management solution for transporting an injured person from the site of a car accident or other emergency.

No stretcher hoist? No way to lift a person safely? Look after your staff and patients by providing them with this inexpensive lifting sheet, designed to protect both patients and health professionals during lifting procedures.

The Lifting Sheet is a lightweight falls management solution, ideal for use when health professionals need to raise someone from the floor in a horizontal position. Designed in conjunction with the Guidelines for Safer Handling during Resuscitation in Healthcare Settings, the lifting sheet is perfect for post-resuscitation transfers where a stretcher hoist is not available to move the patient from the floor onto a bed or trolley.

- Lightweight solution when no stretcher hoist is available
- Robust design supports patients weighing up to 150kg
- Washable at 74°C
- Waterproof
- Fire-retardant
- Can remain under patient

Staebel Safemove

The Safemove has been developed with practitioners and allows radiotherapists to position their patients far more precisely than has been possible until now.

Positioning patients for laser precision has always been difficult and achieving that specific position can be uncomfortable and awkward for the patient and extremely hard, tiring and heavy work for staff. It can result in potential skin damage to patients and injury risk to staff. However, by placing the Safemove under the patient (breast or pelvic areas), staff can use the straight and/or curved handles to much more easily re-position their patients to line up precisely with the lasers.

- No need to manually handle patients
- Reduced risk of tissue damage
- Less stress and discomfort for the patient
- Less strain and risk of injury to the staff

SPECIFICATIONS:

Product Code 210140 (Complete Kit)

SPECIFICATIONS:

Product Code 210138
SWL 150 kg

SPECIFICATIONS:

Product Code 210317
Width 370 mm
Length 620 mm

PATIENT HANDLING



Bed Rail Protector

This Bed Rail is made from a polyurethane material, coupled with a bioelastic incontinence cover, making it the ultimate answer for longevity and effectiveness.

No zips, clips, or flimsy foam means a zero risk of skin tears or wounds caused by the patient/resident hitting the side guards. Simply push it down onto the top rail and it is in place.

No time consuming, technical attachment points that foul up with the side guards.

- Skin friendly
- Anti Microbial safety
- Easy to clean and to disinfect
- Easy to use and install
- 2 x protectors can be used together to cover the complete length of a 4-part side guard or full length side rails.
- Sold singularly



**FOR USE WITH WIBO
BED SIDE GUARD**




SPECIFICATIONS:

Product Code	3983
Width	130 mm
Length	450 mm
Height	800 mm

Positioning And Support Wedge

This foam positioning support wedge is used to place a user in a side position in order to relieve the sacrum, which is an extra sensitive area. For the side position, two support wedges are used together. The support wedge is easy to use and easy to clean. It is supplied with a zipped, liquid resistant and vapour permeable cover. The cover can be wiped or machine washed.

- Soft comfortable foam, stretch vinyl covered



SPECIFICATIONS:

Product Code	3689
Width	260 mm
Length	700 mm
Height	50 - 150 mm

Staebel Crash Mat

- Foldable Falls Impact Reducing Mat.
- Easy to clean
- Made from heavy duty latex free material
- Chamfered edges reduce tripping hazard
- Can be easily folded for storage or transportation

SPECIFICATIONS:


Product Code	9822
Width	910 mm
Length	1760 mm
Thickness	40 mm

Pressure Care

 **LOW RISK**

Stabel Zonecare Mattress	28
Deluxe Foam Cushion	28
Stabel Bariatric Cushion	29

PRESSURE CARE



Stael Zonecare Mattress

The Zonecare foam mattress is a high quality foam mattress featuring a profile-cut pressure relieving foam insert. It provides pressure reduction and assists in the treatment of pressure ulcers, whilst enhancing patient comfort levels in homecare and residential care environments. The Zonecare's specialised inner cover aids in shear reduction on a profiling bed.


- Unique re-reinforced edges designed with longevity in mind



Deluxe Foam Cushion

A comfortable cushion for the low risk user, the deluxe pressure relief foam cushion flows under pressure to provide a fluid like support.

- Suitable for low risk clients
- 40mm top layer of memory foam. 50mm base of standard foam
- 2-way stretch top cover
- Zipped cover included with waterfall zip



SPECIFICATIONS:

Product Code	3728Z
Width	850 mm
Length	1980 mm
Thickness	150 mm
SWL	25 - 150 kg

SIZES AVAILABLE:

Product Code	3729Z
Width	950 mm
Length	1980 mm
SWL	25 - 150 kg

SPECIFICATIONS:

Product Code	3220
Width	450 mm
Length	400 mm
Thickness	90 mm
SWL	100 kg

SIZES AVAILABLE:

Product Code	Width	Length	Thickness
3218	400 mm	400 mm	90 mm
3225	450 mm	450 mm	90 mm
3221	500 mm	450 mm	90 mm
3226	500 mm	500 mm	90 mm
3224	400 mm	510 mm	90 mm
3223	460 mm	560 mm	90 mm
3227	460 mm	610 mm	90 mm

PRESSURE CARE



Staebel Bariatric Cushion

A comfortable bariatric cushion for the low risk user, Helps to eliminate pressure points by transferring support uniformly. 2 foam layers, top layer of 50mm memory foam for comfort and is supported by 100mm of standard foam to ensure no bottoming out. Stretch cover with waterproof zip for easy cleaning.

- Suitable for low risk clients
- Zipped, removable cover included
- 260kg SWL
- Other sizes available upon request

SPECIFICATIONS:

Product Code	3270
Length	660 mm
Width	520 mm
Thickness	150 mm
SWL	260 kg

SIZES AVAILABLE:

Product Code	3271	3272
Length	710 mm	610 mm
Width	550 mm	460 mm
Thickness	150 mm	150 mm
SWL	260 kg	260 kg

WE ARE
Here to Help
TOOWOOMBA

If you would like to see any products in person **book a trial today** at our **brand new Toowoomba showroom!**



11 Grange St
Toowoomba City QLD 4350



07 3518 5974



Seating & Mobility

Staelbel Emperor Wheelchair	31
Staelbel Emperor Transit Wheelchair	31
Staelbel Liberation Ultra - Light Wheelchair	32
Thrifty Self Propelled Wheelchair	32
Staelbel Sovereign Wheelchair	33
Staelbel Patient Transport Chair	34
Staelbel Unilite Forearm Walker	34
Staelbel Standard Bedside Commode	34

SEATING



Staebel Emperor Wheelchair

The Staebel Emperor is a comfortable, lightweight aluminium wheelchair, with a really helpful range of high quality features. With its sleek modern frame and lightweight foldable design, the wheelchair is perfect for travel. All features of this chair have been designed to make life easier for the user.

- Foldable backrest
- Lap belt
- Heel strap
- Height adjustable armrests and footplates
- Removable swing-away footrest
- Available in orange or blue.

SPECIFICATIONS:

Product Code	40411
Seat Width	450 mm
Overall Width	670 mm
Overall Height	950 mm
Product Weight	12 kg
SWL	115 kg

SIZES AVAILABLE:

Product Code	40412
Seat Width	500 mm
Overall Width	720 mm
Overall Height	950 mm
Product Weight	15 kg
SWL	135 kg

ALSO AVAILABLE:

Staebel Emperor Wheelchair - Orange	
Product Code	40416W



Staebel Emperor Transit Wheelchair

The Staebel Emperor is a comfortable, lightweight aluminium wheelchair, with a really helpful range of high quality features. With its sleek modern frame and lightweight foldable design, the wheelchair is perfect for travel. All features of this chair have been designed to make life easier for the user.

- Foldable backrest
- Lap belt
- Heel strap
- Height adjustable armrests and footplates
- Removable swing-away footrest
- Available in orange or blue.

SPECIFICATIONS:

Product Code	40413
Seat Width	450 mm
Overall Width	670 mm
Overall Height	950 mm
Product Weight	12 kg
SWL	115 kg

SIZES AVAILABLE:

Product Code	40414
Seat Width	500 mm
Overall Width	720 mm
Overall Height	950 mm
Product Weight	15 kg
SWL	135 kg

ALSO AVAILABLE:

Staebel Emperor Wheelchair Transit - Orange	
Product Code	40417

SEATING



Staebel Liberation Ultra - Light Wheelchair

The Staebel Liberation Lightweight Wheelchair is a great travel companion. It offers attendant and parking brakes, ultra lightweight frame, folds down in seconds, and has foldable and removable legrests.

- Folding backrest
- Safety Brake and Assistant Brake
- Folding Arm Rest and Foot Rest
- Includes carry case



Thrifty Wheelchair

The Thrifty is a sturdy and reliable standard weight wheelchair, with removable armrests and leg rests, built for individuals, nursing homes and non-prescription requirements. This wheelchair meets Australian Standards for Manual Wheelchairs in Product Performance. All of these features combine to make the Thrifty of unbeatable value in standard wheelchairs.

- Swing away, nonremovable armrests
- Swing away, detachable legrests
- 24" rear wheels with urethane tyres (maintenance free)
- Australian Standard AS/NZS 3695.1:2011 approved



SPECIFICATIONS:

Product Code	40400
Seat Width	450 mm
Overall Width	540 mm
Overall Height	920 mm
Weight	7 kg
SWL	100 kg

SPECIFICATIONS:

Product Code	4246
Seat Width	450 mm
Seat Depth	400 mm
Overall Width	650 mm
Overall Depth	1080 mm
Overall Height	880 mm
Fold Width	300 mm
Fold Length	790 mm
SWL	115 kg

SEATING



Adjustable & removable legrests



Drum Brake



Adjustable armrest height & length



Seat depth, height & angle adjustment

Staebel Sovereign Wheelchair

The Sovereign Wheelchair can be adjusted to suit almost any user. With much consideration in the design from the perspective of OTs, this chair is suited to the vast range of needs required by our clients. Rather than numerous script forms for various users, this chair meets the needs of almost any patient. All wheelchair parts suit any adjustments made - no additional parts required.

- Adjustable & removable leg rest's
- Flip back & removable armrests
- Adjustable armrest height & length
- Seat depth, height & angle adjustment
- Backrest height, tension & angle adjustment

SPECIFICATIONS:

Product Code	4340C - 43	4340C - 46	4340C - 49	4352C - 52	4352C - 55	4352C - 58
Seat Width	430 mm	460 mm	490 mm	520 mm	550 mm	580 mm
Anterior Seat Height	420 mm	445 mm	450 mm	470 mm	475 mm	500 mm
Posterior Seat Height	410 mm	430 mm	450 mm	470 mm	490 mm	510 mm
Seat Depth	400 mm	425 mm	450 mm	475 mm	475 mm	500 mm
Backrest Height (2cm increments)	400 - 460 mm	400 - 460 mm	400 - 460 mm	400 - 460 mm	400 - 460 mm	400 - 460 mm
Backrest Angle (5° increments)	69° - 94°	69° - 94°	69° - 94°	69° - 94°	69° - 94°	69° - 94°
Lower Leg Length (2cm increments)	400 - 480 mm	400 - 480 mm	400 - 480 mm	400 - 480 mm	400 - 480 mm	400 - 480 mm
Armrest Height (2cm increments)	180 - 260 mm	180 - 260 mm	180 - 260 mm	180 - 260 mm	180 - 260 mm	180 - 260 mm
Armrest Length (4cm increments)	330 - 450 mm	330 - 450 mm	330 - 450 mm	330 - 450 mm	330 - 450 mm	330 - 450 mm
Footplate Width	105 - 155 mm	105 - 155 mm	105 - 155 mm	105 - 155 mm	105 - 155 mm	105 - 155 mm
Gross Weight	17 kg	17 kg	17 kg	21 kg	21 kg	
SWL	125 kg	125 kg	125 kg	170 kg	170 kg	205 kg

SEATING



Stael Patient Transport Chair

Stael Patient Transport Chair is a handy, easy to use patient transport chair with flip-up armrests, generous 52cm side open seat and a retractable foot rest.

- 200kg SWL
- PVC document holder
- Adjustable braking system
- Silicone sealed seat for infection control
- HPL footplate
- UV protected Powdercoating
- Dual front castors for easy maneuvering.

SPECIFICATIONS:

Product Code	4750
Width	660 mm
Length	870 mm
Height	1090 mm
Seat Height	620 mm
Seat Width	520 mm
Seat Depth	430 mm
Weight	28.5 kg
SWL	200 kg



Stael Unilite Forearm Walker

The Unilite Forearm Walker is designed to support and tutor the user in an upright gait position. The very stable base gives added confidence and security to the user. Features include padded gutter supports that are adjustable in the forearm length, height and width. The popular folding design allows the Unilite to fold easily for transport and storage. This model is fitted with four large wheels and lockable cable brakes.

- Height adjustable ergonomic shaped hand grips
- Powder coated steel construction
- Lockable cable type brakes
- 200mm castors

SPECIFICATIONS:

Product Code	6742
DVA	BRO213
Seat Width	400 mm
Seat Height	600 mm
Handle Height	890 - 1100 mm
Walker Weight	14 kg
SWL	100 kg



Stael Standard Bedside Commode

The standard bedside commode offers padded armrests for comfort and features adjustable height, a removable padded seat cover and a plastic toilet seat with bowl and lid.

- Value for money
- Stackable design

DVA

SPECIFICATIONS:

Product Code	7952PA
DVA	B1060
Overall Width	575 mm
Overall Height	940 - 1090 mm
Overall Depth	600 mm
Seat Width	440 mm
Seat Height	500 - 650 mm
Seat Depth	440 mm
Back Height	475 mm
Toilet Seat	475 - 625 mm
SWL	125 kg



SEE SOMETHING YOU LIKE? **LET'S HAVE A CHAT!**

Book your consultancy today with one of our friendly team members. See below information for ways you can get in contact....



1800 267 267



Showroom Opening Hours

Monday to Friday | 8:00am to 4:00pm



www.activemedicalsupplies.com.au



Cannon Hill

35 Southgate Ave, Cannon Hill QLD 4170



hello@activemedicalsupplies.com.au



Toowoomba

11 Grange St Toowoomba QLD 4350

activeMEDICAL

S U P P L I E S > S O L U T I O N S > S U P P O R T

P. 1800 267 267 | **F.** 07 3899 8825

E. hello@activemedicalsupplies.com.au

www.activemedicalsupplies.com.au