

Bath, Shower & Toilet Portfolio

Issue 1

Contents

Capability Statement	4
Toilet Surrounds & Over Toilet Frames	6
Shower Chairs & Stools	9
Transfer Benches	13
Bedside Commodes	16
Shower Trolleys & Tables	18
Commodes	23
Juvo Commode Accessories	33
Shower Commode Seats	37
Ancillary Products	39

CAPABILITY STATEMENT

Active Medical Supplies is a pioneer in medical device and equipment manufacturing. The experienced team at Active Medical understand both professionally and personally, the importance of maintaining independence within the comfort of your own home or the community including educational settings. A lot of the people we help are still independent, however, just need assistive solutions as result of an injury, Hospital visit, Healthcare professional's advice or to help loved ones with their daily life.

With over 30 years of experience helping with solutions in community care, we understand sometimes that daily tasks may become difficult. Often things that we have taken for granted like getting out of bed, lifting yourself from a chair or just even moving around the house, inherit more effort than they have before.

For those who come into our showroom, we can replicate the scenarios you would face at home with our daily living set ups. Our experienced and courteous staff gladly provide all the information and assistance you require.

We work extensively with Physiotherapists, Occupational Therapists and Aged Care Clinicians so that we can find the right form of solution for you. Over the course of 30 years, we have provided solutions to community care with diligence and professionalism. Our customers' assurance comes from the fact that we help you find the most suitable solution.


We take great pride in the assembling and manufacturing of all of our products. We culminate to the highest standards in care, comfort for the patient and project cost efficiency.

Procurement & Distribution

Active Medical has developed a network of manufacturers in Australia and around the world – from Germany and Sweden to the United States of America and consults directly with industry professionals to ensure the products we provide responds to unique patient specifications.

The success of our business is supported by our ability to collaborate with manufacturers that not only provide the highest quality of products and services, but are aligned with our mission and commitment to our social, corporate and environmental responsibilities.

Our efficient distribution systems are underpinned by our strong relationships with manufacturers and two centralised warehouses with a capacity to hold up to 50,000 items. This ensures the fast turn around and on time delivery demanded by the health care industry and the needs of people in care. We ensure that all our products are ethically sourced and supplied and meet environmental standards.



*Disclaimer

Active Medical Supplies does not take responsibility for incorrect application of products. This publication is a guide only to assist in selecting products suitable to your clients. Please use in conjunction with professional advice.

CAPABILITY STATEMENT

Did you know...



We sanitize **all products** that leave our warehouses



We offer full product **demonstrations and video trials**



We offer **contact free delivery**



Our staff are **double vaxxed** with boosters




Our showroom is a **controlled environment** that we can also offer trials

Come in & have a chat!

Book a consultancy today!

Active Medical facilitates a showroom on site to replicate the home care environment. Our showrooms are used to replicate a residential care bedroom setting and allows us to offer unique demonstrations and training tailored to your needs.

We are located in Cannon Hill and Toowoomba with wheelchair accessibility. We showcase our extensive range of electric beds, wheelchairs, mobility, patient handling, safety equipment, bathroom solutions and accessories.



Cannon Hill

35 Southgate Ave, Cannon Hill QLD 4170



Toowoomba

11 Grange St Toowoomba City QLD 4350



Showroom Opening Hours

Monday to Friday | 8:00am to 4:00pm



Showroom Opening Hours

Monday to Friday | 8:00am to 4:00pm

Toilet Surrounds & Over Toilet Frames

Toilet Surround	7
Aluminium Over Toilet Frame	7
Folding Aluminium Over Toilet Frame	8
Staebel Extra Care Over Toilet Frame	8

TOILET SURROUNDS & OVER TOILET FRAMES



Toilet Surround

This adjustable height frame fits around the toilet and provides assistance in getting on and off the toilet. The toilet surround is manufactured from mild steel, which is zinc treated and epoxy coated making it both strong and durable.

- Ergonomic shaped side frame makes it easy to stand
- Non-slip feet
- Push button for height adjustment

SPECIFICATIONS:

Product Code	1900
Arm Height	625 - 775mm
Width Between Arms	460 mm
Overall Width	540 mm
Overall Depth	490 mm
SWL	125 kg

SIZES AVAILABLE:

Product Code	1901
Width Between Arms	540 mm
SWL	130 kg



Aluminium Over Toilet Frame

This aluminium over toilet frame is adjustable in height and fits over the toilet. This frame has been designed to assist the user who has difficulty sitting down or standing up by raising the toilet seat height.

- Aluminium construction for long life
- No risk of rust
- Removable splash guard included
- Lift up seat that removes for cleaning
- Ideal for those with hip or knee injuries or who require an increased seat height
- Unique stacking design saves space when storing more than one frame
- Multi D leg adjustment ensures stability on uneven surfaces
- Push button height adjustment
- Non-slip guides on base of legs
- Bowl available as an accessory for commode use

SPECIFICATIONS:

Product Code	1755J
Width Between Arms	450 mm
Seat Height	475 - 625 mm
Overall Width	560 mm
Armrest Height	700 - 850 mm
Footprint (W x D)	560 x 590 mm
SWL	150 kg

TOILET SURROUNDS & OVER TOILET FRAMES



Folding Aluminium Over Toilet Frame

This Over Toilet frame is ideal for those who need to remove or transport it often.

- Powder coated aluminium frame - light to carry
- Seat is able to be removed for replacement and/or cleaning
- Adjustable seat height to suit the individual's requirements



SPECIFICATIONS:

Product Code	1760
Width Between Arms	460 mm
Seat Height	480 - 630 mm
Overall Height	680 mm
Armrest Height	680 - 830 mm
Footprint (WxD)	530 x 510 mm
SWL	90 kg



Staebel Extra Care Over Toilet Frame

This Staebel Extra Care over toilet frame, one of four available sizes, features an Aluminium frame with armrests and plastic seat and splash guard.

- 500mm wide seat width
- Ideal for people with hip or knee injuries or those who need an increased seat height
- Bowl available as an accessory for commode use.
- Splash guard included
- 150kg Safe Weight Limit

SPECIFICATIONS:

Product Code	1750J
Width Between Arms	500 - 600 mm
Seat Width	500 mm
Seat Height	500 mm
Armrest Height	700 - 850 mm
SWL	150 kg

Shower Chairs & Stools

Aluminium Shower Chair	10
Staebel ExtraCare Shower Chair with Arms	10
Shower Chair with Padded Seat	11
Folding Shower Chair	11
Shower Stool with Arms & Back	12
Padded Shower Stool	12
Shower Stool With Plastic Seat	12

SHOWER CHAIRS & STOOLS



Aluminium Shower Chair

This Shower Chair provides a comfortable seat for those unable to stand whilst showering. The armrests provide security when entering or exiting the shower, and assist the user on and off the seat.

- Seat height and angle adjusts to suit the individuals requirements
- Aluminium construction for long life
- One piece moulded supportive seat and backrest
- Small drainage slots allow water to escape and remove risk of entrapment
- Chair design is stackable for space saving

SPECIFICATIONS:

Product Code	1315J
DVA	B4002 - A
Seat Width	440 mm
Seat Depth	430 mm
Seat Height	400 - 550 mm
Bath Rim Clear	460 mm
Seat-Back Height	350 mm
Overall Width	570 mm
Overall Depth	540 mm
SWL	200 kg



Staebel ExtraCare Shower Chair with Arms

The Extracare Shower Chair provides a comfortable seat for those unable to stand whilst showering. The armrests provide security when entering or exiting the shower, and assist the user on and off the seat.

- Strengthened frame for users needing extra room and weight capacity
- Seat height and angle adjusts to suit the individuals requirements
- One piece moulded supportive seat and backrest
- Small drainage slots allow water to escape and remove risk of entrapment

SPECIFICATIONS:

Product Code	1320J
Seat Width	440 mm
Seat Depth	430 mm
Seat Height	400 - 550 mm
Width Between Arms	500 mm
Backrest Height	350 mm
Overall Width	620 mm
Overall Depth	570 mm
SWL	200 kg

SHOWER CHAIRS & STOOLS



Shower Chair with Padded Seat

- 600mm wide seat and reinforced frame with 520mm seat height
- Padded seat and curved backrest provide comfort and postural support
- Powder coated steel frame
- Zinc treated to prevent rusting
- Made in Australia
- 200kg SWL
- Suitable for Bariatric Users

SPECIFICATIONS:

Product Code	1579
Seat Width	600 mm
Seat Height	520 mm
SWL	200 kg



Folding Shower Chair

The Folding Shower Chair, with the ability to fold the frame, is the solution to those needing a shower chair but live in multiple locations or travel. It has a moulded plastic seat and back with a rust proof epoxy coated aluminium frame.

- Non slip textured surface
- Front cut out for personal hygiene
- Adjustable seat height to suit the individual's needs

SPECIFICATIONS:

Product Code	1350
Seat Width	410 mm
Seat Depth	400 mm
Seat Height	455 - 520 mm
Width Between Arms	470 mm
Overall Width	530 mm
Overall Depth	455 mm
SWL	100 kg

SHOWER CHAIRS & STOOLS



Shower Stool with Arms & Back

The shower stool with arms and backrest, provides security when entering or exiting the shower and assist the user on and off the seat.

- Compact size for the smaller shower recess
- Padded seat and backrest for user comfort
- Adjustable seat height for the individual's requirements

SPECIFICATIONS:

Product Code	1569
Seat Width	370 mm
Seat Depth	290 mm
Seat Height	520 - 670 mm
Width Between Arms	640 - 840 mm
Armrest Height	640 - 840 mm
Overall Width	570 mm
Overall Depth	450 mm
SWL	125 kg

Padded Shower Stool

The shower stool with arms, provides security when entering or exiting the shower and assist the user on and off the seat.

- Compact size for the smaller shower recess
- Padded seat for user comfort
- Adjustable seat height for the individual's requirements
- Aluminium frame for longevity

SPECIFICATIONS:

Product Code	1567
Seat Width	370 mm
Seat Depth	290 mm
Seat Height	520 - 670 mm
Width Between Arms	460 mm
Armrest Height	640 - 840 mm
Overall Width	550 mm
Overall Depth	370 mm
SWL	125 kg

Shower Stool with Plastic Seat

The aluminium shower stool with arms, provides security when entering or exiting the shower and assist the user on and off the seat.

- Aluminium construction for long life
- Plastic Seat
- Compact size for the smaller shower recess
- Adjustable seat height for the individual's requirements

SPECIFICATIONS:

Product Code	1565J
Seat Width	370 mm
Seat Depth	290 mm
Seat Height	520 - 670 mm
Width Between Arms	460 mm
Armrest Height	640 - 840 mm
Overall Width	550 mm
Overall Depth	370 mm
SWL	175 kg

OPTIONS AVAILABLE:

Shower Stool with Padded Seat

Product Code	1567
--------------	------

TRANSFER BENCHES



Standard Transfer Bench

This Shower Chair provides a comfortable seat for those unable to stand whilst showering. The armrests provide security when entering or exiting the shower, and assist the user on and off the seat.

- 3 piece padded vinyl seat
- Hand rail for added support while transferring
- Reversible padded backrest to suit every situation
- Welded frame
- Leg Extensions available

SPECIFICATIONS:

Product Code	1018J
DVA	B5910
Seat Width	670 mm
Seat Depth	410 mm
Seat Height	445 - 570 mm
Bath Rim Clear	370 - 495 mm
Back Height	410 mm
Overall Width	790 mm
Overall Depth	540 mm
SWL	140 kg



Padded Transfer Bench

Designed for individuals with restricted mobility, the padded Transfer Bench provides safety and independence while transferring in and out of the bath tub. This durable aluminium transfer bench is adjustable in height.

- Padded seat and curved backrest provide comfort and postural support
- Powder coated steel frame
- Zinc treated to prevent rusting
- Leg Extensions available.

SPECIFICATIONS:

Product Code	1018 - 1J
Seat Width	670 mm
Seat Depth	410 mm
Seat Height	460 - 585 mm
Bath Rim Clear	370 - 495 mm
Back Height	410 mm
Overall Width	790 mm
Overall Depth	540 mm
SWL	140 kg

TRANSFER BENCHES



Commode Transfer Bench

The Commode Transfer Bench has all the features of the standard transfer bench with the extra benefit of padded seat with open front.

- Two piece padded seat allows water to drain away from seating area
- Padded seat for user comfort
- Front opening for easier personal cleaning
- Rack that allows a bowl to slide under for commode use

SPECIFICATIONS:

Product Code	1019
Seat Width	630 mm
Seat Depth	405 mm
Seat Height	460 - 560 mm
Bath Rim Clear	340 - 440 mm
Back Height	380 mm
Overall Width	730 mm
Overall Depth	600 mm
SWL	125 kg



Juvo Heavy Duty Transfer Bench

The Juvo Heavy Duty Transfer bench is both compact and sturdy providing safe and easy access for bathing or showering.

With a weight limit of 200kg this suits a vast range of clients whilst remaining lightweight for easy manoeuvrability. The Heavy Duty bench features height adjustable legs, non slip feet, the option of suction cups and reversible backrest to suit all bath types. With a textured non-slip seat with drainage holes and side rail the transfer bench provides stable seating for our clients.

- Leg Extensions available.
- Suitable for Bariatric Users

SPECIFICATIONS:

Product Code	1785
Seat Width	680 mm
Seat Depth	410 mm
Seat Height	470 - 570 mm
Bath Rim Clear	390 - 490 mm
Back Height	360 mm
Overall Width	780 mm
Overall Depth	550 mm
SWL	200 kg

BEDSIDE COMMODES



Bedside Commode

The standard bedside commode offers padded armrests for comfort and features adjustable height, a removable padded seat cover and a plastic toilet seat with bowl and lid.

- Stackable design
- 125kg SWL

SPECIFICATIONS:

Product Code	7952PA
DVA	B1060
Overall Width	575 mm
Overall Height	940 - 1090 mm
Overall Depth	600 mm
Seat Width	440 mm
Seat Height	500 - 650 mm
Seat Depth	440 mm
Back Height	475 mm
Toilet Seat	475 - 625 mm
SWL	125 kg



Deluxe Bedside Commode

This bedside commode offers a full range of features including a slide back seat with durable toilet seat and pan rack. This comes standard with bowl and lid. A skirt that surrounds the seat conceals the bowl so the chair can be likened to a regular bedroom sitting chair.

- Comfortable upholstered slide back seat
- Comfortable upholstered shallow curve backrest
- Supplied with Plastic toilet seat, bowl and lid
- Padded wide armrests
- Powdercoated steel frame

* Available in various colours in medical grade vinyl upholstery.

SPECIFICATIONS:

Product Code	7955
DVA	B1062WA
Overall Width	600 mm
Overall Height	870 - 1020 mm
Overall Depth	510 mm
Seat Width	460 mm
Seat Height	420 - 570 mm
Seat Depth	440 mm
Back Height	470 mm
SWL	120 kg

SIZES AVAILABLE:

Product Code	Width	SWL
7960	500 mm	135 kg
7961	550 mm	175 kg
7962	600 mm	195 kg
7963	650 mm	210 kg

Shower Trolleys & Tables

Staebel Shower Trolley	19
Staebel Junior Shower Trolley	19
Staebel Bariatric Shower Trolley	20
Staebel Bariatric Electric Shower Trolley	20
Staebel Shower Trolley - Electric with Backrest	21
Iris Profiling Shower Trolley	21
Electric Hi-Lo Change Table	22
Junior Electric Hi-Lo Change Table	22

SHOWER TROLLEYS & TABLES



Stabel Shower Trolley

- Shower trolley with one-section top
- Central braking system
- Easy steering system (for ease of maneuvering trolley's front left castor can be locked in direction of steering)
- Side rails – foldable, height: 185 mm with safety release mechanism (Material – powder coated stainless steel)
- Metal base cover
- Water-resistant shower mattress (height: 270 mm)
- Mattress integrated drain with 645 – 2000 mm long flexible hose
- Water-resistant head pillow – semi – cylinder

SPECIFICATIONS:

Product Code	1209ST
Shower trolley top Width	650 mm
Shower trolley top Length	1900 mm
Shower trolley top Height	540 - 970 mm
Height with Castors	600 - 1030 mm
Trolley top Inclination	0° - 5°
Length	2070 mm
Width	1070 mm
SWL	170 kg



Stabel Junior Shower Trolley

- Shower trolley with one-section trolley top
- Adjustable 5° trolley top inclination
- Hydraulic lowering/raising of trolley top via a console from 54 cm to 97 cm castors with diameter of 150 mm
- Central braking system
- Easy steering system (for ease of maneuvering trolley's front left castor can be locked in direction of steering)
- Side rails – foldable, height: 18,5 cm with safety release mechanism (Material – powder coated stainless steel)
- Metal base cover
- Water-resistant shower mattress (height: 27 cm)
- Mattress integrated drain from 645 – 2000 mm long flexible hose
- Water-resistant head pillow – semi-cylinder
- Trolley weight: 79 kg

SPECIFICATIONS:

Product Code	1214ST
Shower trolley top Width	600 mm
Shower trolley top Length	1500 mm
Trolley top Inclination	0° - 5°
Length	1670 mm
Width	760 mm
SWL	170 kg

SHOWER TROLLEYS & TABLES



Staebel Bariatric Shower Trolley

- Hydraulic lowering/raising of the stretcher top via a console lift
- Lowering/raising range: 540 mm to 970 mm
- Castors with a diameter of 150 mm
- Central braking easy steering system (for ease of maneuvering the trolley's front left castor can be locked in the direction of steering directional wheel)
- Side rails – foldable, height: 185mm with safety release mechanism (Material – powder coated stainless steel)
- Metal base cover
- Water resistant shower mattress (height: 270mm)
- Mattress integrated drain with a 645 up to 2000mm long flexible hose
- Water resistant head pillow – semi – cylinder

Staebel Bariatric Electric Shower Trolley

- Electric lowering/raising of the stretcher top via a console lift
- Lowering/raising range: 540 mm to 970 mm
- Castors with a diameter of 150 mm
- Central braking easy steering system (for ease of maneuvering the trolley's front left castor can be locked in the direction of steering directional wheel)
- Side rails – foldable, height: 185 mm with safety release mechanism (Material – powder coated stainless steel)
- Metal base cover
- Water resistant shower mattress (Height: 270 mm)
- Mattress integrated drain with a 645 up to 2000mm long flexible hose
- Water resistant head pillow semi cylinder
- Trolley weight: 110 kg

SPECIFICATIONS:

Product Code	1216ST
Shower trolley top Width	900 mm
Shower trolley top Length	2000 mm
Shower trolley top Height	540 - 970 mm
Height with Castors	600 - 1030 mm
Trolley top Inclination	0° - 5°
Length	2170 mm
Width	1070 mm
SWL	250 kg

SPECIFICATIONS:

Product Code	1217ST
Shower trolley top Width	900 mm
Shower trolley top Length	2000 mm
Shower trolley top Height	540 - 970 mm
Height with Castors	600 - 1030 mm
Trolley top Inclination	0° - 5°
Length	2170 mm
Width	1060 mm
SWL	250 kg

SHOWER TROLLEYS & TABLES



Staebel Shower Trolley - Electric with Backrest

- Shower trolley with two-section top
- Adjustable 5° trolley top inclination
- Electrically lowering/raising of trolley top via console from 54 cm to 97 cm
- Electrically adjustable head section inclination 45° (head section length: 70 cm)
- Castors with diameter of 150 mm
- Central braking system
- Easy steering system (for ease of maneuvering trolley's front left castor can be locked in direction of steering)
- Side rails – foldable, height: 18,5 cm with safety release mechanism (Material – powder coated stainless steel)
- Metal base cover
- Water-resistant shower mattress (height: 27 cm)
- Mattress integrated drain with 645 – 2000 mm long flexible hose
- Water-resistant head pillow – semi-cylinder
- Trolley weight: 110 kg

SPECIFICATIONS:

Product Code	1211ST
Shower trolley top Width	650 mm
Shower trolley top Length	1900 mm
Trolley top Inclination	0° - 5°
Length	2070 mm
Width	810 mm
SWL	170 kg



Iris Profiling Shower Trolley

Profiling shower trolley has been designed to facilitate shower bathing and hygiene routines, for patients and service users with a wide range of postural and cognitive challenges.

Its three section profiling mattress platform can be adjusted to provide specific postural support in the four key regions of the anatomy, head, back, seat, and legs providing maximum safety, support, and comfort.

Profiling shower trolley is also kind on the care provider thanks to its variable seat and platform height adjustment.

Care workers can maximise their comfort and safety using this function which enhances their well-being. Profiling shower trolley is also a mobile patient transfer system and can be used to transfer patients and service users in the horizontal or seated position.

- Multi-position posture management
- Variable height mattress platform
- Comfortable profiling mattress
- Waterproof mattress with high sides
- Lockable security side rails
- Comfort pillow and body belt
- Retracting stabilising chassis



SPECIFICATIONS:

Product Code	1310ST
Overall Height	580 - 980 mm
Overall Length	1970 mm
Overall Width	550 - 750 mm
SWL	160 kg

SHOWER TROLLEYS & TABLES



Electric Hi-Lo Change Table

Mobile, stable and compact in size, this change table is ideal to fit anywhere and suit any place. The mobile change table features electric hi-low and electric backrest.

- Electric hi-lo backrest raise/lower
- 125mm braking castors
- Designed for hoist access
- Handset operated
- 180kg capacity
- Quality Australian made product

SPECIFICATIONS:

Product Code	7850
Length	1880 mm
Width	7000 mm
Min Electric Hi - Low	620 mm
Max Electric Hi - Low	980 mm
SWL	180 kg



Junior Electric Hi-Lo Change Table

Ideal for children, this compact change table features electric hi-lo, padded top, foot control operation. Optional extras include vertical side rails, lap strap.

- Vinyl upholstered top 700mm wide x 1400mm long
- Electric hi-lo (630mm – 930mm)
- 125mm braking castors
- Designed for hoist access
- Foot control operation
- 90kg capacity
- Quality Australian made product

SPECIFICATIONS:

Product Code	7830
Length	630 mm
Width	930 mm
Min Electric Hi - Low	700 mm
Max Electric Hi - Low	1400 mm
SWL	90 kg

COMMODES



Trix 3 in 1 Shower Chair

- Assists patient with showering and toileting.
- Clear rear open space enabling positioning over the toilet.
- Easy access for hygiene with open armrest design.
- Deeper access front trough for easy cleaning.
- Integral pan carrier glides. Back saver design, assists with emptying the bowl, can be removed half way, self supporting while the bowl/pan is removed.
- Grade 304 stainless steel heavy wall thickness tubing forming robust frame to withstand extensive wet area use.
- Reinforced seat/back frame and armrests – weight bearing.
- Easy removable seat, using the Juvo Registered Design “Klipps” System
- Transportation wheels for effortless manoeuvring.
- Interchangeable feet so chair can be pulled by back or armrest
- Adjustable seat height for improved client positioning.
- Fully adjustable seat depth to assist with client positioning.
- Selection of seamless Polyurethane padded seats & pressure care seats (open or closed fronts).
- Armrests project beyond the seat to assist patient in sitting down and standing up.
- Splayed armrest design to eliminate pinch point to patient when sitting down and standing up.
- Flexible impervious non-woven back straps conforming to client shape for support and comfort. Antibacterial and non-suction designed straps to avoid skin tear when client is wet.
- Curved Backrest design to allow for natural prominence of the gluteal shelf.

SPECIFICATIONS:

Product Code	1783	1779	1780	1781	1782
Seat Width	460 mm	510 mm	610 mm	710 mm	910 mm
Width Between Arms	510 mm	560 mm	660 mm	760 mm	960 mm
Overall Depth	650 - 750 mm	650 - 750 mm	650 - 750 mm	650 - 750 mm	710 - 810 mm
Wall to Rear of Aperture	240 - 340 mm	240 - 340 mm	240 - 340 mm	240 - 340 mm	240 - 340 mm
Toilet Bowl Clearance - Width x Height	390 x 500 mm	440 x 500 mm	540 x 500 mm	640 x 500 mm	840 x 500 mm
To Underside of Pan Carrier - Width x Height	390 x 480 mm	440 x 480 mm	540 x 480 mm	640 x 480 mm	840 x 480 mm
Armrest Height from Seat	210 mm	210 mm	210 mm	210 mm	210 mm
Backrest Height	430 mm	430 mm	430 mm	430 mm	430 mm
Overall Height	920 - 1100 mm	920 - 1100 mm	920 - 1100 mm	920 - 1100 mm	920 - 1100 mm
Overall Width	630 mm	680 mm	780 mm	880 mm	1080 mm
Seat Depth - 5 settings at 25mm increments	470 - 570 mm	470 - 570 mm	470 - 570 mm	470 - 570 mm	495 - 595 mm
Seat Height - 7 settings at 25mm increments	410 - 560 mm	410 - 560 mm	410 - 560 mm	410 - 560 mm	410 - 560 mm
SWL	400 kg	400 kg	400 kg	400 kg	400 kg

COMMUNES



Broda Tilt & Recline Shower Commode

The Broda commode shower chair was designed for use with patients who require a great deal of assistance by a caregiver when being showered or toileted. The chair accommodates easy side transfers with its removable armrests. Front transfers are also facilitated by anterior tilt and a flip up removable footrest. The legrest also has the ability to be elevated. The Broda Commode has a padded commode seat, 12 degrees of posterior tilt, 10 degrees anterior tilt and up to 67 degrees of back recline - numerous positioning options to reduce slumping, sliding and falls while maintaining your residents comfort and safety.

- Stainless steel powdercoated frame for longevity
- Three total lock, one directional lock stainless steel castors
- Padded open front seat upholstered in durable vinyl
- Bed pan & rack (rack slides in and out for easy removal)
- 180 kg SWL

SPECIFICATIONS:

Product Code	10076	10077
Depth	460 mm	560 mm
Width	724 mm	724 mm
Seat Depth	460 mm	460 mm
Seat Height	600 mm	600 mm
Seat Width	460 mm	560 mm
SWL	180 kg	180 kg



Multi Tilt & Recline Mobile Commode


Designed to enhance the user outcome through safety, support and greater interaction with the carer, for a happier toileting and showering experience. The Multi Tilt & Recline mobile shower commode allows the user to be placed in the perfect position for showering thanks to its ability to tilt and recline up to 130 degrees, lift to a seat height of 100cm for best practice manual handling, and support up to 150kg weight.

- Eye to eye contact for carers manual handling and reduced back strain
- Reduced storage space compared to shower trolleys
- More space for transfers compared to shower trolley.
- Reduced risk of patient infection due to better access for perineal care
- Time saved in ease of use and general access
- Greater patient/ resident profile can be used with the Multi due to its unique tilting function
- Automatic leg support function when tilted
- Swing open armrest for sling access
- TGA listed

SPECIFICATIONS:

Product Code	19352
Length	1120 mm
Width	680 mm
Seat Length	430 mm
Seat Height	600 - 1000 mm
Seat Width	540 mm
Weight	53.6 kg
Rating	IP56
SWL	150 kg

COMMODES



Juvo Attendant Propelled Shower Commode with Sliding Footplate

The ultimate in strength, quality and stability this Juvo mobile shower chair is ideal for the Petite client! The wheelbase design prevents forward tipping and has proven to be an effective solution in all areas from home to aged care to hospital.

- Attendant propelled
- Quality stainless steel frame for longevity in a demanding environment
- Flip up armrests for ease of side transfers

- Armrests can be flipped so that they are positioned in front of occupant to reduce risk of slipping forward when being wheeled
- Smooth, slatted backrest design to facilitate easy cleaning
- Generous size sliding stainless steel footplate, weight bearing
- To utilise as a commode you must source the slide-in pan rack

***Please note that the seat for this product is sold separately**

Seat colours can be varied:



Lemon Twist

Lime

Lipstick

Regal

Emerald

DVA **MASS**

SPECIFICATIONS:

Product Code	1910J
DVA	B1020S
Seat Width	460 mm
Seat Depth	435 mm
Width Between Arms	490 mm
Overall Width	540 mm
SWL	200 kg

SIZES AVAILABLE:

1904J	1914J	1912J	1930J	1936J	1971J
	B1020SWA	B1020SWB	B1020SWC	B1020SWD	
410 mm	510 mm	560 mm	610 mm	660 mm	710 mm
435 mm	435 mm	435 mm	435 mm	435 mm	435 mm
440 mm	540 mm	590 mm	640 mm	690 mm	740 mm
490 mm	590 mm	640 mm	690 mm	740 mm	790 mm
200 kg	300 kg	300 kg	300 kg	300 kg	400 kg

COMMODES



Juvo Mobile Shower Commode with Swing Away Legrests

The ultimate in strength, quality and stability the Juvo mobile shower commode is built for heavy use having a 200kg SWL! The long wheelbase stops forward tipping and has proven to be an effective solution in all areas from home to hospital.

- Attendant propelled
- Assists transportation of the patient for showering and toileting
- Can be positioned over the toilet
- Durable total stainless steel frame
- Selection of seamless PU padded seats & pressure care seats (open or closed front)
- Removable seat

- Reinforced seat/back frame
- Swing away, removable and height adjustable legrests with flip up footplates
- Patented design, no protrusions
- Stainless steel 125mm castors (4 braking)
- Swing back independent armrests to assist in side transfer and patient access
- Two position height adjustable armrests
- Safety arms provide user security and safety
- Stainless steel pan rack guides
- This product is tested and certified to Australian Standards

***Please note that the seat for this product is sold separately**

Seat colours can be varied:



Lemon Twist

Lime

Lipstick

Regal

Emerald

SPECIFICATIONS:

Product Code	1915J
Seat Width	460 mm
Anterior Seat Height	590 mm
Posterior Seat Height	565 mm
Width Between Arms	465 mm
Overall Width	540 mm
SWL	200 kg

SIZES AVAILABLE:

1905J	1916J	1917J	1935J
410 mm	510 mm	560 mm	610 mm
595 mm	590 mm	590 mm	590 mm
505 mm	565 mm	565 mm	565 mm
440 mm	520 mm	570 mm	660 mm
490 mm	610 mm	640 mm	710 mm
200 kg	300 kg	300 kg	300 kg

COMMODES



Juvo Self-Propelled Shower Commode with Swing Away Legrests

- Assists transportation of the patient for showering and toileting
 - Can be positioned over the toilet
 - Durable total stainless steel frame
 - Selection of seamless PU padded seats & pressure care seats (open or closed fronts)
 - Removable seat
 - Reinforced seat/back frame
 - Swing away, removable & height adjustable legrests with flip-up footplates.
 - Patented TEARSAFE design, no protrusions
 - Now with braking stainless steel 125mm castors
 - Swing back independent armrests to assist in
 - Side transfer and patient access
 - Two position height adjustable armrests
 - Safety arms provide user security & safety
 - Stainless steel pan rack guides
 - This product is tested and certified to Australian Standards
- *Please note that the seat for this product is sold separately.**

Seat colours can be varied:



Lemon Twist

Lime

Lipstick

Regal

Emerald

SPECIFICATIONS:

Product Code	1925J
Seat Width	460 mm
Anterior Seat Height	590 mm
Posterior Seat Height	565 mm
Width Between Arms	490 mm
SWL	200 kg

SIZES AVAILABLE:

1926J	1927J	1945J	1939J
510 mm	560 mm	610 mm	660 mm
595 mm	595 mm	595 mm	595 mm
565 mm	565 mm	565 mm	565 mm
540 mm	590 mm	640 mm	690 mm
300 kg	300 kg	300 kg	300 kg

COMMODES



Juvo Self Propelled Shower Commode with Sliding Platform

Designed for the bariatric patient this model is the ultimate in strength, quality and stability. The long wheelbase stops forward tipping and has proven to be an effective solution in all areas from home to aged care to hospital.

- Assists transportation of the patient for showering and toileting
- Can be positioned over the toilet
- Durable total stainless steel frame
- Selection of seamless PU padded seats & pressure care seats (open or closed fronts)
- Removable seat
- Reinforced seat/back frame
- Sliding footrest

- Patented TEARSAFE design, no protrusions
- Now with braking stainless steel 125mm castors
- Swing back independent armrests to assist in side transfer and patient access
- Two position height adjustable armrests
- Safety arms provide user security & safety
- Stainless steel pan rack guides
- This product is tested and certified to Australian Standards

*** Please note that the seat for this product is sold separately.**

Seat colours can be varied:



Lemon Twist

Lime

Lipstick

Regal

Emerald

SPECIFICATIONS:

Product Code	1920J
Seat Width	460 mm
Anterior Seat Height	590 mm
Posterior Seat Height	565 mm
Width Between Arms	465 mm
SWL	200 kg

SIZES AVAILABLE:

1921J	1922J	1940J	1938J
510 mm	560 mm	610 mm	660 mm
590 mm	590 mm	590 mm	595 mm
565 mm	565 mm	565 mm	565 mm
530 mm	630 mm	630 mm	680 mm
300 kg	300 kg	300 kg	300 kg

COMMODES



Juvo Height Adjustable Attendant Propelled with Swing Away Legrests

The ultimate in strength, quality and stability the Juvo mobile shower commode is built for heavy use having a 200kg SWL! The long wheelbase stops forward tipping and has proven to be an effective solution in all areas from home to hospital.

- Attendant propelled
- Height adjustable up to 90mm front and rear axle.
- Assists transportation of the patient for showering and toileting
- Can be positioned over the toilet
- Durable total stainless steel frame
- Selection of removable seamless PU padded seats & pressure care seats (NB: Pictured with PU open front option)

- Reinforced seat/back frame
- Swing away, removable and height adjustable legrests with flip up footplates
- Patented design, no protrusions
- Stainless steel 125mm castors (4 braking)
- Swing back independent armrests to assist in side transfer and patient access
- Safety arms provide user security and safety
- Stainless steel pan rack guides
- This product is tested and certified to Australian Standards

***Please note that the seat for this product is sold separately**

SPECIFICATIONS:

Product Code	1906J
Seat Width	410 mm
Overall Width	490 mm
Anterior Seat Height	505 - 595 mm
Posterior Seat Height	475 - 565 mm
Width Between Arms	440 mm
SWL	200 kg

SIZES AVAILABLE:

1913J	1918J	1919J
460 mm	510 mm	560 mm
540 mm	590 mm	640 mm
505 - 595 mm	505 - 595 mm	505 - 595 mm
475 - 565 mm	475 - 565 mm	475 - 565 mm
490 mm	540 mm	590 mm
200 kg	300 kg	300 kg

COMMODOES



Juvo 'Tilting in Space' Shower Commode

- Assists transportation of the patient for showering and toileting
- Gas strut tilt mechanism for smooth and easy tilting
- Multi-position 35° tilt
- Can be positioned over the toilet
- Durable total stainless steel frame
- Selection of seamless PU padded seats & pressure care seats (open or closed front)
- Removable seat
- Reinforced seat/back frame
- Patented TEARSAFE design, no protrusions

- Swing away, removable and height adjustable legrests with flip up footplates
- Patented TEARSAFE design, no protrusions
- Stainless steel 125mm castors (4 braking)
- Swing back independent armrests to assist in side transfer and patient access
- Adjustable seat height
- Two position height adjustable armrests
- Safety arms provide user security and safety
- Stainless steel pan rack guides

***Please note that the seat for this product is sold separately**
Electric kit available for all sizes

Seat colours can be varied:



Lemon Twist Lime Lipstick Regal Emerald

SPECIFICATIONS:

Product Code	1951J
Seat Width	460 mm
Overall Width	560 mm
Armrest Height - High	240 mm
Armrest Height - Low	200 mm
Backrest Height from Seat	800 mm
Width Between Arms	490 mm
SWL	200 kg

SIZES AVAILABLE:

1950J	1952J	1956J	1961J
410 mm	510 mm	560 mm	610 mm
560 mm	610 mm	660 mm	730 mm
240 mm	240 mm	240 mm	240 mm
200 mm	200 mm	200 mm	200 mm
800 mm	800 mm	800 mm	800 mm
440 mm	540 mm	590 mm	610 mm
200 kg	300 kg	300 kg	300 kg

COMMODES



Juvo Folding Self-Propelled with Swing Away Legrests

- Assists transportation of the patient for showering and toileting
- Can be positioned over the toilet
- Durable total stainless steel frame
- Selection of seamless PU padded seats & pressure care seats (open or closed fronts)
- Removable seat
- Reinforced seat/back frame
- Swing away, removable & height adjustable legrests with flip-up footplates.
- Patented design, no protrusions
- Now with braking stainless steel 125mm castors
- Swing back independent armrests to assist in side transfer and patient access
- Two position height adjustable armrests
- Safety arms provide user security & safety
- Stainless steel pan rack guides
- Max User Weight 200kg
- Australian Standards Approved
- Patented folding design into standard suitcase size for easy travel

***Please note that the seat for this product is sold separately**

SPECIFICATIONS:

Product Code	1923J
Seat Width	460 mm
Overall Width	520 mm
Anterior Seat Height	500 - 590 mm
Posterior Seat Height	475 - 565 mm
Width Between Arms	465 mm
SWL	200 kg

SIZES AVAILABLE:

1924J
510 mm
540 mm
590 mm
565 mm
465 mm
200 kg



Juvo Height Adjustable Self Propelled with Swing Away Legrests

The ultimate in strength, quality and stability the Juvo mobile shower chair is built for heavy use having a 200kg SWL! The long wheelbase stops forward tipping and has proven to be an effective solution in all areas from home to aged care to hospital.

- Height adjustable up to 90mm front and rear axle
- Assists transportation of the patient for showering and toileting
- Can be positioned over the toilet
- Durable total stainless steel frame
- Selection of seamless PU padded seats & pressure care seats (open or closed fronts)
- Removable seat
- Reinforced seat/back frame
- Swing away, removable & height adjustable legrests with flip-up footplates.
- Patented TEARSAFE design, no protrusions
- Now with braking stainless steel 125mm castors
- Swing back independent armrests to assist in side transfer and patient access
- Two position height adjustable armrests
- Safety arms provide user security & safety
- Stainless steel pan rack guides
- This product is tested and certified to Australian Standards
- N.B. The legs elevated legrests shown do not come as standard.

***Please note that the seat for this product is sold separately**

SPECIFICATIONS:

Product Code	1928J
Seat Width	460 mm
Overall Width	660 mm
Anterior Seat Height	505 - 595 mm
Posterior Seat Height	475 - 565 mm
Width Between Arms	465 mm
SWL	200 kg

SIZES AVAILABLE:

1929J	1946J
510 mm	560 mm
710 mm	760 mm
505 - 595 mm	505 - 595 mm
475 - 565 mm	475 - 565 mm
540 mm	590 mm
200 kg	300 kg

Juvo Commode Accessories

Padded Armrest	34	Calf Support Strap	35
Arm Trough	34	Contoured Headrest	35
Transfer Bar / Safety Arm	34	Chest / Pelvic Strap	35
One Piece Footplate	34	Juvo Commode IV Pole	35
Elevated Legrest	34	Juvo Oxygen Bottle Carrier	36
Stump Support	34	Juvo Lock Down Armrest Kit	36
Lateral Support	35	Juvo Angle Adjusted Footrest	36
Drop Down Pommel	35	Juvo Angle Adjusted Backrest	36

JUVO COMMODE ACCESSORIES



Padded Armrest

- Padded urethane armrests including hardware to mount
- Sold singularly
- Fixings Included

SPECIFICATIONS:

Product Code 19611



Arm Trough

- The arm trough is of benefit for the shower chair user who has upper limb weakness or paralysis. It attaches to the chair armrest via a sturdy clamp.
- Can be mounted left or right side
 - Sold singularly
 - Fixings Included

SPECIFICATIONS:

Product Code 19613



Transfer Bar / Safety Arm

- Pair
- Handgrip
- Spring and button

SPECIFICATIONS:

Product Code 11950 - T



One Piece Footplate

- One piece footplate for swing-away shower commode models
- Different widths for wider commodes

SPECIFICATIONS:

Product Code 11963 - 46



Elevated Legrest

- Adjustable elevated legrest for Juvo mobile shower commodes. Installed in replacement for normal design legrest on a needs basis.
- Adjustable angle and length to suit the individual user
 - Left and Right side available (specify when ordering)

SPECIFICATIONS:

Product Code	Side
11940 - L	Left
11940 - R	Right



Stump Support

- Adjustable stump support for Juvo mobile shower commodes. Installed in replacement for normal design legrest on a needs basis.
- Adjustable angle and length to suit the individual user
 - Left and Right side available (specify when ordering)

SPECIFICATIONS:

Product Code	Side
11943 - L	Left
11943 - R	Right

JUVO COMMODE ACCESSORIES



Lateral Support

- Quick Release
- Swing-away left and right
- Includes Clamp and PU cover for bolts

SPECIFICATIONS:

Product Code	Side
11999 - L	Left
11999 - R	Right



Dropdown Pommel

- Height and depth adjustable

SPECIFICATIONS:

Product Code	
11985	



Calf Support Strap

Calf Support now has a colour matched, clip-on calf support for a complete, protective, hygienic solutions for adults and children using shower chairs.

- Reduce risks and increase comfort for shower chair users
- 41cm, 46cm and 51cm wide available

SPECIFICATIONS:

Product Code	
1BB216 - 22	



Contoured Headrest

Contoured adjustable headrest to suit Juvo. Comes complete with mounting accessories.

- Fixing Clamp
- Small option available



Chest / Pelvic Strap

Adjustable Chest and Pelvic Strap suitable for Juvo Commodes.

SPECIFICATIONS:

Product Code	
11933	



Juvo Commode IV Pole

Designed for use with shower commodes, this heavy duty, adjustable stainless steel IV pole simply clamps onto the shower commode frame.

SPECIFICATIONS:

Product Code	
11942	
Height (From Floor)	1200 - 1800 mm

SPECIFICATIONS:

Product Code	
19652	

JUVO COMMODE ACCESSORIES



Juvo Oxygen Bottle Carrier

This heavy duty stainless steel oxygen bottle carrier clamps onto the rear of the shower commode frame, and can hold a 130 mm diameter oxygen bottle.



Juvo Lock Down Armrest Kit

The lock down armrest option prevents the armrests from being flipped or pulled back and is useful in situations where this may pose a risk. The armrest is release by simply pulling the release button.

SPECIFICATIONS:

Product Code 11941



Juvo Angle Adjusted Footrest

- Juvo Commode Part - Angle adjustable footrest
- Left and Right side available (specify when ordering)

SPECIFICATIONS:

Product Code	Side
11931 - L	Left
11931 - R	Right

SPECIFICATIONS:

Product Code 19614



Juvo Angle Adjusted Backrest

- Shower Commode Accessory
- Angle Adjustable Backrest
- 3 Straps
- Multiple widths available
- Also available with 4 Straps

SPECIFICATIONS:

Product Code 11991

Shower Commode Seats

Hand Padded Open Front	38
Hand Padded Closed Front	38
Deluxe Hand Padded Open Front	38
Juvo PU Open Front	38
Juvo PU Closed Front	38
Customised Seats	38

SHOWER COMMODE SEATS



Hand Padded Open Front

- Completely sealed for infection control
- No sharp joins or seams
- Padded for increased pressure care
- Front opening

SPECIFICATIONS:

Product Code	Width
196490	410 mm
19648	460 mm
196630	510 mm
196550	560 mm
19668	610 mm



Hand Padded Closed Front

- Completely sealed for infection control
- No sharp joins or seams
- Padded for increased pressure care

SPECIFICATIONS:

Product Code	Width
196490	410 mm
19647	460 mm
19663C	510 mm
19655C	560 mm
19667	610 mm



Deluxe Hand Padded Open Front

- Completely sealed for infection control
- No sharp joins or seams
- Generously padded for high care
- Available with closed front

SPECIFICATIONS:

Product Code	Width
19649D	410 mm
19648D	460 mm
19663D	510 mm
19655D	560 mm
19671D	610 mm
19673D	660 mm



Juvo PU Open Front

- Completely sealed for infection control

SPECIFICATIONS:

Product Code	Width
196430	410 mm
196440	460 mm
196560	510 mm
196560	560 mm
196610	610 mm
196660	660 mm



Juvo PU Closed Front

- Completely sealed for infection control

SPECIFICATIONS:

Product Code	Width
19643C	410 mm
19644C	460 mm
19651C	510 mm
19655C	560 mm
19661C	610 mm
19666C	660 mm



Customised Seats

This is a unique pressure management commode seat cushion combination specially developed for those at high risk of skin breakdown and pressure ulcers.

- 3 foam layers, 79mm combination cushioning to minimise pressure.
- Vinyl fabric

* **Contact us about customised seats.**

SPECIFICATIONS:

Product Code	Width
19780	

ANCILLARY PRODUCTS



Oasis II Bath Board

The Oasis 2 Bath board is a moulded lightweight but strong plastic bath board, with easily adjustable rubber supports. It provides a comfortable and secure platform for showering over a bathtub. Supplied with moulded handle for ease of transfer. Anti-slip rubber supports adjust to fit even the narrowest of baths.

- Adjustable removable grab rail to assist with transfers
- Independently adjustable pillars to securely fit most baths
- Non slip surface with drain holes
- Integrated soap dish holder

SPECIFICATIONS:

Product Code	1020
Width	680 mm
Length	370 mm
Bath Width Min	390 mm
Math Width Max	636 mm
SWL	150 kg



Numo Toilet Safety Arms

With its stylish and clean construction, the Numo arms are the latest addition to the toilet support range. Simply and easily fitted, the Numo arms provide assistance with toilet transfers without the need for a cumbersome frame.

- Very sturdy construction
- Secures to pan with your existing toilet seat bolts
- Unobtrusive stylish design, suitable for public areas
- Australian designed and made



SPECIFICATIONS:

Product Code	1819 - 3
Arm Height Min	200 mm
Arm Height Max	300 mm
Width Between Arms	390 mm
SWL	150 kg



Throne Arms

The Throne is a 3 in 1 toilet support rail system, providing assistance to stand from the toilet. The rails can be positioned in the standard, splayed or fold down position, making this a versatile and adaptable piece of equipment. It is easily attached to the toilet utilising the existing toilet seat fixings.

- Durable powder coated steel construction (White in colour)
- Ergonomic arm design for maximum support
- Rubber lined clamps for stable fitting
- Throne spacer available in 50 or 80mm height if increased toilet seat height is necessary.

SPECIFICATIONS:

Product Code	1830
Arm Height	260 mm
Width Between Arms (Upright)	510 mm
Width Between Arms (Splayed)	640 mm
SWL	110 kg

ANCILLARY PRODUCTS



Bathing Cushion Bath Lift

The Bathing Cushion makes taking a bath easy as it gently lowers you into the bath and, later raises you with confidence up to the top of the bath so that you can simply step out. This new design offers increased stability and comfort and allows you to lie right down at the bottom of the bath.

- Requires no installation
- The option of leaving some air in the product whilst bathing to act as a cushion that moulds around you for additional comfort and support.
- Can be packed away into a small hold for ease of transportation or to store it away when not required, the perfect solution for people who love travelling, weighs only 2kg
- The lumbar support enables you to choose whether to sit up or lie back

SPECIFICATIONS:

Product Code	1649
Width	580 mm
Min Height	20 mm
Max Height	400 mm
Depth	440 mm
Back/Head Length	330 mm
Weight	2 kg
SWL	150 kg



Archimedes Bath Lift

The Archimedes Bath Lift will lower you right down to the bottom of the bath steadily and with confidence, maximising your leg room for an extremely enjoyable bath. The product separates easily into three parts with the heaviest weighing only 5.2 kg (11½ lb.), minimising risk when it is lifted in and out of the bath.

- Great alternative to a walk-in bath
- Simple to fit or remove – separates into three pieces
- Battery operated – no dangerous power cables
- Fully waterproof hand control
- Safety feature: will only lower you into the bath when there is sufficient battery power to raise you back out again
- Lifts up to 150kg (24 stone)
- 2 year guarantee

SPECIFICATIONS:

Product Code	1600
Width	395 mm
Min Height	70 mm
Max Height	440 mm
Depth	445 mm
Back/Head Length	665 mm
Weight	5.2 kg
SWL	150 kg

ANCILLARY PRODUCTS



Standard Bed Pan

- This bedpan is made from polypropylene.
- Bed pan is autoclavable.
- The standard bedpan is heavy duty and ideal for institutional use.

SPECIFICATIONS:

Product Code 1990



Slipper Bed Pan

- This bedpan is made from polypropylene.
- Pan is autoclavable.
- The slipper pan is ideal for use with orthopaedic and frail people.

SPECIFICATIONS:

Product Code 1993



Male Urinal

A clear polypropylene bottle with integrated handle. Calibrated for output recording, bottle is autoclavable.

- Capacity of 1 Litre,
- Graduated in 100ml increments

SPECIFICATIONS:

Product Code 1980



Female Urinal

Clear polypropylene bottle with integral handle, anatomically shaped at opening. Calibrated for output recording and autoclavable.

- Capacity of 1 litre - 100 mL graduations.

SPECIFICATIONS:

Product Code 1983



Bowl & Lid

- Bowl with Lid and Handle
- Can be used with all Shower Commode types

SPECIFICATIONS:

Product Code 11936



Pan/Bowl Carrier

- Pan and Bowl Carrier that can be used for different commodes.
- Different sizing available.

SPECIFICATIONS:

Product Code 19607



Spil-Pruf Urinal

These urinals will not spill or leak in any position. Spil-Pruf urinals help reduces odour problems while increasing confidence.

- Contents are not visible
- Autoclavable

SPECIFICATIONS:

Product Code 1985



Female Adaptor for Spil-Pruf

These urinals will not spill or leak in any position. Spil-Pruf urinals help reduces odour problems while increasing confidence.

- Contents are not visible
- Autoclavable
- Suitable for women

SPECIFICATIONS:

Product Code 1986



Bowl - No Lid

Suitable for the Active Care overtoilet frame and bedside commodes, this oval bowl measures 305mm diameter.

SPECIFICATIONS:

Product Code 1823

ANCILLARY PRODUCTS



Raised Toilet Seat - With Lid

- Lid included for attractive appearance
- Increases seat height of existing toilet simply and easily
- Secured with two screw in brackets with rubber grips
- Lip at front for increased stability
- 136kg weight capacity

SPECIFICATIONS:

Product Code 1726



Gel Commode Cushion - Closed

The polymer gel pad is made from Akton polymer to provide pressure management on commode and toilet seats.

- Easy to clean
- Hygienic
- Includes velcro attachments to prevent sliding

SPECIFICATIONS:

Product Code 1971



Gel Commode Cushion - Open

The polymer gel pad is made from Akton polymer to provide pressure management on commode and toilet seats.

- Easy to clean
- Hygienic
- Includes velcro attachments to prevent sliding

SPECIFICATIONS:

Product Code 1973



Infl for Closed Front Seat

- Fills in the gap in Shower Commode
- For Closed front seat
- Fits Juvo PU seat

SPECIFICATIONS:

Product Code 11998



Infl for Open Front Seat

- Fills in the gap in Shower Commode
- For open front seats
- Fits Juvo PU seat

SPECIFICATIONS:

Product Code 11997



Raised Toilet Seat - 50mm

- Increases seat height of existing toilet simply and easily
- Secured with two screw in brackets with rubber grips
- Lip at front for increased stability
- 136kg weight capacity

SPECIFICATIONS:

Product Code 1737



SEE SOMETHING YOU LIKE? **LET'S HAVE A CHAT!**

Book your consultancy today with one of our friendly team members. See below information for ways you can get in contact....



1800 267 267



Showroom Opening Hours

Monday to Friday | 8:00am to 4:00pm



www.activemedicalsupplies.com.au



Cannon Hill

35 Southgate Ave, Cannon Hill QLD 4170



hello@activemedicalsupplies.com.au



Toowoomba

11 Grange St Toowoomba QLD 4350

activeMEDICAL

S U P P L I E S > S O L U T I O N S > S U P P O R T

P. 1800 267 267 | **F.** 07 3899 8825

E. hello@activemedicalsupplies.com.au

www.activemedicalsupplies.com.au