

# NUM's Quick Guide For Essentials

Issue 1



# Contents

Capability Statement	4
Bath, Toilet & Shower	6
Patient Handling	11
Pressure Care	21
Wheelchairs	29
Seating	33
Mobility Solutions	38
Bedroom Solutions	41
Patient Safety	48
Treatment Tables	51

## CAPABILITY STATEMENT

---

**Active Medical Supplies is a pioneer in medical device and equipment manufacturing.** The experienced team at Active Medical understand both professionally and personally, the importance of maintaining independence within the comfort of your own home or the community including educational settings. A lot of the people we help are still independent, however, just need assistive solutions as result of an injury, Hospital visit, Healthcare professional's advice or to help loved ones with their daily life.

With over 30 years of experience helping with solutions in community care, we understand sometimes that daily tasks may become difficult. Often things that we have taken for granted like getting out of bed, lifting yourself from a chair or just even moving around the house, inherit more effort than they have before.

For those who come into our showroom, we can replicate the scenarios you would face at home with our daily living set ups. Our experienced and courteous staff gladly provide all the information and assistance you require.

We work extensively with Physiotherapists, Occupational Therapists and Aged Care Clinicians so that we can find the right form of solution for you. Over the course of 30 years, we have provided solutions to community care with diligence and professionalism. Our customers' assurance comes from the fact that we help you find the most suitable solution.


We take great pride in the assembling and manufacturing of all of our products. We culminate to the highest standards in care, comfort for the patient and project cost efficiency.

# Procurement & Distribution

Active Medical has developed a network of manufacturers in Australia and around the world – from Germany and Sweden to the United States of America and consults directly with industry professionals to ensure the products we provide responds to unique patient specifications.

The success of our business is supported by our ability to collaborate with manufacturers that not only provide the highest quality of products and services, but are aligned with our mission and commitment to our social, corporate and environmental responsibilities.

Our efficient distribution systems are underpinned by our strong relationships with manufacturers and two centralised warehouses with a capacity to hold up to 50,000 items. This ensures the fast turn around and on time delivery demanded by the health care industry and the needs of people in care. We ensure that all our products are ethically sourced and supplied and meet environmental standards.



### \*Disclaimer

Active Medical Supplies does not take responsibility for incorrect application of products. This publication is a guide only to assist in selecting products suitable to your clients. Please use in conjunction with professional advice.



## CAPABILITY STATEMENT

# Did you know...



We sanitize **all products** that leave our warehouses



We offer full product **demonstrations and video trials**



We offer **contact free delivery**



Our staff are **double vaxxed** with boosters




Our showroom is a **controlled environment** that we can also offer trials

# Come in & have a chat!

## Book a consultancy today!

Active Medical facilitates a showroom on site to replicate the home care environment. Our showrooms are used to replicate a residential care bedroom setting and allows us to offer unique demonstrations and training tailored to your needs.

We are located in Cannon Hill and Toowoomba with wheelchair accessibility. We showcase our extensive range of electric beds, wheelchairs, mobility, patient handling, safety equipment, bathroom solutions and accessories.



### Cannon Hill

35 Southgate Ave, Cannon Hill QLD 4170



### Toowoomba

11 Grange St Toowoomba City QLD 4350



### Showroom Opening Hours

Monday to Friday | 8:00am to 4:00pm



### Showroom Opening Hours

Monday to Friday | 8:00am to 4:00pm

\*Products and information listed in this catalogue are subject to change. Website and digital information supercedes all print. For the most up to date information, please contact us directly via the details listed above.



## BATH, TOILET & SHOWER



### Aluminium Shower Chair

This Shower Chair provides a comfortable seat for those unable to stand whilst showering. The armrests provide security when entering or exiting the shower, and assist the user on and off the seat.

- 3 piece padded vinyl seat
- Hand rail for added support while transferring
- Reversible padded backrest to suit every situation
- Welded frame

### SPECIFICATIONS:

<b>Product Code</b>	1315J
<b>DVA</b>	B4002 - A
<b>Seat Width</b>	440 mm
<b>Seat Depth</b>	430 mm
<b>Seat Height</b>	400 - 550 mm
<b>Bath Rim Clear</b>	460 mm
<b>Seat-Back Height</b>	350 mm
<b>Overall Width</b>	570 mm
<b>Overall Depth</b>	540 mm
<b>SWL</b>	200 kg



### Aluminium Over Toilet Frame

This aluminium over toilet frame is adjustable in height and fits over the toilet. This frame has been designed to assist the user who has difficulty sitting down or standing up by raising the toilet seat height.

- Aluminium construction for long life
- No risk of rust
- Removable splash guard included
- Lift up seat that removes for cleaning
- Ideal for those with hip or knee injuries or who require an increased seat height
- Unique stacking design saves space when storing more than one frame
- Multi D leg adjustment ensures stability on uneven surfaces
- Push button height adjustment
- Non slip guides on base of legs

### SPECIFICATIONS:

<b>Product Code</b>	1755J
<b>Width Between Arms</b>	450 mm
<b>Seat Height</b>	475 - 625 mm
<b>Overall Width</b>	560 mm
<b>Armrest Height</b>	700 - 850 mm
<b>Footprint (W x D)</b>	560 x 590 mm
<b>SWL</b>	150 kg

## BATH, TOILET & SHOWER



### Juvo PU Open Front

- Completely sealed for infection control
- Ideal for patients that won't spend a long time in the shower chair
- Open front for better access
- Easy removable for cleaning and replacement



### Standard Bed Pan

- This bedpan is made from polypropylene.
- Also autoclavable.
- The standard bedpan is heavy duty and ideal for institutional use.



### Pan/Bowl Carrier

- Pan and Bowl Carrier that can be used for different commodes. Different sizing available.

### SPECIFICATIONS:

Product Code	Size
196430	410 mm
196440	460 mm
196560	510 mm
196560	560 mm
196610	610 mm
196660	660 mm

### SPECIFICATIONS:

Product Code	Size
1990	

### SPECIFICATIONS:

Product Code	Size
19607	

## BATH, TOILET & SHOWER



### Trix 3 in 1 Shower Chair

- Assists patient with showering and toileting.
- Clear rear open space enabling positioning over the toilet.
- Easy access for hygiene with open armrest design.
- Deeper access front trough for easy cleaning.
- Integral pan carrier glides. Back saver design, assists with emptying the bowl, can be removed half way, self supporting while the bowl/pan is removed.
- Grade 304 stainless steel heavy wall thickness tubing forming robust frame to withstand extensive wet area use.
- Reinforced seat/back frame and armrests – weight bearing.
- Easy removable seat, using the Juvo Registered Design “Klipps” System
- Transportation wheels for effortless manoeuvring.
- Interchangeable so wheelchair can be pulled by back or armrest.
- Adjustable seat height for improved client positioning.
- Fully adjustable seat depth to assist with client positioning.
- Selection of seamless Polyurethane padded seats & pressure care seats (open or closed fronts).
- Armrests project beyond the seat to assist patient in sitting down and standing up.
- Splayed armrest design to eliminate pinch point to patient when sitting down and standing up.
- Flexible impervious non-woven back straps conforming to client shape for support and comfort. Antibacterial and non-suction designed straps to avoid skin tear when client is wet.
- Curved Backrest design to allow for natural prominence of the gluteal shelf.

### SPECIFICATIONS:

Product Code	1783	1779	1780	1781	1782
Seat Width	460 mm	510 mm	610 mm	710 mm	910 mm
Width Between Arms	510 mm	560 mm	660 mm	760 mm	960 mm
Overall Depth	650 - 750 mm	650 - 750 mm	650 - 750 mm	650 - 750 mm	710 - 810 mm
Wall to Rear of Aperture	240 - 340 mm	240 - 340 mm	240 - 340 mm	240 - 340 mm	240 - 340 mm
Toilet Bowl Clearance - Width x Height	390 x 500 mm	440 x 500 mm	540 x 500 mm	640 x 500 mm	840 x 500 mm
To Underside of Pan Carrier - Width x Height	390 x 480 mm	440 x 480 mm	540 x 480 mm	640 x 480 mm	840 x 480 mm
Armrest Height from Seat	210 mm	210 mm	210 mm	210 mm	210 mm
Backrest Height	430 mm	430 mm	430 mm	430 mm	430 mm
Overall Height	920 - 1100 mm	920 - 1100 mm	920 - 1100 mm	920 - 1100 mm	920 - 1100 mm
Overall Width	630 mm	680 mm	780 mm	880 mm	1080 mm
Seat Depth - 5 settings at 25mm increments	470 - 570 mm	470 - 570 mm	470 - 570 mm	470 - 570 mm	495 - 595 mm
Seat Height - 7 settings at 25mm increments	410 - 560 mm	410 - 560 mm	410 - 560 mm	410 - 560 mm	410 - 560 mm
SWL	400 kg	400 kg	400 kg	400 kg	400 kg

## BATH, TOILET & SHOWER



### Juvo Attendant Propelled Shower Commode with Sliding Footplate

The ultimate in strength, quality and stability this Juvo mobile shower chair is ideal for the Petite client! The wheelbase design prevents forward tipping and has proven to be an effective solution in all areas from home to aged care to hospital.

- Attendant propelled
- Quality stainless steel frame for longevity in a demanding environment
- Flip up armrests for ease of side transfers

- Armrests can be flipped so that they are positioned in front of occupant to reduce risk of slipping forward when being wheeled
- Smooth, slatted backrest design to facilitate easy cleaning
- Generous size sliding stainless steel footplate, weight bearing
- To utilise as a commode you must source the slide-in pan rack

**\*Please note that the seat for this product is sold separately**

Seat colours can be varied:



Lemon Twist

Lime

Lipstick

Regal

Emerald

DVA MASS

#### SPECIFICATIONS:





Product Code	1904J
DVA	
Seat Width	410 mm
Seat Depth	435 mm
Width between arms	440 mm
Overall Width	490 mm
SWL	200 kg

#### SIZES AVAILABLE:

1910J	1914J	1912J	1930J	1936J	1971J
B1020S	B1020SWA	B1020SWB	B1020SWC	B1020SWD	
460 mm	510 mm	560 mm	610 mm	660 mm	710 mm
435 mm	435 mm	435 mm	435 mm	435 mm	435 mm
490 mm	540 mm	590 mm	640 mm	690 mm	740 mm
540 mm	590 mm	640 mm	690 mm	740 mm	790 mm
200 kg	300 kg	300 kg	300 kg	300 kg	400 kg



## PATIENT HANDLING



### Alto 200 Patient Hoist

Alto Lift 200: A medium sized heavy duty mobile patient lifting hoist with a safe working load of 200kg and a compact ergonomic design. The Alto Lift 200, whilst being compact, provides superior patient lift space through the use of SPACE LINK™ technology.

The Alto Lift 200 is the lifter of choice for institutional use or home care where a high capacity lift is required. Designed and tested to comply with: AS ISO 10535-2002 Hoists for transfer of disabled persons - Requirements and test methods.

### SPECIFICATIONS:

Product Code	2891
Base Length	1150 mm
External Base Width	645 mm
Maximum Internal Length	860 mm
Internal Base Width Closed	578 mm
Internal Base Width Open	1130 mm
Base Width Adjustment	Electric
Under Bed Clearance	125 mm
Minimum Height	420 mm
Spreader Bar	300 mm
Maximum Height	870 mm
Spreader Bar	1920 mm
Hoisting Range*	1620 mm

Minimum Storage Height	1350 mm
Under Base Height	50 mm
Maximum Hoist Reach	1060 mm
Minimum Hoist Reach	610 mm
Max Hoist Reach at	41 kg
Max Height	4400 mm
Min Hoist Reach at Min Height	430 mm
Handle Height	1000 mm
Turning Radius	1210 mm
Mass of Heaviest Part (Base Assembly)	17.5 kg
Lifter Weight	42 kg
SWL	200 kg



## PATIENT HANDLING



### Salsa 200 Standing Lifter

Designed as an active standing aid allowing a single carer to raise a patient into a comfortable secure standing position. The Salsa 200 is ideal for transferring from bed to chair, or chair to commode. Allegro Concepts standing slings are manufactured to suit the Salsa, allowing superior access for dressing and personal hygiene.

With a small footprint and low overall weight, the Salsa 200 is highly maneuverable. A high capacity safe working load of 200kg means the Salsa is equally at home in an institutional or home environment.

Designed and tested to comply with:  
AS ISO 10535-2002 Hoists for the transfer of disabled persons – Requirements and test methods.

WANT TO LEARN MORE ABOUT OUR  
**SALSA 200 STANDING LIFTER?**

WATCH OUR VIDEOS ONLINE TODAY!

FIND US ON  
**YOUTUBE**   
**activemedicalsupplies**





### SPECIFICATIONS:

<b>Product Code</b>	27700
<b>Base Length</b>	1010 mm
<b>External Base Width Closed</b>	665 mm
<b>External Base Width Open</b>	1030 mm
<b>Internal Base Width Closed</b>	565 mm
<b>Internal Base Width Open</b>	930 mm
<b>Base Width Adjustment</b>	Electric
<b>Under Bed Clearance</b>	130 mm

<b>Minimum Height Lift Arm</b>	880 mm
<b>Maximum Height Lift Arm</b>	1685 mm
<b>Under Base Height</b>	55 mm
<b>Lifter Weight</b>	51 kg
<b>Minimum Storage Height</b>	1190 mm
<b>Footplate</b>	Removable
<b>Kneepad</b>	Height Adjustable
<b>SWL</b>	200 kg

## PATIENT HANDLING



### Pivot Sling

The Pivot sling is designed for use with the pivot/tilt frame. It attaches with the keyhole tags engineered to the highest standards in strength and size tolerances and feature auto clip on action (no safety straps required). It

also features sling centre marking and positioning handles to assist correct positioning and woven label that will not wash out or fade.

**\*Available in Allegro Mesh and Fabric**

### SPECIFICATIONS:


#### Allegro Mesh

Product Code	27P02XS	27P02S	27P02M	27P02L	27P02XL	27P02XXL
Size	Extra Small	Small	Medium	Large	Extra Large	Extra Extra Large
A	720 mm	870 mm	1000 mm	1080 mm	1170 mm	1320 mm
B	190 mm	230 mm	270 mm	260 mm	270 mm	290 mm
C	580 mm	830 mm	920 mm	1010 mm	1090 mm	1180 mm
D	870 mm	1220 mm	1310 mm	1410 mm	1540 mm	1710 mm

#### Allegro Fabric

Product Code	27P01XS	27P01S	27P01M	27P01L	27P01XL	27P01XXL
Size	Extra Small	Small	Medium	Large	Extra Large	Extra Extra Large
A	720 mm	870 mm	1000 mm	1080 mm	1170 mm	1320 mm
B	190 mm	230 mm	270 mm	260 mm	270 mm	290 mm
C	580 mm	830 mm	920 mm	1010 mm	1090 mm	1180 mm
D	870 mm	1220 mm	1310 mm	1410 mm	1540 mm	1710 mm

## PATIENT HANDLING



### General Purpose Head Support Sling

The General Purpose Head Support sling, as the name suggests, is a full body support style sling and is recommended for general lifting and in situations of lifting from the floor, and for the less secure patient. The loop style attachments have three lengths at each connection point and are used on

a patient lifter with a yoke style attachment. Available in sizes from x-small through to XX-large, providing a popular choice for individuals and care facilities. Slings are manufactured from polyester fabric or mesh (suitable for bathing).

### SPECIFICATIONS:


#### Allegro Mesh

Product Code	27GPHS02XS	27GPHS02S	27GPHS02M	27GPHS02L	27GPHS02XL	27GPHS02XXL
Size	Extra Small	Small	Medium	Large	Extra Large	Extra Extra Large
A	820 mm	790 mm	920 mm	1110 mm	1270 mm	1360 mm
B	210 mm	230 mm	240 mm	260 mm	270 mm	290 mm
C	750 mm	780 mm	920 mm	1000 mm	1020 mm	1050 mm
D	1310 mm	1280 mm	1410 mm	1630 mm	1680 mm	1780 mm

#### Allegro Fabric

Product Code	27GPHS01XS	27GPHS01S	27GPHS01M	27GPHS01L	27GPHS01XL	27GPHS01XXL
Size	Extra Small	Small	Medium	Large	Extra Large	Extra Extra Large
A	820 mm	790 mm	920 mm	1110 mm	1270 mm	1360 mm
B	210 mm	230 mm	240 mm	260 mm	270 mm	290 mm
C	750 mm	780 mm	920 mm	1000 mm	1020 mm	1050 mm
D	1310 mm	1280 mm	1410 mm	1630 mm	1680 mm	1780 mm


## PATIENT HANDLING



### Salsa Standing Sling

The Salsa sling is for use with the Salsa stand up hoist. It features optimum back support and under arm padding for a safe and comfortable lift and also allows for the easy removal of clothing for toilet transfers.

The Allegro Salsa sling features a torso band for security and stability with Buckle (b) or Velcro (v) fastening.



### SPECIFICATIONS:

#### Allegro Mesh

Product Code	27702	27703	27704
Size	Medium	Large	Extra Large
A	880 mm	1060 mm	1220 mm
B	290 mm	370 mm	410 mm
C (b)	960 mm	1020 mm	1220 mm
C (v)	1180 mm	1320 mm	1440 mm

#### Velcro Fastening

Product Code	27705	27706	27707
Size	Medium	Large	Extra Large
A	880 mm	1060 mm	1220 mm
B	290 mm	370 mm	410 mm
C (b)	960 mm	1020 mm	1220 mm
C (v)	1180 mm	1320 mm	1440 mm



### Staebel Patient Transfer

The Staebel Patient Transfer Device is the perfect solution for safe transfer of a patient or resident within the home or facility. This is your all-in-one solution without the need for additional lifting straps and belts, and promoting independence. The width adjustable base allows good access to most seating, whilst also designed for good access beneath beds etc.



### SPECIFICATIONS:

Product Code	212145
Seat Width	450 mm
Overall Width	640 mm
Overall Length	930 mm
SWL	150 kg

## PATIENT HANDLING



### Stabel Transfer Platform

Stabel Transfer Platform is a manually-operated device designed as a platform for use by healthcare professionals to assist in the physical moving (transfer) of a patient in an upright (standing) position over a short distance. The device consists of a wheeled platform, a surrounding framework, kneepad and handlebar. It can include safety belt, to hold the person being transferred in the standing position.



#### SPECIFICATIONS:

<b>Product Code</b>	212150
<b>Width</b>	620 mm
<b>Height</b>	1180 mm
<b>SWL</b>	150 kg



### Return Transfer Platform

The ReTurn Transfer Platform, an assistance device for transferring clients, is specially designed to make different sit-to-stand and transfer situations easier for both clients and care givers. It is very easy to manoeuvre. It has been developed to facilitate work in confined spaces and it requires less effort than conventional lifting. In addition, instead of remaining passive, the patient is activated, which is a great advantage in comparison with lifting. Since the patient's own capacity is utilised and the device is easy to manoeuvre, the ReTurn is ideal for both home care and institutional settings.

The design of the ReTurn allows the client to get close to and to rise to standing position using natural movements. All of this adds up to a high degree of utilisation. It is specially designed for short transfers in confined spaces. The knee supports are height adjustable which makes it easy to adapt the device for any patient. Many patients can use the Return from the same seat height. The ReTurn's base plate is low and stable. With a increase in height of only 4 cm, the patient does not need to lift his or her feet very high. Foot markers on the base plate simplify and ensure safe placement. ReTurn has two sections. This, together with its low weight, makes it ideal in a transport situation and the ReTurn is easy to keep clean and to store.

#### SPECIFICATIONS:

<b>Product Code</b>	217500I
<b>Size</b>	Adult
<b>Base Width</b>	685 mm
<b>Base Length</b>	570 mm
<b>Height</b>	1150 mm
<b>SWL</b>	150 kg

#### SIZES AVAILABLE:

<b>Product Code</b>	217400
<b>Size</b>	Small
<b>Base Width</b>	685 mm
<b>Base Length</b>	570 mm
<b>Height</b>	995 mm
<b>SWL</b>	150 kg



## PATIENT HANDLING



### Return Belt

ReTurnBelt is an excellent complement to ReTurn when more support and stability is needed, for instance when the user's weight-bearing ability is questionable or varied. ReTurnBelt can also be used stand alone for sit-to-stand, transfer and gait training procedures.

- Anti-slip inner surface
- PU-coated Nylon



### Walk Belt

The Walk Belt is designed to assist with the safe transfer and ambulation of the client. The construction of the Walk Belt is of nylon and cotton with the inside lining a slip resistant material to minimise the risk of the belt riding up during use.

The belts are 25cm wide and fasten with Velcro. Available in Small, Medium, Large and X-Large.



### SPECIFICATIONS:

<b>Product Code</b>	216034M
<b>Size</b>	Medium
<b>Width</b>	1050 mm
<b>Length</b>	1250 mm
<b>SWL</b>	205 kg

### SIZES AVAILABLE:

Product Code	216033S	216035L	216036XL	216037XXL
<b>Size</b>	Small	Large	Extra Large	Extra Extra Large
<b>Width</b>	850 mm	1200 mm	1350 mm	1550 mm
<b>Length</b>	1050 mm	1450 mm	1550 mm	1750 mm

### SPECIFICATIONS:

<b>Product Code</b>	2501S
<b>Width</b>	250 mm
<b>Size</b>	Small
<b>Length</b>	690 - 930 mm

### SIZES AVAILABLE:

Product Code	2501M	2501L	2501XL
<b>Size</b>	Medium	Large	Extra Large
<b>Width</b>	250 mm	250 mm	250 mm
<b>Length</b>	850 - 1100 mm	900 - 1170 mm	1080 - 1330 mm

## PATIENT HANDLING



### Staebel Chair Scale

The Staebel Chair Scale with its patented design maintains reliability and uncompromising precision. All of the user's care needs, including safety, comfort, care-giver's use and maintenance, are met. The compact design ensures the most cost efficient shipping method and reduces the risk of carer strain particularly when weighing an obese/ bariatric patient.

Unique base frame design gives care staff assurance that scale is compatible with other equipment utilised in care setting

- Twin wheel castors are standard and do not create excessive noise when wheeling scale along a corridor
- One piece footplate flips up and has been designed with skin tear prone patients in mind
- Armrests flip back to allow scale to be utilise in bariatric care setting
- Base frame comes standard with a level indicator
- BMI function/ Tare/Zero/Hold/Print/kg/lb
- 5 digit ¾" LCD screen that can be tilted to accommodate different care staff heights
- Automatic re-weigh
- Up-datable firmware
- Ports: RS232 and USB
- Amazing lightweight yet durable construction - only weighs 17.5 kg
- 24 month warranty

#### SPECIFICATIONS:

<b>Product Code</b>	2957
<b>Seat Width</b>	430 mm
<b>Seat Depth</b>	420 mm
<b>Overall Width</b>	630 mm
<b>Overall Height</b>	965 mm
<b>SWL</b>	300 kg



### Premium Slidesheet

The ideal solution when transferring or repositioning patients across, up or down the bed. Made from a strong, non flammable silicone coated polyester taffeta these premium quality sheets retain their slip for over 150 machine washes.

Large identification label sewn into edge to allow care facilities to allocate by unit and identify age.

- Can be machine washed to 75°C
- Slip consistently in all directions/both sides
- TGA Registered product

#### SPECIFICATIONS:

<b>Product Code</b>	2433
<b>Width</b>	2000 mm
<b>Length</b>	1450 mm

#### SIZES AVAILABLE:

<b>Product Code</b>	2441R
<b>Width</b>	2000 mm
<b>Length</b>	1500 mm



### EasyGlide Transfer Board

The flexible three-part sliding board is easy to place in front of the wheels of a wheelchair to protect against soiling during a sliding transfer. The smooth top surface makes it easy to move and provides for easy transfers while the double anti-slip tape on the underside prevents the board from sliding during transfer. It is especially suitable for levelling out height differences or for bridging short distances.

- Anti-slip protection on the underside provides for safety and security.
- The long sides can be folded to protect clothing from the wheelchair's wheels.

#### SPECIFICATIONS:

<b>Product Code</b>	215020
<b>Width</b>	330 mm
<b>Length</b>	600 mm
<b>Thickness</b>	5 mm
<b>SWL</b>	135 kg

#### SIZES AVAILABLE:

<b>Product Code</b>	215050
<b>Width</b>	330 mm
<b>Length</b>	750 mm
<b>Thickness</b>	6 mm

## PATIENT HANDLING



### Standard Yoke

The yoke is the standard attachment for general purpose slings. A unique sling retainer has been designed into the hook of the yoke to prevent the sling loop from coming detached during lifting.

#### SPECIFICATIONS:

Product Code 2882



### Pivot Attachment with Integrated Weigh Scales

Complete as one attachment, this allows for standard lifting without losing any critical lift height with separate attachment. Weighing scale Safe Weight Limit 300kg.

- Suitable for allegro Sonata & Alto Hoists.

#### SPECIFICATIONS:

Product Code 2895



### Four Point Yoke

The Allegro four point yoke can also be used with the general purpose slings and is particularly useful in the care and comfort of larger patients. A unique sling retainer has been designed into the hook of the yoke to prevent the sling loop from coming detached during lifting.

#### SPECIFICATIONS:

Product Code 2886



### Pivot Frame

The pivot frame is used with the pivot sling. The frame spreads the sling attachment points to ensure the patients hips are not flexed excessively, avoiding discomfort some patients experience with standard slings. The ability to position a person in a sitting or lying position, assisting in providing correct seating and positioning.

#### SPECIFICATIONS:

Product Code 2884



### Standard Yoke with Weigh Scales

Complete as one attachment, this allows for standard lifting without losing any critical lift height with separate attachment. Weighing scale Safe Weight Limit 300kg. Suitable for allegro Sonata & Alto Hoists.

#### SPECIFICATIONS:

Product Code 2894



### Lifting Frame with Gliders

The Lifting Frame is an essential piece of equipment for the safe lifting and handling of patients requiring spinal or orthopaedic immobility.




The specialised frame made from a durable metal comes complete with 12 removable nylon gliders and can be dismantled for ease of storage. Strap kit sold separately. \*Pictured being utilised with 4-point yoke. SWL 200kg.

#### SPECIFICATIONS:

Product Code 227201



# Pressure Care

Positioning Roll	22	Relax Easy Cushion	25
Molly Heel Protection	22	Static Air Cushion	25
Positioning & Support Wedge	22	 <b>HIGH RISK</b>	
Multi Pillow	23	Optimal 5zon Mattress	26
Reusable Polymer Elbow/ Heel Protector Pair	23	CuroCell® SAM BAS PRO Air Mattress	26
Action Polymer Elbow/ Heel Protector Pair	23	CuroCell® Area Zone Mattress	27
 <b>LOW RISK</b>		CuroCell® A4 CX20 Air Mattress	27
Relax Trio Cushion	24	CuroCell® UNO Air Mattress	28
 <b>MEDIUM RISK</b>		Action Gel Segment Bed Pad	28
Action Gel Pilot Cushion	24		

# PRESSURE CARE



## Positioning Roll

The positioning roll can give support to the resident in sitting and lying positions. It can be placed in various ways to provide support and relief.

- Supports shoulders, spine, hips, buttocks and heels
- All-round support in seated and reclining positions
- Relieves the knees and supports feet
- Supports in lateral position



### SPECIFICATIONS:

<b>Product Code</b>	32307
<b>Width</b>	150 mm
<b>Length</b>	600 mm

### SIZES AVAILABLE:

<b>Product Code</b>	32308	32309
<b>Width</b>	200 mm	200 mm
<b>Length</b>	1250 mm	2250 mm

## Molly Heel Protection

- Easy to use
- Equalizes pressure
- Faster healing
- Improves circulation
- Makes it easier to freely move the legs
- Prevents crossing the legs after hip surgery
- Soft and comfortable
- Wipe with cleaning agent and/or disinfect
- Machine wash max. 95°C
- Easy to clean
- Liquid resistant
- Long durability



Example of positioning



### SPECIFICATIONS:

<b>Product Code</b>	32311
<b>Leg Length</b>	780 - 850 mm
<b>Thickness</b>	75 mm

### SIZES AVAILABLE:

<b>Product Code</b>	32312	32313
<b>Leg Length</b>	860 - 940 mm	1020 - 1030 mm
<b>Thickness</b>	90 mm	110 mm

## Positioning & Support Wedge

The Support Wedge is used to place a user in side position in order to relieve pressure on the sacrum, which is an extra sensitive area. For the side position, two support wedges are used together.

The support wedge is easy to use and easy to clean. It is delivered with a liquid resistant and vapour permeable cover. The cover can be wiped off with an alcohol-based cleaning agent, as well as machine washed.

- Machine wash max. 95°C, tumble dry.
- Wipe off with cleaning agent, disinfect with alcohol-based cleaning agent.
- 100% HR-foam
- Polyurethane 2-way stretch polyester

### Examples of different uses:

- When the user has problems to lie completely on the side
- To prevent pressure ulcers
- When the user has problems with position changes
- For other positioning



Example of positionings

### SPECIFICATIONS:

<b>Product Code</b>	3689
<b>Width</b>	260 mm
<b>Length</b>	700 mm
<b>Thickness</b>	50 - 150 mm

## PRESSURE CARE



### Multi Pillow

The Multi Pillow is a versatile pillow that can be placed between the knees or under legs and heels for support and relief.

- Supports the feet, counters heel pressure against surfaces
- Locked position in connection with spasticity
- Supports in lateral position
- Relief for heels
- Supports under the knees



### SPECIFICATIONS:

Product Code	32302
Width	450 mm
Length	450 mm

### Reusable Polymer Elbow/ Heel Protector Pair

A comfortable, lightweight, and effective solution for zero heel pressure, the Pedisafe is designed for heel ulcer prevention and treatment. It also helps prevent footdrop and provides pressure relief. It comes reusable and is suitable in all healthcare environments.

- Lightweight design
- Helps prevent footdrop and relieves pressure on the heels
- Provides patient comfort and protects extremely fragile skin
- Open hole design for air & moisture ventilation
- Three adjustable straps allow custom fit



### SPECIFICATIONS:

Product Code	3131
--------------	------

### Action Polymer Elbow/ Heel Protector Pair

These gel protectors are for use in long term care settings and for assistance in wound management. They allow clients taking medication that affects skin integrity or who are prone to skin tears or abrasions to continue daily activities or to sleep comfortably.

They need only a wipe over with antiseptic solution or mild detergent to be kept fresh and odour free.

- Unique polymer consistency (not liquid gel) so bottoming out cannot occur
- Helps prevent footdrop



### SPECIFICATIONS:

Product Code	3170
Gel Thickness	48 mm

### SIZES AVAILABLE:

Product Code	3171
Gel Thickness	95 mm

# PRESSURE CARE




## Relax Trio Cushion


The Relax Trio is an ultralight cushion that benefits from four sections of FloGel being positioned into a polyurethane base, with a viscoelastic foam upper layer. It is designed to conform to the shape of the user, increasing surface area and relieving pressure.

- Viscoelastic slow
- Recovery Foam
- Fast Recovery Foam
- Includes outer cover (not pictured)
- Pressure Relieving FloGel (4 sections)
- Internal clear protective PU film coating for easy clean

\*Pressure Mapping Reports and successful Case Studies Available on request



Pressure diagram with Relax Trio cushion



Pressure diagram without wheelchair cushion



## SPECIFICATIONS:

Product Code	35332
Width	450 mm
Length	400 mm
Thickness	38 mm
SWL	80 kg

## SIZES AVAILABLE:

Product Code	35330	35333
Width	400 mm	450 mm
Length	400 mm	450 mm
Thickness	38 mm	38 mm
SWL	80 kg	80 kg



## Action Gel Pilot Cushion

The Pilot is a low-profile cushion that provides exceptional pressure relief for most users with a pressure ulcer history. 25mm of Akton polymer provides immersion for bony prominences and protects clients with fragile skin integrity. Incontinence covers are an optional extra

- Not restricted by a weight limit
- Easy to care for
- Quality, durable Pilot cover included - can be machine washed

\*Custom sizes also available.




## SPECIFICATIONS:

Product Code	3531
Width	450 mm
Length	400 mm
Thickness	25 mm
SWL	Unlimited

## SIZES AVAILABLE:

Product Code	3530	3532	3533
Width	400 mm	450 mm	500 mm
Length	400 mm	450 mm	450 mm
Thickness	25 mm	25 mm	25 mm

## PRESSURE CARE




### Relax Easy Cushion


This cushion from the Relax range is easy to use and is an ideal pressure relieving cushion for use in the nursing home environment. This combination of pressure relieving FloGel and viscoelastic foam facilitates effective pressure re-distribution while allowing the patient to immerse themselves comfortably into the cushion. The cushion has no specific front or back, top or bottom, making it simple to use.

- Conforms to shape of user, increasing surface area and relieving pressure
- Includes waterproof wipeable PU stretch outer cover
- Outer cover machine washable to 60° C
- Pressure Relieving FloGel - Segmented For Added Conformity

\*Pressure Mapping Reports and successful Case Studies Available on request



Pressure diagram with Relax Easy cushion



Pressure diagram without wheelchair cushion



### SPECIFICATIONS:

Product Code	35312
Width	450 mm
Length	400 mm
Thickness	50 mm
SWL	130 kg

### SIZES AVAILABLE:

Product Code	35100	35313	35315	353120
Width	400 mm	450 mm	500 mm	500 mm
Length	400 mm	450 mm	450 mm	500 mm
Thickness	50 mm	50 mm	50 mm	50 mm
SWL	130 kg	130 kg	130 kg	130 kg



### Static Air Cushion

CurroCell® S.A.C. is a static air seat cushion for comfort, prevention and as an aid in the treatment of pressure injury up to category 2.

The cushion pre-filled and very easy to handle. The durable soft material for long term use provides a superior combination of comfort and support.

- Wipe with cleaning agent and/or disinfectants
- Air holes for venting heat and moisture
- Pre-inflated low profile design
- Fits all chairs
- Durable soft material for long term use
- Patient ready

\*Pressure Mapping Reports and successful Case Studies available on request



### SPECIFICATIONS:

Product Code	3417	3418
Width	450 mm	500 mm
Length	450 mm	500 mm
SWL	160 kg	160 kg

## PRESSURE CARE



### Optimal 5zon Mattress

The Optimal 5zon® is ideal for all healthcare environments and is used as an aid in the prevention and treatment of pressure ulcers up to and including category 2, and for pain therapy. It is a soft and flexible mattress replacement that is placed directly on the bed and replaces both the standard mattress and mattress pad.

- Highly elastic, high-quality material and a unique zone and module system gives effective pressure equalisation.
- Specially designed to minimise shear forces
- Simple and effective cleaning, cover has a very well documented hygiene standard
- The hygienic cover is vapour permeable, which reduces the risk of skin maceration
- The patented sit function DuoCore® is designed to give support when sitting and also to provide good comfort when lying
- Heel Protection provides comfort for relief of heel. The user can move the leg and make position changes. Comes with a removable liquid resistant hygiene cover.



#### SPECIFICATIONS:

Product Code	3771
Width	850 mm
Length	2000 mm
Thickness	140 mm
SWL	180 kg

#### SIZES AVAILABLE:

Product Code	3778	3775
Width	1000 mm	1050 mm
Length	2000 mm	2000 mm
SWL	180 kg	180 kg



### CuroCell® SAM BAS PRO Air Mattress

SAM BAS PRO is used as an aid in the prevention and treatment of pressure ulcers up to category 3, and for pain therapy. It is a mattress replacement system with an integrated foam core. The mattress is made of very soft and durable materials. The included support pump is easy to use and very quiet. The mattress comes with CPR (rapid deflation of air) and has a built-in heel function.

- This mattress is supplied with a four way stretch hygienic outer cover which conforms to relevant standards for hygiene, quality and safety.
- Silent, energy-saving pump
- Soft, elastic, vapour permeable hygienic cover
- Easy to handle and maintain
- Integrated heel function

\*Clinical Test Report Available on request

For more sizes please contact the Active Medical team on 1800 267 267



#### SPECIFICATIONS:

Product Code	3683
Width	850 mm
Length	2000 mm
Thickness	160 mm
SWL	250 kg

#### SIZES AVAILABLE:

Product Code	3684	3685
Width	1050 mm	1000 mm
Length	2000 mm	2000 mm
SWL	250 kg	250 kg



## PRESSURE CARE



### CuroCell® Area Zone Mattress

The AREA Zone is a mattress that combines air with foam for effective pressure redistribution and great comfort. The mattress is self-adjusting and non-powered. It is a flexible mattress replacement system that can be used in all kinds of healthcare settings, easy to use and handle. The mattress system is used as an aid in prevention and treatment of pressure injuries up to category 3.

- Zone castellations in top layer foam provides reduction of shearing forces
- Handle for transportation
- Flexible and self-adjusting internal air cells, no pump required
- Integrated heel function
- Reinforced cover with welded seams
- Non-powered mattress replacement system

\*Clinical Test Report Available on request



### CuroCell® A4 CX20 Air Mattress

The CuroCell® A4 CX20 provides simplicity together with user comfort while effectively aiding in prevention and treatment of pressure injuries up to and including category 4, and unclassified pressure ulcers. It provides features for excellent user comfort, reduction of shear forces and improvement of the microclimate.

The CX20 mattress is 20 cm thick and the alternation of the cells is every third cell down ensuring 75% of the body is being supported at one time. The pump is intelligent and fully automatic without need for manual settings. It identifies the user and automatically manages all settings to optimise the pressure according to the user's needs.

All mattresses are supplied with a four way stretch hygienic cover which meets standards for hygiene, quality and safety. All covers are manufactured from materials which reduce the risk of skin maceration.

- Unique CPR-function with transport function
- Heel Function
- Audible and visual notifications
- Silent running and vibration free
- Three program modes: Pulsating, Alternating and CLP (constant low pressure)
- Individually replaceable air cells with safety mattress
- Smart cable management

For more sizes please contact the Active Medical team on 1800 267 267



#### SPECIFICATIONS:

Product Code	3678Z
Width	850 mm
Length	2000 mm
Thickness	180 mm
SWL	230 kg

#### SIZES AVAILABLE:

Product Code	3686Z
Width	1000 mm
Length	2000 mm
SWL	230 kg



#### SPECIFICATIONS:

Product Code	3641
Width	850 mm
Length	2000 mm
Thickness	200 mm
SWL	250 kg

#### SIZES AVAILABLE:

Product Code	3642	3643
Width	1000 mm	1050 mm
Length	2000 mm	2000 mm
SWL	250 kg	250 kg

## PRESSURE CARE



### CuroCell® UNO Air Mattress

The CuroCell UNO® is used as an aid in the prevention and treatment of pressure injuries up to and including category 4. It is an alternating pressure mattress. It is equipped with a large touch screen panel with a user-friendly interface. The fully automated pump offers automatic adjustment, providing optimal pressure relief and comfort for each individual user while aiding simple use. The mattress has a four way stretch hygienic outer cover which conforms to relevant standards for hygiene, quality and safety. All covers are manufactured from modern, recyclable materials and are vapour permeable, which reduce the risk of skin maceration.

- Fully automated pump that adapts individually with large colour touch display
- Cell on Cell construction with safety mattress
- Unique design for easy handling and usage
- Silent running and vibration free
- Heel Function
- Smart cable management
- 3 modes: Alternating, Pulsating and Static



#### SPECIFICATIONS:

Product Code	3651
Width	850 mm
Length	2000 mm
Thickness	200 mm
SWL	250 kg

#### SIZES AVAILABLE:

Product Code	3652
Width	1050 mm
Length	2000 mm
SWL	250 kg



### Action Gel Segment Bed Pad

The Action bed overlays are made of Akton polymer which mimics human fatty tissue in consistency, does not migrate like liquid gel and provides equally distributed pressure management. With varying sizes available, these overlays can be used to cover the entire surface or isolated areas of an existing standard mattress. Clients can use one, two or three segments in combination on their areas of need e.g shoulders and hips, feet and hips or full body etc.

- Easily cleaned, no maintenance
- Not restricted by any weight limit
- 3 Segments cover length of single bed



#### SPECIFICATIONS:

Product Code	3853
Width	630 mm
Length	690 mm
Thickness	180 mm
SWL	250 kg





## WHEELCHAIRS



### Thrifty Wheelchair

The Thrifty is a sturdy and reliable standard weight wheelchair, with removable armrests and leg rests, built for individuals, nursing homes and non-prescription requirements. This wheelchair meets Australian Standards for Manual Wheelchairs in Product Performance. All of these features combine to make the Thrifty of unbeatable value in standard wheelchairs.

- Swing away, nonremovable armrests
- Swing away, detachable legrests
- 24" rear wheels with urethane tyres ( maintenance free )
- Australian Standard AS/NZS 3695.1:2011 approved

#### SPECIFICATIONS:

Product Code	4246
Seat Width	450 mm
Seat Depth	400 mm
Overall Width	650 mm
Overall Depth	1080 mm
Overall Height	880 mm
Fold Width	300 mm
Fold Length	790 mm
SWL	115 kg



### Thrifty Transit Wheelchair

The Thrifty Transit is an excellent chair designed to suit the person who cannot propel themselves. This folding chair has removable armrests, removable swing away leg rests and 12.5" rear wheels. Ideal for transporting in the boot of a car.

- Removable armrests
- Swing away, detachable legrests
- 12½" rear wheels with urethane tyres
- Powdercoated steel frame
- Padded nylon upholstery

#### SPECIFICATIONS:

Product Code	4146
Seat Width	460 mm
Seat Depth	400 mm
Seat Height	450 mm
Backrest Height	430 mm
Armrest Height	230 mm
Overall Width	630 mm
Overall Depth	980 mm
Overall Height	905 mm
Fold Width	280 mm
Net / Gross Weight	12 - 16 kg
SWL	100 kg

## WHEELCHAIRS




## Deluxe Alloy Wheelchair

The Deluxe Ultralight Alloy Wheelchairs have been designed to accurately meet the needs of our clients. They are an ultralight weight wheelchair with the benefit of many additional features. The improvement compared to previous designs (ranging from 25 - 30 kg gross weight) is considerable. Much input from clients and nurses has assisted in the development of these superior quality wheelchairs.

- Flip back armrests
- Swingaway, detachable legrests
- Fold down backrest
- Quick release rear wheels
- 24" rear wheels with urethane tyres
- Aluminium frame
- Padded nylon upholstery
- Complies with AS3695:2008



Height Adjustable armrests



Urethane Tyres



No Upholstery Screws



Flip back Armrests

## SPECIFICATIONS:

<b>Product Code</b>	4629
<b>Seat Width</b>	460 mm
<b>Seat Depth</b>	400 mm
<b>Seat Height</b>	480 mm
<b>Backrest Height</b>	430 mm
<b>Armrest Height</b>	230 mm
<b>Overall Width</b>	650 mm
<b>Overall Depth</b>	1060 mm
<b>Overall Height</b>	960 mm
<b>Fold Width</b>	320 mm
<b>Net/Gross Weight</b>	11 / 14 kg
<b>SWL</b>	125 kg

## SIZES AVAILABLE:

<b>Product Code</b>	45129	45629
<b>Seat Width</b>	510 mm	560 mm
<b>Seat Depth</b>	400 mm	400 mm
<b>Seat Height</b>	490 mm	490 mm
<b>Backrest Height</b>	430 mm	430 mm
<b>Armrest Height</b>	210 - 280 mm	210 - 280 mm
<b>Overall Width</b>	700 mm	750 mm
<b>Overall Depth</b>	1000 mm	1000 mm
<b>Overall Height</b>	960 mm	960 mm
<b>Fold Width</b>	300 mm	300 mm
<b>Net/Gross Weight</b>	14.5 / 17 kg	15 / 17.5 kg
<b>SWL</b>	150 kg	150 kg

## WHEELCHAIRS



### Hospital Ward Alloy Wheelchair

The Hospital Ward Alloy Wheelchair has been designed to accurately meet the needs of our patients. Much input from hospital staff has assisted in the development of these superior quality wheelchairs.

- Swingaway, detachable legrests (if non-removable legrests are required this can be done at time of order)
- Height adjustable, flip back armrests - this assists in side transfers of the client
- IV Pole and oxygen bottle carrier
- Anti tip bars
- 12" rear wheels with urethane tyres - no flat tyres
- Aluminium construction - less force to push
- Padded nylon upholstery - can be wiped over to clean
- Complies with AS3695:2008

#### SPECIFICATIONS:

<b>Product Code</b>	45128WD
<b>Seat Width</b>	510 mm
<b>Seat Depth</b>	400 mm
<b>Seat Height</b>	490 mm
<b>Backrest Height</b>	430 mm
<b>Armrest Height</b>	210 - 280 mm
<b>Overall Width</b>	700 mm
<b>Overall Depth</b>	1000 mm
<b>Overall Height</b>	960 mm
<b>Fold Width</b>	350 mm
<b>Net/Gross Weight</b>	15 / 17 kg
<b>SWL</b>	150 kg



### Gusto Heavy Duty Wheelchair

The Gusto heavy duty wheelchair caters for the bariatric client. Featuring padded vinyl upholstery, heavy duty folding frame, 24 inch rear wheels and inset foot plates for client comfort, this wheelchair is ideal in the aged care and hospital environment. It comes standard with urethane tyres, eliminating the need to maintain tyre pressure. The swing away leg rests are removable and also adjust in length. The armrests will remove to enable side transfers if necessary. Footplates are alloy for longevity and ease of cleaning.

- Removable, swingaway, adjustable length legrests
- Removable armrests
- Padded vinyl upholstery
- 24" rear mag type wheels, urethane tyres

#### SPECIFICATIONS:

<b>Product Code</b>	46108
<b>Seat Width</b>	610 mm
<b>Seat Depth</b>	460 mm
<b>Seat Height</b>	520 mm
<b>Backrest Height</b>	430 mm
<b>Armrest Height</b>	220 mm
<b>Overall Width</b>	850 mm
<b>Overall Depth</b>	1100 mm
<b>Overall Height</b>	930 mm
<b>Fold Width</b>	350 mm
<b>Net/Gross Weight</b>	26 / 30 kg
<b>SWL</b>	240 kg



## SEATING



### Regency Care Chair

The 5050G Regency Care Chair is designed to give the client comfort by combining the original tilt system with an independent multi-positional back rest. Its many features take into account the needs of both carer and user of the chair. The Regency 5050G derives its status from a proven and successful track record across a broad spectrum of the industry.

- **Other sizes available.** Contact the team today to find out more.
- **Drop down armrests and side wings:** Regency armrests/wings fold right away to allow chair to be positioned close to bed for supine transfers. Also making sling insertion much easier.
- **Manoeuvrability:** Regency chairs have two front directional lock castors, two locking rear castors and a horizontal push handle reducing the effort to manoeuvre.
- **Gel cushion:** The guesswork has been removed ensuring the seat has the correct pressure to support the resident with the highest level of comfort.
- **Leg and footrest foldaway feature:** The uniquely designed non removable footplates prevents persons from sliding forward lowering the number of repositions required.
- **Manoeuvrable castor base:** 175mm castors reduces effort to move the chair.

Colours available:



### SPECIFICATIONS:

<b>Product Code</b>	5050G
<b>Seat Width</b>	460 mm
<b>Seat Depth</b>	520 mm
<b>Seat Height</b>	600 mm
<b>Backrest Length</b>	700 mm
<b>Armrest Height</b>	175 mm

<b>Legrest Length</b>	480 - 650 mm
<b>Overall Width</b>	660 mm
<b>Overall Depth</b>	1040 mm
<b>Overall Height</b>	1310 mm
<b>SWL</b>	180 kg



## SEATING




### Ashley Mobile Recliner - Standard

Economically priced, the Ashley features a waterfall style backrest, sculptured armrests, side pockets and padded leg chaise. Australian made. Cascading roll back with isolated head rest and lumbar roll to support both these areas.

Available in either a fabric or vinyl upholstery in a pattern or plain. Contact our team to discuss! Vinyl available in Austex Studio Encore colour range.

- Allow +/- 20mm tolerances

**Other colours available:**



Magma Wine    Glamour Scarlett    Deluxe Chambray    Topsuede Denim



Lovely Jade    Regis Ginger    Mazza Charcoal    Jazz Walnut

Available in either a fabric or vinyl upholstery in a pattern or plain.  
Contact our team to discuss!

### SPECIFICATIONS:

Product Code	5521
Size	Medium
Seat Width	480 mm
Seat Depth	500 mm
Seat Height	480 mm
Backrest Height	640 mm
Overall Width	780 mm
Overall Height	1050 mm
SWL	120 kg



### Oasis Mobile Recliner with Tray

Generously padded all over, and with a gas cylinder mechanism allowing infinite recline positioning, the Oasis recliner is a treat to relax in. Equally important to patient comfort is ease of operation for the care-giver. When developing the Oasis Recliner we gave extra consideration to this, which is evidenced in the features below.

- Ergonomically placed handles enable easier patient transfer
- Heavy duty design throughout for extra strength and product longevity
- Arm cap designed to assist patient ingress and egress
- Infinite recline positioning
- Side lever permits both patient and caregiver to operate leg extension
- Geometrically balanced mechanism substantially reduces effort required to operate leg extensions



### SPECIFICATIONS:

Product Code	5813
Seat Width	584 mm
Recline Depth	1880 mm
Seat Height	520 mm
Arm Height	692 mm
Overall Width	902 mm
Overall Height	1219 mm
SWL	225 kg

## SEATING



### Embark Deluxe Extracare Day Chair

The Extra Care Embark Deluxe has an impressive 250kg weight capacity and extra wide 60cm width, making it ideal for the heavier client. It still provides all the benefits of the Embark Deluxe, e.g. adjustable seat depth and height, and comes standard with tilt and move wheels.

- High backrest incorporating lumbar and lateral support
- Padded armrests
- Comfortable vinyl upholstered seat and back
- Adjustable seat height and depth
- Powder-coated steel frame with reinforcement for increased load

#### SPECIFICATIONS:

<b>Product Code</b>	5931
<b>Seat Width</b>	600 mm
<b>Seat Depth</b>	450 - 530 mm
<b>Seat Height</b>	470 - 620 mm
<b>Backrest Height</b>	765 mm
<b>Arm Height</b>	220 mm
<b>Overall Width</b>	800 mm
<b>Overall Depth</b>	740 - 820 mm
<b>Overall Height</b>	1180 - 1330 mm
<b>SWL</b>	250 kg



### Utilite Extra Care Day Chair

The Extra Care Utilite is a heavy duty, multi purpose adjustable seat height and angle chair. Ideal for use in an area where adjustable seating is required, for example, lounge, dining room, bedroom, activity room, etc. This low back chair provides comfort, stability and practicality at an affordable price.

- Extra wide frame – 550mm seat width
- Shallow curved low backrest
- Padded armrests
- Comfortable vinyl upholstered seat and back
- Adjustable seat height to accommodate any user
- Powder-coated steel frame

#### SPECIFICATIONS:

<b>Product Code</b>	5911
<b>Seat Width</b>	550 mm
<b>Seat Depth</b>	430 mm
<b>Seat Height</b>	450 - 630 mm
<b>Backrest Height</b>	450 mm
<b>Arm Height</b>	230 mm
<b>Overall Width</b>	670 mm
<b>Overall Depth</b>	420 mm
<b>Overall Height</b>	850 - 1030 mm
<b>SWL</b>	200 kg



## SEATING



### Embark Superior Day Chair

The Embark Superior Day Chair was designed and developed by a team of clinicians. It has a strong frame and a high weight capacity of 200 kg. The chair has low resilient foam and 2-way stretch upholstery on seating surfaces. For increased longevity in demanding environments it has 6cm tilt and move wheels and other small improvements.

There have been instances in the past, whereby patients, staff and visitors have sustained cuts and injuries from square tube frames. As a result, the Embark was not only designed for comfort, but to avoid these injuries.

- Low back option also available.
- Comfortable upholstered seat and back
- Padded armrests
- Adjustable seat height and depth
- Powder-coated steel frame
- High backrest incorporating lumbar and lateral support
- Available in Grey and Periwinkle



Periwinkle



Grey

#### SPECIFICATIONS:

<b>Product Code</b>	5933
<b>Seat Width</b>	500 mm
<b>Seat Depth</b>	450 - 530 mm
<b>Seat Height</b>	450 - 600 mm
<b>Backrest Height</b>	790 mm
<b>Arm Height</b>	230 mm
<b>Overall Width</b>	700 mm
<b>Overall Depth</b>	730 - 800 mm
<b>Overall Height</b>	1200 - 1350 mm
<b>SWL</b>	200 kg



### Utilite Day Chair

The Utilite is a multi purpose adjustable seat height and angle chair. Ideal for use in an area where adjustable seating is required, for example, lounge, dining room, bedroom, activity room, outdoor area etc. This low back chair provides comfort, stability and practicality at an affordable price.

- Shallow curved low backrest
- Padded armrests
- Comfortable, durable upholstered seat, back and armrests
- Powdercoated steel frame



Blue



Grey



Green Slate



Champagne

#### SPECIFICATIONS:

<b>Product Code</b>	5910
<b>Seat Width</b>	460 mm
<b>Seat Depth</b>	460 mm
<b>Seat Height</b>	485 - 640 mm
<b>Backrest Height</b>	455 mm
<b>Arm Height</b>	230 mm
<b>Overall Width</b>	610 mm
<b>Overall Depth</b>	550 mm
<b>Overall Height</b>	900 - 1075 mm
<b>SWL</b>	120 kg



## MOBILITY SOLUTIONS



### Crutches - Bariatric Height Adjustable

The Bariatric Crutches offer push button height adjustment, providing support, strength and comfort for consumers up to 250kg.

- 3 sizes available. (Codes: 6120, 6121, 6122)
- Supplied as a pair

#### SPECIFICATIONS:

Product Code	6020	6021	6022
<b>Crutch Height</b>	940 - 1150 mm	1130 - 1340 mm	1330 - 1540 mm
<b>Patient Height</b>	1219 - 1372 mm	1372 - 1575 mm	1575 - 1778 mm
<b>SWL</b>	250 kg	250 kg	250 kg



### Crutches - Anatomic Grip

The Anatomic Elbow Crutches feature contoured moulded handles that evenly distribute pressure across the palm of the hand, contributing to the users ease of operation and comfort. Added to this is the ability to adjust the crutch length from cuff to hand grip and also from hand grip to ground, enabling the user to adjust the crutch to their individual requirements. Manufactured from a larger tubing giving increased stability and weight capacity.

- Supplied as a pair.
- 3 sizes available (Codes: 6150, 6151J, 6152J)

#### SPECIFICATIONS:

<b>Product Code</b>	6150
<b>DVA</b>	410L
<b>Cuff to Hand Grip</b>	790 - 945mm
<b>Hand Grip Height</b>	1219 - 1372 mm
<b>SWL</b>	100 kg

#### SIZES AVAILABLE:

Product Code	6151J	6152J
<b>DVA</b>	410L	410L
<b>Cuff to Hand Grip</b>	930 - 1130 mm	1335 - 1540 mm
<b>Hand Grip Height</b>	1372 - 1575 mm	1777 - 1981 mm
<b>SWL</b>	180 kg	180 kg



### Staebel Unilite Forearm Walker

The Unilite Forearm Walker is designed to support and tutor the user in an upright gait position. The very stable base gives added confidence and security to the user. Features include padded gutter supports that are adjustable in the forearm length, height and width. The popular folding design allows the Unilite to fold easily for transport and storage. This model is fitted with four large wheels and lockable cable brakes.

- Height adjustable ergonomic shaped hand grips
- Powder coated steel construction
- Lockable cable type brakes
- 200mm castors

#### SPECIFICATIONS:

<b>Product Code</b>	6742
<b>DVA</b>	BRO213
<b>Seat Width</b>	400 mm
<b>Seat Height</b>	600 mm
<b>Handle Height</b>	890 - 1100 mm
<b>Walker Weight</b>	14 kg
<b>SWL</b>	100 kg

## MOBILITY SOLUTIONS



### Stimulator Rehab Forearm Support Walker

“The best rehab program and the best rollator will guarantee you the best results,” – is this what you’re looking for? Our Stimulator Rehab Forearm Walker has been specifically designed for the needs of the hospital environment.

- Very lightweight, easy for the frail, weak patient to use
- High quality castors with super smooth rolling and minimal friction even on carpet, taking away any shuddering or wobbling
- Easy, simple and very gentle lever height adjustment
- Soft PU multi adjustable arm pads for any patient with any specific position requirement
- Multi adjustable hand grips, for the arthritic, ABI or stroke patients.
- Four swivelling castors for ultimate manoeuvrability
- Individual brakes on each castor, as well as hand brakes.

#### SPECIFICATIONS:

<b>Product Code</b>	6755
<b>Arm Support Height</b>	905 - 1310 mm
<b>Width Between Arms</b>	310 - 470 mm
<b>Overall Width</b>	700 mm
<b>Length</b>	810 mm
<b>SWL</b>	200 kg



### Cruiser Folding Walker

The Cruiser is a popular seat walker manufactured from durable lightweight aluminium. 6” castors make this an ideal outdoor walker. But is also compact enough for use indoors. Features include a padded seat, bag under seat, lockable brakes and padded, folding backrest for comfort while seated. The adjustable height handles are fitted with orthopaedic hand grips, therefore providing the user with excellent support. The walker also folds flat for transport.

- Height adjustable ergonomic shaped handgrips
- Lightweight aluminium construction
- Lockable handbrakes
- 152mm castors

#### SPECIFICATIONS:

<b>Product Code</b>	6730
<b>DVA</b>	2907
<b>Seat Width</b>	440 mm
<b>Seat Height</b>	560 mm
<b>Handle Height</b>	820 - 930 mm
<b>Walker Weight</b>	6 kg
<b>SWL</b>	125 kg



### Ligero Seat Walker

The Ligero is a popular seat walker manufactured from durable lightweight aluminium. The larger castors make this an ideal outdoor walker, however is also compact enough for use indoors. Features include a padded seat and padded folding backrest for comfort while seated. The adjustable height handles are fitted with Orthopaedic hand grips therefore providing the user with excellent support. The walker also folds for storage and transport.

- Height adjustable ergonomic shaped hand grips
- Lightweight aluminium construction
- Lockable handbrakes
- 200mm (8”) castors

#### SPECIFICATIONS:

<b>Product Code</b>	6865
<b>DVA</b>	2967
<b>Seat Width</b>	440 mm
<b>Seat Height</b>	550 mm
<b>Handle Height</b>	790 - 910 mm
<b>Walker Weight</b>	7.5 kg
<b>SWL</b>	120 kg

#### SIZES AVAILABLE:

<b>Ligero Seat Walker - 510 mm</b>	
<b>Product Code</b>	6869



## BEDROOM SOLUTIONS



### Eleganza 1 Bed

With full electrical adjustment, the Eleganza 1 bed is designed for the standard hospital ward. The bed has a simple construction, a wide range of features and functions and at a low cost, the Eleganza 1 provides customers with great value for money. Its sophisticated electric functions include the safety concept of the 'GO' button, protecting the patient against an injury caused by any unintentional movement of the bed; nurse control panel with pre-programmed cardiac chair, CPR and trendelenberg settings; ergonomic patient handset.

Side guard protection options include standard collapsible side rails, or 4-part side guards which also have the option of the patient control integrated.

The Eleganza 1 uses the Ergoframe, which enables the mattress platform to widen the space under the patients buttocks for optimum comfort and posture whilst eliminating the negative effects of compression, shear and friction

forces which can lead to the development of pressure ulcers. The bed can be equipped with a number of optional functions and useful versatile features and will meet the specific needs of every hospital department.

Nosocomial infection prevention is achieved through joint-free construction and easy cleaning. Safety is paramount, with the Eleganza 1 bed meeting all of the requirements of the latest standard IEC 60601-2-52:2010.

- 4-section mattress platform
- Height range 395 – 768mm
- Range of colours to co-ordinate with your interior
- European construction

### SPECIFICATIONS:

<b>Product Code</b>	73014
<b>External Dimensions</b>	990 x 2180 mm
<b>Bed Extension</b>	150 mm
<b>Mattress Platform</b>	4 section electric adjustable
<b>Mattress Surface</b>	900 x 2000 x 140 mm
<b>Electric Height Adjustment</b> (125mm Castors)	370 - 735 mm
<b>Electric Height Adjustment</b> (150mm Castors)	395 - 768 mm
<b>Electric Backrest Adjustment</b>	0 - 70°

<b>Electric Thighrest</b>	0 - 39°
<b>TR / ATR</b>	15° - 15°
<b>Ergoframe®</b>	100 + 60 mm
<b>Bed Weight</b>	135 kg
<b>Standard</b>	IEC:60601-2-52
<b>Cover Against Water and Dust</b>	IP54
<b>Clearance for Server</b> (125mm Castors)	160 mm
<b>Clearance for Server</b> (150mm Castors)	1880 mm
<b>SWL</b>	250 kg



## BEDROOM SOLUTIONS





### Sentida 6 Bed

Known for its ambient non-institutional appearance the Sentida's features include:

- Compliance with standard IEC 60601-2-52
- Unique Multi-position Mobilisation Aid which have 4 positions and very useful for effective mobilisation as they are easy to grab and clasp. They can be used as a safe standing-up aid and verifiably support the resident's mobility. Easy to adjust and are hidden behind a side board when lowered.
- Great movability with 100mm easy-roll castors at both bed ends
- Foot-end bar brake with 4 castor central braking system
- Minimal gap between edge of mattress surface and side guard
- Double retracting sleeping surface to reduce friction & shearing forces. Also will aid in relief of pressure in sacral area
- Unique comfort seating position to make patient breathing easier, as well as digestion
- Vertical raise therefore no risk of holes in walls behind the bed
- Underbed clearance for patient lifters at any position
- Handset lockout function so care staff have control in dementia care setting

- Unique hand control with 2 button operation reducing risk of accidental bed adjustment
  - TGA approved
  - Havana Cherry Decor (pictured), other colours available
  - 36 month electrical warranty
  - 10 month frame warranty
  - Quality German made construction
  - Low-low 3 stop bed height
  - Comfort and Relax Trendelenburg features
- \* Also available in King Single and LE version



(Picture for illustration purposes only)

### SPECIFICATIONS:

<b>External Dimensions</b>	2080 x 1060 mm
<b>Sleeping Surface</b>	2000 x 850 mm
<b>Bed Extensions</b>	200 mm
<b>Maximum Sleeping Surface Height</b>	800 mm
<b>Minimum Sleeping Surface Height</b>	270 mm
<b>SWL</b>	270 kg

### SIZES AVAILABLE:

**Sentida 6 Bed, Low Height, King Single**

<b>Product Code</b>	77713C
<b>SWL</b>	270 kg

### USE WITH:

<b>Product Code</b>	77607
---------------------	-------

## BEDROOM SOLUTIONS



### Staebel Floorline Bed

Introducing the first ever universal bed that suits low care, high care and bariatric care. Plus, this innovative three-in-one bed has a low height of 100mm and NO TOWERS! Suitable for care of frail elderly or dementia patients, through to acute and bariatric settings.

- Single Width
- Electric Hi/Low
- Auto Regression Backrest
- Electric Kneebend
- Electric Forward & Reverse Trendelenburg
- Standard Head and Foot Boards - Fini Blonde
- 75mm Twin Wheel Castors
- SWL 200kg
- Low Position - 120mm | High Position - 720mm
- Underbed Clearance - Clearview
- Angles Back - 62° Knee - 35°
- Trendelenburg - 15° Reverse - 15°
- Overall Dimensions 2080mm x 905mm
- Sleeping Surface 2000mm x 890mm
- 10 year frame / 5 year electric warranty / 1 year Handset
- Various colours and styles available
- Also Available in King Single



### SPECIFICATIONS:

<b>Product Code</b>	7770S	<b>Back Rest Retraction</b>	115 mm
<b>Sleeping Surface Width</b>	870 mm	<b>Undercarriage Clearance</b>	138 mm
<b>Sleeping Surface Length</b>	2000 mm	<b>External Width</b>	1060 mm
<b>Min Sleeping Surface Height</b>	270 mm	<b>External Length</b>	2080 mm
<b>Max Sleeping Surface Height</b>	800 mm	<b>Max User Weight</b>	215 kg
<b>Thigh Rest Retraction</b>	60 mm	<b>SWL</b>	250 kg

## BEDROOM SOLUTIONS



### Image 3 Bed

Image 3, electrically adjustable bed respects the growing demands of the patient and nursing staff. Its attractive design will find universal use in a whole spectrum of hospital care. Image 3 deals with the patient's safety and mobilisation in a unique way: combining the SafeFree® siderails and reduced low height.

Image 3 was designed with focusing on 100% safe recovery of patient in the bed. Safety is the key characteristics of Image 3. Reduced height combined with telescopic SafeFree® siderails create absolutely safe environment for the patient, that can be adjusted according to their specific needs.

Also included is a complex system of supportive automatic features called 3-STOP STRATEGY. It contains functions for safe sleeping, independent mobilisation for quicker recovery, convenient care and cleaning procedures for hospital staff. This bed contains extended features which prevent pressure injuries and contribute to their easier treatment. Mattress platform profiling

called Ergoframe® provides needed comfort for the patient.

Furthermore, this bed is designed with emphasis to maximally ease everyday activities such as nursing procedures, transfers and other routine tasks. Ergonomics, intuitive controls, smooth manipulation with bed and easy cleaning bring relief to hospital staff.

Image 3 is elegant and its simple design with clear lines impresses with pleasant, home looking atmosphere. Though it is very functional, it brings ergonomics and user-friendly handling. This way, the Image 3 can be easily used in any hospital area of care. Thanks to the functions and its timeless design it can be implemented to the nursing care as well. Numerous options of equipment and accessories will adjust to the specific needs of its owner.

An extra wide version, the Image 3 Bariatric is also available, with a 1m wide patient lying surface.

### SPECIFICATIONS:

<b>Product Code</b>	74800
<b>External Dimensions</b>	1040 × 2200 mm
<b>Clearance Height</b>	130 / 155 mm
<b>Mattress Dimensions</b>	2000 x 860 x 900 mm
<b>Mattress Platform Height Adjustment</b>	280 - 800 mm / 300 - 820 mm
<b>Bed Extension</b>	100, 220, 300 mm
<b>Height of Siderails (Split)</b>	400 mm
<b>Height of Siderails (Single Collapsible)</b>	390 mm

<b>Maximum Mattress Height (Split Siderails)</b>	180 mm
<b>Maximum Mattress Height (Single Collapsible Siderails)</b>	170 mm
<b>Bed Weight</b>	145 kg
<b>Backrest Angle</b>	0 - 70°
<b>Thighrest Angle</b>	0 - 34°
<b>TR / ATR Angle</b>	0 - 14°
<b>SWL</b>	265 kg



## BEDROOM SOLUTIONS



### Viva Mobile Overbed Table

The Viva Overbed Table offers the latest in smart and modern design, with an easy clean, seamless PE top to promote ease of use and satisfy infection control requirements. Featuring smooth lines, no joins or hard to reach contours, this innovative design makes previously hard to maintain items easier and quicker to clean.

- Easy gas lift height adjustment, non-tilt
- PE easy clean top, 34mm deep, max 20kg load.
- Base height 98mm
- Four castors, including two with brakes
- Supplied fully assembled

### SPECIFICATIONS:

<b>Product Code</b>	4210
<b>Top Width</b>	900 mm
<b>Top Depth</b>	430 mm
<b>Minimun Height</b>	725 mm
<b>Maximun Height</b>	1050 mm
<b>SWL</b>	20 kg



### Mobile Overbed Table

A very popular mobile table with a number of useful features, making it an ideal choice in many situations. With its easy lever operated height adjustment and smooth curved U-design, it fulfils all needs.

- Spring assisted lever height adjustment – lever stays with top eliminating need to bend to adjust table height
- Powdercoated U-shape base, round tube minimises injury risk
- Swivel castors, two are lockable

**\* Table top is sold seperately. There are a range of colours on offer**

### SPECIFICATIONS:

<b>Product Code</b>	4215 - B00
<b>Width</b>	450 mm
<b>Length</b>	750 mm

## BEDROOM SOLUTIONS



### Deluxe Commode Bedside

This bedside commode offers a full range of features including a slide back seat with durable toilet seat and pan rack. This comes standard with bowl and lid. A skirt that surrounds the seat conceals the bowl etc. so the chair can be likened to a regular bedroom sitting chair

- Adjustable seat height
- Comfortable upholstered slide back seat
- Comfortable upholstered shallow curve backrest
- Plastic toilet seat, bowl and lid
- Padded armrests
- Powdercoated steel frame
- Vinyl upholstery

### SPECIFICATIONS:

<b>Product Code</b>	7955
<b>Overall Width</b>	600 mm
<b>Overall Height</b>	870 - 1020 mm
<b>Overall Depth</b>	510 mm
<b>Seat Width</b>	460 mm
<b>Seat Height</b>	420 - 570 mm
<b>Seat Depth</b>	440 mm
<b>Back Height</b>	470 mm
<b>SWL</b>	120 kg

### SIZES AVAILABLE:

Product Code	Width	SWL
7960	500 mm	135 kg
7961	550 mm	175 kg
7962	600 mm	195 kg
7963	650 mm	210 kg



### Triangle With Retractable Belt

The triangle can assist patients when in bed. It is a unique retractable design to avoid the triangle hitting clients head and/or restricting clear vision when bed backrest has been raised.

For use in conjunction with Overhead pole.

### SPECIFICATIONS:

<b>Product Code</b>	77605
---------------------	-------


### USE WITH:

<b>Product Code</b>	77692
---------------------	-------





## PATIENT SAFETY



### Cura1 Cordless Floor Mat

- 1 year life expectancy
- Simple operation and maintenance
- Easy to clean
- For use with cura1 Premium Falls Monitor
- Virtually no false alarms over the lifetime of the pad
- All bed/chair pads made from non-allergenic soft PVC
- Vent holes in the neck reduce the risk of moisture ingress



### Staebel Socks

- Staebel socks provide safety and prevention of falls.
- Soft terry cloth upper. Latex-free.
  - Colour helps care staff identify high fall risk patients
  - Grip on top and bottom
  - Mesh band on large socks comfortably fits swollen ankles and calves
  - Standard size fits up to adult size 13

Available in Green, Orange, Purple, Red or Yellow



### Floor Impact Cushion - Bevelled Edge

- EVA Foam with bevelled edge on three sides
- Helps reduce impact forces by up to 85%
- Low profile edge reduces risk of tripping
- Tested according to ASTM test protocol for shock absorbency
- Wipe clean vinyl cover
- Non-Skid bottom
- Tri-fold design with convenient carrying handle
- Generous size helps provide full coverage of the bed exit area
- Quality product made in USA

### SPECIFICATIONS:

<b>Product Code</b>	92663
<b>Width</b>	610 mm
<b>Length</b>	1220 mm
<b>SWL</b>	200 kg

### SPECIFICATIONS:

<b>Product Code</b>	9900	
<b>Size</b>	Small	
<b>Min Foot Size</b>	215 mm	
<b>Max Foot Size</b>	255 mm	


### SIZES AVAILABLE:

<b>Product Code</b>	9902	9904
<b>Size</b>	Medium	Large
<b>Min Foot Size</b>	255 mm	280 mm
<b>Max Foot Size</b>	280 mm	310 mm

### SPECIFICATIONS:

<b>Product Code</b>	9PP58229
<b>Width</b>	660 mm
<b>Length</b>	1700 mm
<b>Thickness</b>	40 mm

# PATIENT SAFETY



## Cura1 Six Channel Cordless Monitor

- Monitor up to 7 pads or mats simultaneously
- Comes with soft touch removable protective cover, fixing strap and batteries
- Adjustable alarm tone and volume
- Selectable 0, 2sec, 15 min delay setting
- Can be connected to a nurse call system
- Alarm memory and automatic reset
- Pause feature for caregiver to attend patient without setting off the alarm

## Cura1 Cordless Chair Pad

- 1 year life expectancy
- Simple operation and maintenance
- Easy to clean
- For use with cura1 Premium Falls Monitor
- Virtually no false alarms over the lifetime of the pad
- All bed/chair pads made from non-allergenic soft PVC
- Vent holes in the neck reduce the risk of moisture ingress

## Cura1 Cordless Bed Pad

- 1 year life expectancy
- Simple operation and maintenance
- Easy to clean
- For use with cura1 cordless Falls Monitor
- Virtually no false alarms over the lifetime of the pad
- All bed/chair pads made from non-allergenic soft PVC
- Vent holes in the neck reduce the risk of moisture ingress

### SPECIFICATIONS:

<b>Product Code</b>	92700
<b>Width</b>	300 mm
<b>Height</b>	770 mm
<b>Length</b>	1220 mm

### SPECIFICATIONS:

<b>Product Code</b>	92661
<b>Width</b>	255 mm
<b>Length</b>	380 mm
<b>SWL</b>	200 kg

### SPECIFICATIONS:

<b>Product Code</b>	92662
<b>Width</b>	510 mm
<b>Length</b>	760 mm
<b>SWL</b>	200 kg



## TREATMENT TABLES



### Fortress 2 Section Table - Electric Hi-Lo

- 2 section support surface with electrical height adjustment
- Foot control operation (can be hand control if specified at time of order)
- Fitted with 4 castors and central lifting system
- Head section has ergonomically designed removable cushion and facial opening
- Nicely curved cushions with premium padding
- Durable, anti-bacterial, washable and flame retardant upholstery
- Excellent stability due to streamlined base design

#### SPECIFICATIONS:

<b>Product Code</b>	111002
<b>Table Width</b>	700 mm
<b>Head Tilt Angle</b>	-35° - 87°
<b>Min Height</b>	460 mm
<b>Max Height</b>	910 mm
<b>Head Section Length</b>	685 mm
<b>Leg Section Length</b>	1210 mm
<b>Overall Width</b>	700 mm
<b>Overall Length</b>	1930 mm
<b>SWL</b>	225 kg



### Fortress 3 Section Table - Electric Hi-Lo

- 3 section support surface with electrical height adjustment
- Foot control operation (can be hand control if specified at time of order)
- Fitted with 4 castors and central lifting system
- Head section with ergonomically designed removable cushion and facial opening
- Nicely curved cushions with premium padding
- Durable, anti-bacterial, washable and flame retardant upholstery
- Excellent stability due to streamlined base design
- 7 year frame and 3 year electronics warranty

#### SPECIFICATIONS:

<b>Product Code</b>	111003U - WN
<b>Table Width</b>	700 mm
<b>Head Tilt Angle</b>	-35° - 87°
<b>Min Height</b>	430 mm
<b>Max Height</b>	910 mm
<b>Head Section Length</b>	685 mm
<b>Mid Section Length</b>	565 mm
<b>Leg Section Length</b>	680 mm
<b>Overall Width</b>	700 mm
<b>Overall Length</b>	1950 mm
<b>SWL</b>	250 kg



### Mobile Therapy Stool

- This high quality gas stool features a comfortable padded seat
- Five castor base offers greater stability and easy manoeuvrability
- Height adjustable for various treatment positions

#### SPECIFICATIONS:

<b>Product Code</b>	111305
<b>Min Height</b>	490 mm
<b>Max Height</b>	630 mm
<b>Seat Diameter</b>	350 mm
<b>SWL</b>	150 kg













# SEE SOMETHING YOU LIKE? **LET'S HAVE A CHAT!**

**Book your consultancy today with one of our friendly team members.** See below information for ways you can get in contact....



**1800 267 267**



**Showroom Opening Hours**

Monday to Friday | 8:00am to 4:00pm



**[www.activemedicalsupplies.com.au](http://www.activemedicalsupplies.com.au)**



**Cannon Hill**

35 Southgate Ave, Cannon Hill QLD 4170



**[hello@activemedicalsupplies.com.au](mailto:hello@activemedicalsupplies.com.au)**



**Toowoomba**

11 Grange St Toowoomba QLD 4350

**activeMEDICAL**

S U P P L I E S > S O L U T I O N S > S U P P O R T

**P.** 1800 267 267 | **F.** 07 3899 8825

**E.** [hello@activemedicalsupplies.com.au](mailto:hello@activemedicalsupplies.com.au)

**[www.activemedicalsupplies.com.au](http://www.activemedicalsupplies.com.au)**