

Positioning Product Portfolio

Issue 1

Contents

Capability Statement	4
Curera Positioning Neck Horseshoe Pillow	6
Curera Positioning Pillow	6
Curera Positioning Triangle Pillow	7
Back Pillow	7
Positioning Multi Pillow	8
Molly Heel Protection	8
Curera Positioning Curved Pillow	8
Curera Positioning S Pillow	9
Curera Positioning Torso Support Pillow	9
Positioning and Support Wedge	9

CAPABILITY STATEMENT

Active Medical Supplies is a pioneer in medical device and equipment manufacturing. The experienced team at Active Medical understand both professionally and personally, the importance of maintaining independence within the comfort of your own home or the community including educational settings. A lot of the people we help are still independent, however, just need assistive solutions as result of an injury, Hospital visit, Healthcare professional's advice or to help loved ones with their daily life.

With over 30 years of experience helping with solutions in community care, we understand sometimes that daily tasks may become difficult. Often things that we have taken for granted like getting out of bed, lifting yourself from a chair or just even moving around the house, inherit more effort than they have before.

For those who come into our showroom, we can replicate the scenarios you would face at home with our daily living set ups. Our experienced and courteous staff gladly provide all the information and assistance you require.

We work extensively with Physiotherapists, Occupational Therapists and Aged Care Clinicians so that we can find the right form of solution for you. Over the course of 30 years, we have provided solutions to community care with diligence and professionalism. Our customers' assurance comes from the fact that we help you find the most suitable solution.


We take great pride in the assembling and manufacturing of all of our products. We culminate to the highest standards in care, comfort for the patient and project cost efficiency.

Procurement & Distribution

Active Medical has developed a network of manufacturers in Australia and around the world – from Germany and Sweden to the United States of America and consults directly with industry professionals to ensure the products we provide responds to unique patient specifications.

The success of our business is supported by our ability to collaborate with manufacturers that not only provide the highest quality of products and services, but are aligned with our mission and commitment to our social, corporate and environmental responsibilities.

Our efficient distribution systems are underpinned by our strong relationships with manufacturers and two centralised warehouses with a capacity to hold up to 50,000 items. This ensures the fast turn around and on time delivery demanded by the health care industry and the needs of people in care. We ensure that all our products are ethically sourced and supplied and meet environmental standards.



*Disclaimer

Active Medical Supplies does not take responsibility for incorrect application of products. This publication is a guide only to assist in selecting products suitable to your clients. Please use in conjunction with professional advice.

CAPABILITY STATEMENT

Did you know...



We sanitize **all products** that leave our warehouses



We offer full product **demonstrations and video trials**



We offer **contact free delivery**



Our staff are **double vaxxed** with boosters




Our showroom is a **controlled environment** that we can also offer trials

Come in & have a chat!

Book a consultancy today!

Active Medical facilitates a showroom on site to replicate the home care environment. Our showrooms are used to replicate a residential care bedroom setting and allows us to offer unique demonstrations and training tailored to your needs.

We are located in Cannon Hill and Toowoomba with wheelchair accessibility. We showcase our extensive range of electric beds, wheelchairs, mobility, patient handling, safety equipment, bathroom solutions and accessories.



Cannon Hill

35 Southgate Ave, Cannon Hill QLD 4170



Toowoomba

11 Grange St Toowoomba City QLD 4350



Showroom Opening Hours

Monday to Friday | 8:00am to 4:00pm



Showroom Opening Hours

Monday to Friday | 8:00am to 4:00pm

POSITIONING PRODUCTS



Curera Positioning Neck Horseshoe Pillow

If a person makes too few spontaneous movements, there is an increased risk of pressure ulcers. Regular changes of position reduce the risk of pressure ulcers. Position as required in both reclining and seated position.

- Curera Positioning Neck Pillow offers support and pressure relief of head / neck.
- Machine washable



SPECIFICATIONS:

Product Code	32303
Width	380 mm
Length	500 mm



Curera Positioning Pillow

Pillows that provide support in multiple positions. If a person makes too few spontaneous movements, there is an increased risk of pressure ulcers. Regular changes of position reduce the risk of pressure ulcers. Position as required in both reclining and seated position.

Curera® Positioning Pillow can be used as an aid in the facilitation of lung expansion for users with breathing and/or ventilation difficulties. They can also be used to support side lying, Support for neck and arms, Support for users affected by stroke and Positioning for legs and feet.



SPECIFICATIONS:

Product Code	32305
Width	500 mm
Length	600 mm

SIZES AVAILABLE:

Product Code	32304	32306
Width	500 mm	500 mm
Length	250 mm	850 mm

POSITIONING PRODUCTS



Curera Positioning Triangle Pillow

A Triangle Pillow under the knees can stabilize the pelvis to prevent sliding in combination with an elevated backrest.

The Triangle Pillow is also suitable for user's with shortened hip flexors and/or hamstrings. The increased flexion in hips and knees can relieve lower back strain. In some cases, the Triangle Pillow can aid in alleviating pressure from under the heels.



SPECIFICATIONS:

Product Code	32318
Width	280 mm
Length	470 mm
Thickness	280 mm

SIZES AVAILABLE:

Product Code	32319
Width	350 mm
Length	620 mm
Thickness	350 mm



Back Pillow

The Back Pillow is a winged-pillow for back relief. The pillow is placed behind the back, and the wing rests under the user to set the position. It can be used for support in side lying, for beds without side guards. Use a back Support Pillow with a flap placed either under the user or the mattress overlay.



SPECIFICATIONS:

Product Code	32301
Width	500 mm
Length	850 mm

POSITIONING PRODUCTS



Positioning Multi Pillow

Avoid ulcers with knee & heel relief. When bed bound individuals don't move enough, it increases the risk of pressure ulcers. With regular position changes, this risk can be reduced and the comfort level increased.

Multiple uses for the Multi Pillow:

- An aid to prevent "drop-foot"
- Reduce pressure under heels
- Support for the knees

Molly Heel Protection

- Easy to use
- Equalizes pressure
- Faster healing
- Improves circulation
- Makes it easier to freely move the legs
- Prevents crossing the legs after hip surgery
- Soft and comfortable
- Wipe with cleaning agent and /or disinfect
- Machine wash max. 95°C
- Easy to clean
- Liquid resistant
- Long durability

Curera Positioning Curved Pillow

If a person makes too few spontaneous movements, there is an increased risk of pressure ulcers. Regular changes of position reduce the risk of pressure ulcers. Position as required in both reclining and seated position.

- Curera Positioning Neck Pillow offers support and pressure relief of head/neck.
- Machine washable



Example of positioning



SPECIFICATIONS:

Product Code	32310
Width	600 mm
Length	700 mm

SIZES AVAILABLE:

Product Code	32302
Width	450 mm
Length	450 mm

SPECIFICATIONS:

Product Code	32311
Leg Length	780 - 850 mm
Thickness	75 mm

SIZES AVAILABLE:

Product Code	32312	32313
Leg Length	860 - 940 mm	1020 - 1030 mm
Thickness	90 mm	110 mm

SPECIFICATIONS:

Product Code	32321
Width	750 mm
Length	1200 mm
Thickness	30 mm

POSITIONING PRODUCTS



Curera Positioning S Pillow

Pillows that provide support in multiple positions. If a person makes too few spontaneous movements, there is an increased risk of pressure ulcers. Regular changes of position reduce the risk of pressure ulcers. Position as required in both reclining and seated position.

Curera® Positioning Pillow can be used as an aid in the facilitation of lung expansion for users with breathing and/or ventilation difficulties.

They can also be used to support side lying, support for neck and arms, Support for users affected by stroke and Positioning for legs and feet.



SPECIFICATIONS:

Product Code	32322
Width	280 mm
Length	1700 mm



Curera Positioning Torso Support Pillow

The Torso Support Pillow can be placed either behind the user for support when seated on the bedside or in front of the user, preferably in combination with a Curved Pillow behind.

The Trunk Support Pillow can also be used as a support in a half sitting position when in bed or while seated in a chair or wheelchair.



SPECIFICATIONS:

Product Code	32324
Width	650 mm
Length	850 mm
Thickness	300 mm



Positioning and Support Wedge

This foam positioning support wedge is used to place a user in a side position in order to relieve the sacrum, which is an extra sensitive area. For the side position, two support wedges are used together.

The support wedge is easy to use and easy to clean. It is supplied with a zipped, liquid resistant and vapour permeable cover. The cover can be wiped or machine washed.

- Soft comfortable foam, stretch vinyl covered




SPECIFICATIONS:

Product Code	3689
Width	260 mm
Length	700 mm
Thickness	50 - 150 mm

HAVE YOU SEEN OUR OTHER CATALOGUES?

REQUEST A CATALOGUE TODAY



Please email hello@activemedicalsupplies.com.au to request your copy of any of our catalogues!



SEE SOMETHING YOU LIKE? **LET'S HAVE A CHAT!**

Book your consultancy today with one of our friendly team members. See below information for ways you can get in contact....



1800 267 267



Showroom Opening Hours

Monday to Friday | 8:00am to 4:00pm



www.activemedicalsupplies.com.au



Cannon Hill

35 Southgate Ave, Cannon Hill QLD 4170



hello@activemedicalsupplies.com.au



Toowoomba

11 Grange St Toowoomba QLD 4350

activeMEDICAL

S U P P L I E S > S O L U T I O N S > S U P P O R T

P. 1800 267 267 | **F.** 07 3899 8825

E. hello@activemedicalsupplies.com.au

www.activemedicalsupplies.com.au