


activeMEDICAL
SUPPLIES > SOLUTIONS > SUPPORT

JUVO
solutions

Trolleys Portfolio

Issue 1





Contents

01. Instrument & Dressing Trolleys	4
02. Specialist Trolleys	8
03. General Purpose Trolleys	11
04. General Purpose Products	13
05. Accessories	17

01.

Instrument & Dressing Trolleys

Instrument Trolley	5
Standard Dressing Trolley	5
Dressing Trolley with 1 Drawer	6
Dressing Trolley with 2 Drawer	6
Dressing Trolley with 3 Drawer	7
Dressing Trolley with 4 Drawer	7



Instrument Trolley

- 75mm twin braking castors medical design
- Rounded corners —minimise impact damage
- Material grade 304 stainless steel
- Anti-bacterial, Fully welded joints leaving no crevasse to harbour infection
- Soft close drawer runners

Specifications:	
Product Code	8i050
Width	500 mm
Length	500 mm
Height	900 mm

Also Available:			
Product Code	8i070	8i090	8i120
Width	700 mm	900 mm	1100 mm
Length	500 mm	500 mm	500 mm
Height	900 mm	900 mm	900 mm



Standard Dressing Trolley

- Rails on shelves 3 sides
- Material grade 304 stainless steel
- Anti-bacterial, Fully welded joints leaving no crevasse to harbour infection
- Soft close drawer runners

Specifications:	
Product Code	8050
Width	500 mm
Length	500 mm
Height	900 mm

Also Available:			
Product Code	8070	8090	80120
Width	700 mm	900 mm	1100 mm
Length	500 mm	500 mm	500 mm
Height	900 mm	900 mm	900 mm



Dressing Trolley with 1 Drawer

- Rounded corners —minimise impact damage
- One drawer
- Material grade 304 stainless steel
- Anti-bacterial, Fully welded joints leaving no crevasse to harbour infection
- Soft close drawer runners

Specifications:

Product Code	81050
Width	500 mm
Length	500 mm
Height	900 mm

Also Available:

Product Code	81070	81090	810120
Width	700 mm	900 mm	1100 mm
Length	500 mm	500 mm	500 mm
Height	900 mm	900 mm	900 mm



Dressing Trolley with 2 Drawer

- Rounded corners —minimise impact damage
- Two drawers
- Material grade 304 stainless steel
- Anti-bacterial, Fully welded joints leaving no crevasse to harbour infection
- Soft close drawer runners

Specifications:

Product Code	82050
Width	500 mm
Length	500 mm
Height	900 mm

Also Available:

Product Code	82070	82090	82120
Width	700 mm	900 mm	1100 mm
Length	500 mm	500 mm	500 mm
Height	900 mm	900 mm	900 mm



Dressing Trolley with 3 Drawer

- Rounded corners —minimise impact damage
- Three drawers
- Material grade 304 stainless steel
- Anti-bacterial, Fully welded joints leaving no crevasse to harbour infection
- Soft close drawer runners

Specifications:

Product Code	83050
Width	500 mm
Length	500 mm
Height	900 mm

Also Available:

Product Code	83070	83090	83120
Width	700 mm	900 mm	1100 mm
Length	500 mm	500 mm	500 mm
Height	900 mm	900 mm	900 mm



Dressing Trolley with 4 Drawer

- Rounded corners —minimise impact damage
- Four drawers
- Material grade 304 stainless steel
- Anti-bacterial, Fully welded joints leaving no crevasse to harbour infection
- Soft close drawer runners

Specifications:

Product Code	84050
Width	500 mm
Length	500 mm
Height	900 mm

Also Available:

Product Code	84070	84090
Width	700 mm	900 mm
Length	500 mm	500 mm
Height	900 mm	900 mm

02.

Specialist Trolleys

Plaster Trolley	9
Anaesthesia Trolley	9
Resuscitation Trolley	10
Flow Chart Table	10



Plaster Trolley

- Large handle for easy access
- Plaster roll holder
- Fold down shelf with bowl hole. Anti rattle design for silent transport. Bowl included.
- 6 drawers or 3 drawers

Specifications:

Product Code	83063
Height to Top Shelf	900 mm
Depth	500 mm
Width / Width with Handles	600 / 1000 mm
Shelf Fold Down	500 x 400 mm
Bowl Diameter Hole	350 mm

Also Available:

Product Code	83093
Height to Top Shelf	900 mm
Depth	600 mm
Width / Width with Handles	900 / 1700 mm
Shelf Fold Down	600 x 400 mm
Bowl Diameter Hole	350 mm



Anaesthesia Trolley

- Large handle for easy access
- Plaster roll holder
- Fold down shelf with bowl hole. Anti rattle design for silent transport. Bowl included.
- 6 drawers or 3 drawers

Specifications:

Product Code	84064
Height to Top Shelf	900 mm
Depth	500 mm
Width	600 mm



Resuscitation Trolley

- Anti-Bacterial – Fully welded joints leaving no crevasse to harbour infection
- 4 Stainless steel castors 125mm diameter – to assist with manoeuvring – including 2 total locking
- Accessories Rails as standard
- Soft close drawer runners
- 2 x 2 Drawers
- One Shelf
- 3 Side Rails on top shelf
- Material Grade 304 Stainless Steel



Flow Chart Table

- Full width pen holder and bottom of board
- Positive tilt mechanism with safety lock for horizontal position
- 2 clips at top of board
- Gas strut assisted height adjustment
- Adjustable Height

Specifications:

Product Code	84065
Height to Top Shelf	900 mm
Depth	500 mm
Width	600 mm

Specifications:

Product Code	80075
Width	750 mm
Length	660 mm
Height	850 - 1150 mm
Table	750 x 650 mm
Shelves Removable	350 x 520 mm

03.

General Purpose Trolleys

Bassinette Trolley

12

Chart Trolley

12



Bassinette Trolley

- Rounded corners – minimise impact damage
- Angle adjustment to acrylic tub with flat bar support rails
- Clear acrylic tub and mattress included



Chart Trolley

- Removable hanging rails with removable divider — convert to all purpose trolley/cart
- Holds A4 charts
- Capacity 12 or 24 chart

Specifications:

Product Code	80086
Width	420 mm
Length	860 mm
Height	1080 mm
Mattress Base Height	820 mm

Specifications:

Product Code	80039
Width / Width with Handles	390 / 460 mm
Length	500 mm
Height	890 mm

Also Available:

Product Code	80082
Width / Width with Handles	820 / 960 mm
Length	500 mm
Height	890 mm

04.

General Purpose Products

IV Stand	14
Bed Pan & Bottle Holder Rack	14
Oxygen Bottle Holder & Clamp	15
IV Pole & Clamp for Instrument	15
IV Hook Accessories	16



Oxygen Bottle Holder & Clamp

- Heavy duty cast alloy 5 leg base
- 75mm twin braking castors
- Adjustable height
- 23mm diameter pole
- A pull handle as an accessory
- Removable head — 8mm stainless steel wire
- One stand suits 5 head types (see below)
- Comes standard with 4 prong head

Specifications:

Product Code	8302
--------------	------



IV Pole & Clamp Instrument Trolley

- Hold 4 or 6 bottles
- Antibacterial fully welded design

Specifications:

Product Code	8301
--------------	------



IV Stand

- Heavy duty cast alloy 5 leg base
- 75mm twin braking castors
- Adjustable height
- 23mm diameter pole
- A pull handle as an accessory
- Removable head — 8mm stainless steel wire
- One stand suits 5 head types (see below)
- Comes standard with 4 prong head

Specifications:

Product Code	8008
Width	580 mm
Height	1430 – 2430 mm



Bed Pan & Bottle Holder Rack


- Hold 4 or 6 bottles
- Antibacterial fully welded design

Specifications:

Product Code	8009
Width	550 mm
Height	560 mm

Also Available:

Product Code	8010
Width	800 mm
Height	560 mm



4 Prong Head

Specifications:

Product Code 8008B



Outside Hook

Specifications:

Product Code 8008C



Inside Hook

Specifications:

Product Code 8008D



Carol Head

Specifications:

Product Code 8008E



Ramead

Specifications:

Product Code 8008F



Pull Handle

Specifications:

Product Code 8008H

05.

Accessories

- Additional Rails 18
- Trolley Handle 18
- Plaster Roll Holder 18
- Fold Down Shelf 19
- Fold Down Shelf 19
- SS Dia Bowl 19



Additional Rails

- Order with any Juvo trolley. Suits fixing centres of accessories from wide range of suppliers.
- Welded to Side of trolley



Trolley Handle

- Add to any trolley with Additional rails



Plaster Roll Holder

- Add to any trolley with Additional rails
- Removable



Fold Down Shelf

- Add to any trolley with Additional rails
- NO bucket hole



Fold Down Shelf

- Add to any trolley with Additional rails
- With bucket hole



SS Dia Bowl

- 320 mm

Specifications:

Product Code 8005

Specifications:

Product Code 8001

Specifications:

Product Code 8002

Specifications:

Product Code 8003

Specifications:

Product Code 8004

Specifications:

Product Code 8006

activeMEDICAL
SUPPLIES > SOLUTIONS > SUPPORT

P. 1800 267 267 | **F.** 07 3899 8825
E. hello@activemedicalsupplies.com.au

www.activemedicalsupplies.com.au