

activeMEDICAL

Seating Portfolio

Issue 2





Contents

01. General Wheelchairs	5
02. Hospital Focus	10
03. Embark Chairs	15
04. Recliners	19
05. Regency Chairs	30
06. Sertain Chairs	34
07. Bariatric	39
08. Cushion Matrix	43

We are Here to Help

For over 30 years, Active Medical Supplies has delivered safe and innovative solutions for the care of clients in the Aged Care, Hospital, School and Home Care sectors. We strive to provide the highest quality equipment for all healthcare environments which will withstand the test of time.

We partner with experts in the field of design and manufacturing to provide the very best resources to ensure that your loved ones are provided with the correct equipment to receive optimum care.

We have a commitment to lower our carbon footprint and protect our environment. We utilise our own delivery vehicles and use sustainable materials for packaging to reduce environmental damage. We are also proud to provide beds and equipment with world leading, energy-saving electronics.

For your peace of mind our product is manufactured to very high quality and we have a secure supply chain. 90% of our products are sourced from Australia, Europe, USA, Canada, United Kingdom, Spain, India, Sweden, Vietnam and Taiwan.

Procurement & Distribution

Active Medical has developed a network of manufacturers in Australia and around the world – from Germany and Sweden to the United States of America and consults directly with industry professionals to ensure the products we provide responds to unique patient specifications.

The success of our business is supported by our ability to collaborate with manufacturers that not only provide the highest quality of products and services, but are aligned with our mission and commitment to our social,

corporate and environmental responsibilities.

Our efficient distribution systems are underpinned by our strong relationships with manufacturers and two centralised warehouses with a capacity to hold up to 50,000 items. This ensures the fast turn around and on time delivery demanded by the health care industry and the needs of people in care. We ensure that all our products are ethically sourced and supplied and meet environmental standards.



Did you know...



We sanitize **all products** that leave our warehouses



We offer full product **demonstrations and video trials**



We offer **contact free delivery**



Our staff are **double vaxxed** with boosters



Our showroom is a **controlled environment** that we can also offer trials

*Disclaimer

Active Medical Supplies does not take responsibility for incorrect application of products. This publication is a guide only to assist in selecting products suitable to your clients. Please use in conjunction with professional advice. Products and information listed in this catalogue are subject to change. Website and digital information supercedes all print. For the most up to date information, please contact us directly.

01.

General Wheelchairs

Stael Emperor Wheelchair	6
Stael Transit Emperor Wheelchair	6
Stael Thrifty Wheelchair	7
Stael Liberation Ultra - Light Wheelchair	7
Stael Deluxe Alloy Wheelchair	8
Stael Sovereign Wheelchair	9
Wheelchair Accessories	10

NDIS

NDIS



Stael Emperor Wheelchair

The Stael Emperor is a comfortable, lightweight aluminium wheelchair, with a really helpful range of high quality features. With its sleek modern frame and lightweight foldable design, the wheelchair is perfect for travel. All features of this chair have been designed to make life easier for the user.

- Foldable backrest
- Lap belt
- Heel strap
- Height adjustable armrests and footplates
- Removable swing-away footrest
- Available in orange or blue.



Specifications:

Product Code	40411
Seat Width	450 mm
Overall Width	670 mm
Overall Height	950 mm
Product Weight	12 kg
SWL	115 kg

Also Available:

Product Code	40412
Seat Width	500 mm
Overall Width	720 mm
Overall Height	950 mm
Product Weight	15 kg
SWL	135 kg

Stael Transit Emperor Wheelchair

The Stael Emperor is a comfortable, lightweight aluminium wheelchair, with a really helpful range of high quality features. With its sleek modern frame and lightweight foldable design, the wheelchair is perfect for travel. All features of this chair have been designed to make life easier for the user.

- Foldable backrest
- Lap belt
- Heel strap
- Height adjustable armrests and footplates
- Removable swing-away footrest
- Available in orange or blue.



Specifications:

Product Code	40413
Seat Width	450 mm
Overall Width	670 mm
Overall Height	950 mm
Product Weight	12 kg
SWL	115 kg

Also Available:

Product Code	40414
Seat Width	500 mm
Overall Width	720 mm
Overall Height	950 mm
Product Weight	15 kg
SWL	135 kg



Staebel Thrifty Wheelchair

The Thrifty is a sturdy and reliable standard weight wheelchair, with removable armrests and leg rests, built for individuals, nursing homes and non-prescription requirements. This wheelchair meets Australian Standards for Manual Wheelchairs in Product Performance. All of these features combine to make the Thrifty of unbeatable value in standard wheelchairs.

- Swing away, nonremovable armrests
- Swing away, detachable legrests
- 24" rear wheels with urethane tyres (maintenance free)
- Australian Standard AS/NZS 3695.1:2011 approved



Specifications:

Product Code	4246
Seat Width	450 mm
Seat Depth	400 mm
Overall Width	650 mm
Overall Depth	1080 mm
Overall Height	880 mm
Fold Width	300 mm
Fold Length	790 mm
SWL	115 kg



Staebel Liberation Utra Lightweight Wheelchair

The Staebel Liberation Lightweight Wheelchair is a great travel companion. It offers attendant and parking brakes, ultra lightweight frame, folds down in seconds, and has foldable and removable legrests.

- Folding backrest
- Safety Brake and Assistant Brake
- Folding Arm Rest and Foot Rest
- Includes carry case



Specifications:

Product Code	40400
Seat Width	450 mm
Overall Width	540 mm
Overall Height	920 mm
Weight	7 kg
SWL	100 kg



NDIS



Staebel Deluxe Alloy Wheelchair

The Deluxe Ultralight Alloy Wheelchairs have been designed to accurately meet the needs of our clients. They are an ultralight weight wheelchair with the benefit of many additional features. The improvement compared to previous designs (ranging from 25 - 30 kg gross weight) is considerable. Much input from clients and nurses has assisted in the development of these superior quality wheelchairs.

- Flip back armrests
- Swingaway, detachable legrests

- Fold down backrest
- Quick release rear wheels
- 24" rear wheels with urethane tyres
- Aluminium frame
- Padded nylon upholstery
- Complies with AS3695:2008



Specifications:		Sizes Available:	
Product Code	4629	Product Code	45129 45629
Seat Width	460 mm	Seat Width	510 mm 560 mm
Seat Depth	400 mm	Seat Depth	400 mm
Seat Height	480 mm	Seat Height	490 mm
Backrest Height	430 mm	Backrest Height	430 mm
Armrest Height	230 mm	Armrest Height	210 - 280 mm 210 - 280 mm
Overall Width	650 mm	Overall Width	700 mm 750 mm
Overall Depth	1060 mm	Overall Depth	1000 mm 1000 mm
Overall Height	960 mm	Overall Height	960 mm 960 mm
Fold Width	320 mm	Fold Width	300 mm 300 mm
Net/Gross Weight	11 / 14 kg	Net/Gross Weight	14.5 / 17 kg 15 / 17.5 kg
SWL	125 kg	SWL	150 kg 150 kg

NDIS
MASS



Staebel Sovereign Wheelchair

The Sovereign wheelchair can be adjusted to suit almost any user. With much consideration in the design from the perspective of OTs, this chair is suited to the vast range of needs required by our clients. Rather than numerous script forms for various users, this chair meets the needs of almost any patient. All wheelchair parts suit any adjustments made - no additional parts required.

- Adjustable & removable leg rest's
- Flip back & removable armrests
- Adjustable armrest height & length
- Seat depth, height & angle adjustment
- Backrest height, tension & angle adjustment
- N.B. Seat Cushion is not included with Wheelchair.



Specifications:

Product Code	4340C - 43	4340C - 46	4340C - 49	4352C - 52	4352C - 55	4352C - 58
Seat Width	430 mm	460 mm	490 mm	520 mm	550 mm	580 mm
Anterior Seat Height	420 mm	445 mm	450 mm	470 mm	475 mm	500 mm
Posterior Seat Height	410 mm	430 mm	450 mm	470 mm	490 mm	510 mm
Seat Depth	400 mm	425 mm	450 mm	475 mm	475 mm	500 mm
Backrest Height (2cm increments)	400 - 460 mm	400 - 460 mm	400 - 460 mm	400 - 460 mm	400 - 460 mm	400 - 460 mm
Backrest Angle (5° increments)	69° - 94°	69° - 94°	69° - 94°	69° - 94°	69° - 94°	69° - 94°
Lower Leg Length (2cm increments)	400 - 480 mm	400 - 480 mm	400 - 480 mm	400 - 480 mm	400 - 480 mm	400 - 480 mm
Armrest Height (2cm increments)	180 - 260 mm	180 - 260 mm	180 - 260 mm	180 - 260 mm	180 - 260 mm	180 - 260 mm
Armrest Length (4cm increments)	330 - 450 mm	330 - 450 mm	330 - 450 mm	330 - 450 mm	330 - 450 mm	330 - 450 mm
Footplate Width	105 - 155 mm	105 - 155 mm	105 - 155 mm	105 - 155 mm	105 - 155 mm	105 - 155 mm
Net/Gross Weight	17 kg	17 kg	17 kg	21 kg	21 kg	
SWL	125 kg	125 kg	125 kg	170 kg	170 kg	170 kg



Stael Anti-Tip Bar

The anti-tip bars are fitted to the wheelchair to prevent the wheelchair tipping backwards. It is particularly beneficial for amputees or active users.

- Attaches to wheelchair tilt bar
- Includes small plastic rollers

Product codes may vary depending on wheelchair. Please call 1800 267 267 for the correct code.



Stael Swivel Arm Trough

The swivel arm trough is of benefit for the wheelchair user who has upper limb weakness or paralysis. It attaches to the wheelchair armrest via a swivelling clamp.

- One piece moulded polyurethane trough
- Supplied with hardware
- Can be mounted left or right side



Stael Elevating Legrest

The elevating legrest gives the wheelchair user the ability to support their leg in the elevated position. It attaches to the wheelchair in place of the existing legrest.

- Adjustable angle and length to suit the individual user
- Configured left or right

Product codes may vary depending on wheelchair. Please call 1800 267 267 for the correct code.

Specifications:

Product Code 4982



Stael Calf Strap

This padded calf strap provides comfortable support and keeps the feet safely on the footplates. The Velcro straps simply and easily secure to the legrests of the wheelchair.

- Generous size (suitable for 430-490mm chair width)
- Calf section measures 460mm x 89mm

Specifications:

Product Code 4986



Stael Oxygen Bottle Carrier

The oxygen bottle carrier is designed to be attached to the wheelchair to enable the oxygen dependent user to be transported in the wheelchair.

- Suitable for type C oxygen cylinders
- Removes easily for storage and transportation
- Slides onto existing wheelchair tilt bar
- Suitable for Stael® wheelchairs and some other brands

Specifications:

Product Code 4980 - L

Product Code 4980 - R



Stael I.V. Pole

The IV pole is attached to the wheelchair and allows the intravenous drip dependent client to be transported in a wheelchair.

- Easy height adjustment with triwheel
- Slides onto existing wheelchair tilt bar

Specifications:

Product Code 4979

Specifications:

Product Code 4983

Specifications:

Product Code 4984



Seat Belt Metal Buckle - Screw On

This seat belt is a simple device for safety and positioning of the wheel chair user. It features a push button buckle and is adjustable in length. It has metal eyelets to enable attachment via existing seat upholstery screws on the users wheelchair.

- Adjustable in length
- Simply and easily fitted to most wheelchairs
- 62-154cm

Staebel Elevating Legrest

The elevating legrest gives the wheelchair user the ability to support their leg in the elevated position. It attaches to the wheelchair in place of the existing legrest.

- Adjustable angle and length to suit the individual user
- Configured left or right - specify when ordering

*Only suitable for Deluxe Alloy range of wheelchairs

Contact us for the codes for other models

Staebel Stump Support

The stump support is used on a wheelchair in place of the regular legrest, and is useful for amputees with a below knee amputation.

- Padded stump support
- Adjustable pad height
- Left or right side configured - specify when ordering

*Only suitable for Deluxe Alloy range of wheelchairs

Contact us for the codes for other models

Specifications:

Product Code 4974

Specifications:


Left Leg 4980 - L

Right Leg 4980 - R

Specifications:

Left Leg 4981 - L

Right Leg 4981 - R



Staebel Anti - Fold/ Theft Device

The anti-fold/ anti-theft device is used on a wheelchair to prevent the frame being folded, and to minimise the risk of theft. It is easily fitted to the frame between the rear wheels.

- Adjustable in width

Specifications:

Product Code 4978

Length 400 - 490 mm

Sizes Available:

Product Code 4978L 4978XL

Size Large Extra Large

Length 510 - 560 mm 610 - 710 mm

Staebel Universal Tray

Universal Tray fit to most chairs and wheelchairs.

- Cutout allows closer fitting for the convenience of the user
- Raised lip edge helps prevent items falling of tray
- Moulded cup holder recess
- Fits securely by adjustable velcro straps
- Lightweight and durable
- Suits most wheelchairs

Specifications:

Product Code 4972

Width 600 mm

Staebel Handy Bar

Transfer Bar wheelchair/car transfer bar.

- Firm secure handle
- Simple to use
- Easy to install and remove
- Can be used to break glass in emergency
- Inbuilt seat belt cutter
- Fits most vehicles - no modification needed
- Recommended by Doctors and Health Professionals

Specifications:

Product Code 4998



Stael Sovereign Arm Trough One Piece

The arm trough is of benefit for the wheelchair/ shower chair user who has upper limb weakness or paralysis. It attaches to the chair armrest via a clamp.

- One piece moulded polyurethane trough
- Supplied with hardware
- Can be mounted left or right side

Specifications:

Product Code	4985
---------------------	------



Stael Sovereign Stump Support

The stump support is used on a wheelchair in place of the regular legrest, and is useful for amputees with a below knee amputation.

- Adjustable pad height
- Left or right side configured - specify when ordering

Specifications:

Product Code	43417 - L
Product Code	43417 - R



Stael Sovereign One Arm Drive Kit

One arm drive system with double push rim (Must specify left or right at time of order)

Specifications:

Product Code	43425
---------------------	-------



Stael Sovereign Seat Belt

This seat belt is a simple device for safety and positioning. It has a push button buckle and is adjustable in length. It has metal eyelets to enable attachment via existing seat upholstery screws on the wheelchair.

- Adjustable in length
- Simply and easily fitted to most wheelchairs

Specifications:

Product Code	4974
---------------------	------



Stael Sovereign Elevating Legrest

This seat belt is a simple device for safety and positioning of the wheel chair user. It features a push button buckle and has metal eyelets to enable attachment via existing seat upholstery screws on the wheelchair.

- Adjustable in length

Specifications:

Product Code	43416 - L
Product Code	43416 - R



Stael Sovereign Transit Wheel Kit

Transit wheel kit, to suit Sovereign Wheelchair
Not compatible with 52, 55, 58 cm

Specifications:

Product Code	43426
---------------------	-------



Stael Sovereign Drum Brake Kit

- Drum Brake Conversion Kit for the Sovereign Wheelchair
- Includes Brakes, Cables & Wheels

Specifications:

Product Code	43410
---------------------	-------



Stael Sovereign Pelvic Strap

- Plastic side release buckle belt, 2-piece design
- Loop tie attachment around wheelchair/shower chair frame
- Length adjustable 80-155 cm

Specifications:

Product Code	19636
---------------------	-------



Stael Sovereign Head Support

- Head Support to suit Sovereign Wheelchair
- Height Adjustable

Specifications:

Product Code	49608
---------------------	-------

02.

Hospital Focus

Staebel Patient Transport Chair	14
Staebel Gynecology Chair	14
Staebel Blood Donor Chair	15
Pura Multifunctional Chair with Arms	15
Hospital Ward Wheelchair	16
EasyGlide Transfer Board	16
Turn Table	17
EasyGlide Oval Board	17
Staebel Mobile Therapy Stool	18
Staebel Saddle Treatment Stool	18
Staebel Stool with Back - Saddle Design	18

NDIS

NDIS



Stael Patient Transport Chair

Stael Patient Transport Chair is a handy, easy to use patient transport chair with flip-up armrests, generous 52cm side open seat and a retractable foot rest.

- 200kg SWL
- PVC document holder
- Adjustable braking system
- Silicone sealed seat for infection control
- HPL footplate
- UV protected Powdercoating
- Dual front castors for easy maneuvering.

Stael Gynecology Chair

- Three section chair top
- Head inclination 12° (head length 30 cm)
- Back part inclination electrically adjustable from - 10° to +65° (length 60 cm)
- Seat inclination electrically adjustable from 0° to +30°
- Footrest manual adjustable
- Arm supports manual adjustable
- Height adjustment electrically adjustable - column lift from 56 to 94 cm
- Stable base on the leveling feets
- Seamless padding from artificial leather for medical purposes (5 cm) wear-resistant and disinfectant tolerant material
- Range of upholstery colours available



Specifications:

Product Code	4750
Width	660 mm
Length	870 mm
Height	1090 mm
Seat Height	620 mm
Seat Width	520 mm
Seat Depth	430 mm
Weight	28.5 kg
SWL	200 kg



Specifications:

Product Code	71160
Width	700 mm
Length	1380 mm
Weight	74.5 kg
SWL	150 kg



Staebel Blood Donor Chair

An ergonomically designed chair, which enables a relaxed and comfortable sitting or lying during the procedure. The armrests are movable into all directions. Armrests suitable for blood donation procedures. All parts have quality upholstery with rounded corners with antibacterial artificial leather-skai, special for medical purposes.

- Patient maximum weight: 170 kg;
- Electrically adjustable seat section height
- Seamless padding 5 cm.
- Swivel double wheel casters with diameter of 75mm with brake
- Metal frame powder-coated with antibacterial paint
- The base is made of powder coated metal elements
- Backrest inclination adjustable via gas spring from -10° to $+80^{\circ}$
- Seat section fixed $2,5^{\circ}$ seat section inclination
- Inclination adjustable via gas spring from -90° to $+10^{\circ}$ footrest
- Upholstery colour by agreement. Contact our team for more colour options.



Specifications:

Product Code	5201
Seat Width	600 mm
Seat Length	550 mm
Overall Width	600 mm
Overall Length	1000 mm
Overall Length (Lying Position)	1850 mm
Seat Height Adjustment	550 - 750 mm
Backrest (Width x Length)	600 x 850 mm
Footrest (Width x Length)	450 x 400 mm
SWL	170 kg

Pura Multifunctional Chair with Arms

The Pura was designed to improve the patient flow with advanced innovations during long-term procedures such as day surgery, diagnostic procedures, dialysis, oncology treatment or blood donation.

- Three section chair top
- Head inclination 12° (head length 30 cm)
- Back part inclination electrically adjustable from -10° to $+65^{\circ}$ (length 60 cm)
- Seat inclination electrically adjustable from 0° to $+30^{\circ}$
- Footrest manual adjustable
- Arm supports manual adjustable
- Height adjustment electrically adjustable - column lift from 56 to 94 cm
- Stable base on the leveling feets
- Seamless padding from artificial leather for medical purposes (5 cm) wear-resistant and disinfectant tolerant material
- Range of upholstery colours available



Specifications:

Product Code	5400
Seat Width	590 mm
Seat Depth	500 mm
Backrest Height	950 mm
Legrest Length	650 mm
Weight	104 kg
SWL	200 kg

Also Available:

Pura Multifunctional Chair with Soft Drop Siderail	
Product Code	5405

NDIS

NDIS



Hospital Ward Wheelchair

The Hospital Ward Alloy Wheelchair has been designed to accurately meet the needs of our patients. Much input from hospital staff has assisted in the development of these superior quality wheelchairs.

- Swingaway, detachable legrests (if non-removable legrests are required this can be done at time of order)
- Height adjustable, flip back armrests - this assists in side transfers of the client
- IV Pole and oxygen bottle carrier
- Anti tip bars
- 12" rear wheels with urethane tyres - no flat tyres
- Aluminium construction - less force to push
- Padded nylon upholstery - can be wiped over to clean
- Complies with AS3695:2008

EasyGlide Transfer Board

The flexible three-part sliding board is easy to place in front of the wheels of a wheelchair to protect against soiling during a sliding transfer. The smooth top surface makes it easy to move and provides for easy transfers while the double anti-slip tape on the underside prevents the board from sliding during transfer. It is especially suitable for levelling out height differences or for bridging short distances. EasyGlide is easy to transport and convenient to always have close at hands.

EasyGlide is made of strong polyethylene and is available in different sizes. The longer versions are especially suitable for transfers from wheelchair to car.

- Anti-slip protection on the underside provides for safety and security.
- The long sides can be folded to protect clothing from the wheelchair's wheels.



Specifications:

Product Code	45128WD
Seat Width	510 mm
Seat Depth	400 mm
Seat Height	490 mm
Backrest Height	430 mm
Armrest Height	210 - 280 mm
Overall Width	700 mm
Overall Depth	1000 mm
Overall Height	960 mm
Fold Width	350 mm
Net/Gross Weight	15 / 17 kg
SWL	150 kg



Specifications:

Product Code	215020
Width	330 mm
Length	600 mm
Thickness	5 mm
SWL	135 kg

Sizes Available:

Product Code	215050
Width	330 mm
Length	750 mm
Thickness	6 mm



Turn Table

Many people have enough strength to stand on their legs but not the ability to move their feet. In these circumstances, TurnTable, which is placed on the floor under the feet, is an outstanding, easy-to-use aid for the more skilled carer and for transferring a user between chair and wheelchair, wheelchair and toilet or from wheelchair to bed. The carer can guide and slow the transfer by keeping half a foot on the TurnTable.

For maximum safety in transfers, without the risk of twisting ankles and knees, we recommend using EasyBelt or FlexiBelt, alternatively FlexiMove, together with TurnTable.

TurnTable is available in three sizes with or without handle. TurnTable is made of durable polyethylene with a rubbery anti-slip surface.

- The surface provides a good grip against both feet and floor.
- The two plates rotate smoothly over each other.
- A practical handle facilitates handling.
- Convenient size and low weight.
- Easy to clean with an alcohol or soap solution.

EasyGlide Oval Board

Another popular product from the RoMedic range the EasyGlide black oval boards can be used to facilitate placing of leg supports during application of lifting slings with divided leg supports in a seated position. Due to their unique design and ability to flex these also eliminate the friction under the user and make it very easy to position the user further back in the chair.



Specifications:

Product Code	215040
Width	850 mm
Length	390 mm
Thickness	3 mm
SWL	135 kg

Also Available:

Product Code	215045
Width	450 mm
Length	200 mm

Specifications:

Product Code	215085
SWL	380 mm

Also Available:

Product Code	215084
Length	450 mm

NDIS

NDIS

NDIS



Stael Mobile Therapy Stool

- This high quality gas stool features a comfortable padded seat
- Five castor base offers greater stability and easy manoeuvrability
- Height adjustable for various positions



Stael Saddle Stool

This ergonomically designed saddle stool helps to reduce the discomfort of sitting for long periods and allows the user free range of upper body movement.

- High quality gas stool features a deluxe comfortable padded seat
- Five castor base offers greater stability and easy manoeuvrability
- Height adjustable for various treatment positions
- Anti-bacterial upholstery
- Foam thickness: 6cm



Stael Stool with Back - Saddle Design

This ergonomically designed saddle stool helps to reduce the discomfort of sitting for long periods and allows the user free range of upper body movement.

- High quality gas stool features a deluxe comfortable padded seat
- Anti-bacterial upholstery
- Five castor base offers greater stability and easy manoeuvrability
- Height adjustable for various positions



Specifications:

Product Code	111305
Min Height	490 mm
Max Height	630 mm
Seat Diameter	350 mm
SWL	150 kg



Specifications:

Product Code	111311
Min Height	560 mm
Max Height	630 mm
Cushion Diameter	350 mm
SWL	160 kg



Specifications:

Product Code	111313
Min Height	550 mm
Max Height	760 mm
SWL	150 kg

03.

Embark Chairs

Utilite Day Chair	20
Embark Deluxe Day Chair	20
Embark Superior Day Chair	21
Blood Collection Chair	21
Embark Leg Rest	22
Chair / Bed Raisers Biscuit	22
Chair / Bed Raiser	22
Footrest Stool - Non - Slip	22

NDIS



NDIS



Utilite Day Chair

The Utilite is a multi purpose adjustable seat height and angle chair. Ideal for use in an area where adjustable seating is required, for example, lounge, dining room, bedroom, activity room, outdoor area etc. This low back chair provides comfort, stability and practicality at an affordable price.

- Shallow curved low backrest
- Padded armrests
- Comfortable, durable upholstered seat, back and armrests
- Powdercoated steel frame

Embark Deluxe Day Chair

An exceptional rehabilitation chair, the Embark Deluxe is available in sturdy quality vinyls for ease of cleaning. The adjustable seat height, seat depth and angle ensures Embark Deluxe can adapt to the widest range of people, making this chair ideal for hospital, home and nursing home. The padded seat and back rest have been designed to enhance the posture, whilst providing maximum comfort.

- High backrest incorporating lumbar and lateral support
- Padded armrests
- Comfortable vinyl upholstered seat and back
- Adjustable seat height and depth
- Powdercoated steel frame
- Standard with tilt and move wheels for ease of manoeuvrability



Specifications:

Product Code	5910
Seat Width	460 mm
Seat Depth	460 mm
Seat Height	485 - 640 mm
Backrest Height	455 mm
Arm Height	230 mm
Overall Width	610 mm
Overall Depth	550 mm
Overall Height	900 - 1075 mm
SWL	120 kg

Also Available:

Embark Day Chair Legrest

Product Code	5958
---------------------	------



Specifications:

Product Code	5930
Seat Width	500 mm
Seat Depth	450 - 530 mm
Seat Height	450 - 600 mm
Backrest Height	790 mm
Arm Height	230 mm
Overall Width	700 mm
Overall Depth	730 - 800 mm
Overall Height	1200 - 1350 mm
SWL	200 kg

Also Available:

Embark Day Chair Legrest

Product Code	5958
---------------------	------

NDIS



NDIS



Embank Superior Day Chair

The Embark Superior Day Chair was designed and developed by a team of clinicians. It has a strong frame and a high weight capacity of 200 kg. The chair has low resilient foam and 2-way stretch upholstery on seating surfaces. For increased longevity in demanding environments it has 6cm tilt and move wheels and other small improvements.

There have been instances in the past, whereby patients, staff and visitors have sustained cuts and injuries from square tube frames. As a result, the Embark was not only designed for comfort, but to avoid these injuries.

- Low back option also available.
- Comfortable upholstered seat and back
- Padded armrests
- Adjustable seat height and depth
- Powder-coated steel frame
- High backrest incorporating lumbar and lateral support

Blood Collection Chair

The Blood Collection chair has been designed in conjunction with leading pathology companies to ensure that sample collection is safe for both staff and patient, whilst affording the patient a reasonable degree of comfort. The sloping armrests are positioned to present the arm at the optimum angle and the step is fitted with a non slip surface.

- Elevated Seat height for ease of access
- Wide padded armrests
- Non-slip step
- Reinforced armrest support
- Range of hospital grade vinyl coverings
- A range of colours available, contact our team to discuss.



Scan the QR code to watch the **video demonstration!**

Specifications:

Product Code	5933
Seat Width	500 mm
Seat Depth	450 - 530 mm
Seat Height	450 - 600 mm
Backrest Height	790 mm
Arm Height	230 mm
Overall Width	700 mm
Overall Depth	730 - 800 mm
Overall Height	1200 - 1350 mm
SWL	200 kg



Specifications:

Product Code	5202G
Seat Width	550 mm
Width (With Armsrests)	920 mm
Length	650 mm
Length (With Footrests)	840 mm
Seat Height	640 mm
SWL	130 kg



NDIS

Embark Legrest

This vinyl covered legrest provides comfort for resting legs, with its generously padded top set on a 10 degree angle. Height and angle adjustable using simple push button system. Fitted with non-slip rubber feet.

- Sturdy powdercoated steel frame with vinyl upholstery
- Complements the Embark chair range



Specifications:

Product Code	5958
Height	435 - 585 mm
Width	550 mm
Length	430 mm
SWL	200 kg

DVA

NDIS



Chair / Bed Raiser

Lightweight but extremely strong with a broad base and a slip resistant pad in a recessed top for a chair or bed leg. Can be stacked for greater height increase and storage. Supports up to 250kg.

- Can be stacked/combined to achieve greater height
- Round cavity to support chair leg or bed castor, up to 65mm in diameter



Specifications:

Product Code	5990 / 5991
Hole Diameter	65 mm
Height (5990)	40 mm
Height (5991)	100 mm
SWL	250 kg



NDIS

Chair / Bed Raisers Biscuit

This plastic disc is used with the plastic bed raiser, when the bed leg is too large to fit into the raiser. The biscuit is attached with a screw to the base of the bed leg. The biscuit then can then be slotted into the centre top indent of a plastic 40mm or 100mm height raiser.

Specifications:

Product Code	5992
--------------	------

NDIS



Footrest Stool - Non - Slip

Standard footrest provides support for the feet and assist in the correction of hip angle. The rubber non-slip surface provides added grip.

- Chrome steel frame

Specifications:

Product Code	5950
Width	360 mm
Length	420 mm
Height	150 mm

Also Available:

Product Code	5951
Width	360 mm
Length	420 mm
Height	225 mm

04.

Recliners

Ashley Lift & Recline Chair - One Motor	24
Ashley Lift & Recline Chair - Two Motor	24
Ashley Mobile Recliner Luxor Optima	25
Ashley Mobile Recliner - Standard	25
Ashley Luxor Optima Lift & Recline Chair - One Motor	26
Ashley Luxor Optima Lift & Recline Chair - Two Motor	26
Comflex Windsor Chair	27
Shoalhaven Aircell Recliner	27
Verity Mobile Recliner	28
Carson Recliner Chair	28
Oasis Mobile Recliner with Tray	29
Regency Tilt & Recline	29
Regency Result3 Electric Recliner	30
Regency Bahama Comfort Chair	30
Cocoon G1 Chair	31
Cocoon G2 Chair	31
Configura Chair	32

NDIS

NDIS



Ashley Lift & Recline Chair - One Motor

The Ashley range of chairs have long been the favourite choice with their comfortable seating and range of options. Economically priced and upholstered in standard fabric, the Ashley features a waterfall style backrest, sculptured armrests, side pockets and padded leg chaise. The Ashley chairs can also be ordered in a variety of sizes to suit the individual. This quality powered rise and recline chair is both comfortable and durable, and provides the user assistance with hassle free push button operation.

This model features a single motor allowing the occupant to recline and be raised to the standing position.

Available in either a fabric or vinyl upholstery in a pattern or plain. Contact our team to discuss! Vinyl available in Austex Studio Encore colour range.

- Allow +/- 20mm tolerances



Ashley Lift & Recline Chair - Two Motor

The Ashley range of chairs have long been the favourite choice with their comfortable seating and range of options. Economically priced, the Ashley features a waterfall style backrest, sculptured armrests, side pockets and padded leg chaise. Australian made, The Ashley chairs can also be ordered in a variety of sizes to suit the individual. This quality powered rise and recline chair is both comfortable and durable, and provides the user assistance with hassle free push button operation. The two motor action enables independent control of backrest and leg rest operation for ultimate flexibility of chair position.

Available in either a fabric or vinyl upholstery in a pattern or plain. Contact our team to discuss! Vinyl available in Austex Studio Encore colour range.

- Allow +/- 20mm tolerances



Specifications:

Product Code	5524	5524MI	5524MX
Size	Medium	Mini	Maxi
Seat Width	480 mm	480 mm	530 mm
Seat Depth	500 mm	490 mm	520 mm
Seat Height	500 mm	470 mm	530 mm
Backrest Height	640 mm	600 mm	680 mm
Overall Width	780 mm	780 mm	830 mm
Overall Height	1070 mm	1020 mm	1170 mm
SWL	120 kg	120 kg	120 kg



Specifications:

Product Code	5525	5525MI	5525MX
Size	Medium	Mini	Maxi
Seat Width	480 mm	480 mm	530 mm
Seat Depth	500 mm	490 mm	520 mm
Seat Height	500 mm	470 mm	530 mm
Backrest Height	640 mm	600 mm	680 mm
Overall Width	780 mm	780 mm	830 mm
Overall Height	1070 mm	1020 mm	1170 mm
SWL	120 kg	120 kg	120 kg

NDIS



NDIS



Ashley Mobile Recliner Luxor Optima

The Ashley Luxor range of chairs have been specially developed for those with pressure care needs or those needing increased comfort when seated. This is achieved with layered pressure reducing foam and stretch PU upholstered on the seating surfaces. All the favourite features of the standard Ashley are still included - a waterfall style backrest, sculptured armrests, side pockets and padded leg chaise. Cascading roll back with isolated head rest and lumbar roll to support both these areas

- Upholstered leg rest with calf support
- Manual operated recliner action with lever
- Castors and push handle for mobility
- No sag springs in seat and lumbar for long lasting support
- Quality Australian made
- 12 month warranty on fabric, 24 month warranty on frame
- Adjustable head rest cushion
- 150 mm castors, 3 x braking, 1 x directional lock
- Allow +/- 20mm tolerances

Ashley Mobile Recliner - Standard

Economically priced, the Ashley features a waterfall style backrest, sculptured armrests, side pockets and padded leg chaise. Australian made. Cascading roll back with isolated head rest and lumbar roll to support both these areas.

Available in either a fabric or vinyl upholstery in a pattern or plain. Contact our team to discuss! Vinyl available in Austex Studio Encore colour range.

- Allow +/- 20mm tolerances



Specifications:

Product Code	5519MI	5519	5519MX
Size	Mini	Medium	Maxi
Seat Width	480 mm	480 mm	530 mm
Seat Depth	470 mm	500 mm	520 mm
Seat Height	460 mm	480 mm	500 mm
Backrest Height	600 mm	640 mm	680 mm
Overall Width	780 mm	780 mm	830 mm
Overall Height	1020 mm	1050 mm	1150 mm
SWL	120 kg	120 kg	120 kg

Specifications:

Product Code	5521
Size	Medium
Seat Width	480 mm
Seat Depth	500 mm
Seat Height	480 mm
Backrest Height	640 mm
Overall Width	780 mm
Overall Height	1050 mm
SWL	120 kg

NDIS



NDIS



Ashley Luxor Optima Lift & Recline Chair - One Motor

The Ashley Luxor Optima range of chairs have been specially developed for those with pressure care needs or those needing increased comfort when seated. This is achieved with layered pressure reducing foam and a specialised 2-way stretch upholstery on the seating surfaces. All the favourite features of the standard Ashley are still included – a waterfall style backrest, sculptured armrests, side pockets and padded leg chaise. This model which features a single motor reclines to enable the occupant to rest and will also raise to the standing position.

- Upholstered leg rest with calf support
- One motor wall saver action – chair can be positioned close to within 100mm from wall and still recline
- Backrest reclines to half lying down position
- Chair raises to assist standing transfer
- No sag springs in seat and lumbar for long lasting support
- Tasmanian Oak hardwood frame for durability
- High quality Dunlop Enduro foams used
- Battery back up capability
- Small castors enable chair to be wheeled (not patient transport)
- Quality Australian made product
- 2 year electronic warranty, 5 year frame warranty
- Allow +/- 20mm tolerances



Specifications:

Product Code	5522MI	5522	5522MX
Size	Mini	Medium	Maxi
Seat Width	480 mm	480 mm	530 mm
Seat Depth	490 mm	500 mm	520 mm
Seat Height	470 mm	500 mm	530 mm
Backrest Height	600 mm	640 mm	680 mm
Overall Width	780 mm	780 mm	830 mm
Overall Height	1020 mm	1070 mm	1170 mm
SWL	120 kg	120 kg	120 kg

Ashley Luxor Optima Lift & Recline Chair - Two Motor

The Ashley Luxor Optima range of chairs have been specially developed for those with pressure care needs or those needing increased comfort when seated. This is achieved with layered pressure reducing foam and stretch PU upholstered on the seating surfaces. All the favourite features of the standard Ashley are still included – a waterfall style backrest, sculptured armrests, side pockets and padded leg chaise. The Optima chairs can be ordered in a variety of sizes to suit the individual.

The two motor action enables independent control of backrest and leg rest operation for ultimate flexibility of chair position.

- Upholstered leg rest with calf support
- Backrest reclines to flat position
- Chair raises to assist standing transfer
- No sag springs in seat and lumbar for long lasting support
- High quality Dunlop Enduro foams used
- Battery back up capability
- Small castors enable chair to be wheeled (not patient transport)
- Quality Australian made product
- Battery back up system
- Allow +/- 20mm tolerances



Specifications:

Product Code	5532MI	5532	5532MX
Size	Mini	Medium	Maxi
Seat Width	480 mm	480 mm	530 mm
Seat Depth	490 mm	500 mm	520 mm
Seat Height	470 mm	500 mm	530 mm
Backrest Height	600 mm	640 mm	680 mm
Overall Width	780 mm	780 mm	830 mm
Overall Height	1020 mm	1070 mm	1170 mm
SWL	120 kg	120 kg	120 kg



Scan the QR code to watch the video demonstration!

NDIS



Comflex Windsor Chair

This comfortable recliner is packed with features for both the patient and carer. The framework has been specifically engineered around stability for anti-tipping and it has a one piece backrest to keep cleaning easy. Choose to add a mobile base and slide-out footplate to increase functionality.

- Solid Tasmanian Oak frame
- Available as electric chair or manual (flush trigger pull operation)
- Unique design to limit tipping
- One piece backrest for easy clean & infection control
- Padded leg-rest chaise
- Available in a huge range of Healthcare grade upholstery colours
- Australian made



Specifications:

Product Code	5251M / 5251E
Seat Width	560 mm
Seat Depth	560 mm
Seat Height	550 mm
Backrest Length	520 mm
Overall Width	820 mm
Overall Depth	860 mm
Overall Height	1070 mm
SWL	120 kg

NDIS



Shoalhaven Aircell Recliner

Following on the success of the Shoalhaven comes the Shoalhaven Aircell. This has been engineered to ensure maximum comfort for those confined to a chair, whilst allowing the carer to operate the chair with maximum safety and care. Soft, stretch fabric, quality pressure sensitive foam and a shaped backrest giving lumbar and lateral support provide ultimate comfort.

- Non maintenance, air cell cushioning combined with slow release memory foam, for excellent pressure care.
- Adjustable / removable wings
- Compact postural design with tilt in space recline action, and high legrest elevation for reduced fluid retention
- Resilient, multiple stretch, water proof fabric, offering low skin shear, anti bacterial and non allergenic properties
- Designed for use with standing lifters & transfer hoists
- Low seat height for fallout safety, and easier stand assist transfers
- Hospital quality 100mm locking castors
- Full length padded armrest with bumpers
- Head rest cushion included
- Manufactured in Australia



Specifications:

Product Code	5230
Seat Width	520 mm
Seat Depth	500 mm
Seat Height	500 mm
Backrest Height	750 mm
Arm Height	220 mm
Legrest Length	400 mm
Overall Width	780 mm
Overall Depth	1000 mm
Overall Height	1170 mm
SWL	130 kg

NDIS



Verity Mobile Recliner

The Verity Recliner offers healthcare-grade durability and exudes an inviting, homelike aesthetic that brings residential warmth to clinical environments.

- 4 locking castors
- 3 position recline with easy to operate mechanism
- Ergonomically-positioned push handles for patient transport
- Stitching method reduces seams for increased infection resistance
- CAL TB133 compliance
- Moisture Barrier
- Optional arm caps available in polyurethane
- Contrasting fabrics are obtainable

NDIS



Carson Recliner Chair

Carson Recliner provides a clean, modern aesthetic and is thoughtfully designed for healthcare environments, with a high level of comfort and ease-of-use for patients and caregivers.

- Arm caps available in wood, polyurethane, rigid polyester or solid surface
- Central locking caster system
- Trendelenburg feature
- Side-mount tablet with 360° swivel
- Flip-up side tablet
- Moisture barrier
- Field-replaceable upholstery covers
- Back construction features 14-gauge inner steel structure
- Lever positioning permits both caregiver and patient to operate footrest
- Ergonomically-positioned handles for easy transport

3-Position Patient Recliner



Specifications:

Product Code	5266
Width	740 mm
Depth	740 mm
Height	1120 mm
Seat Height	480 mm
Seat Width	584.2 mm
Arm Height	660 mm
SWL	225 kg



Specifications:

Product Code	5850
Seat Width	584 mm
Recline Depth	1880 mm
Seat Height	520 mm
Arm Height	660 mm
Overall Width	736 mm
Overall Height	1219 mm
SWL	225 kg



Oasis Mobile Recliner with Tray

Generously padded all over, and with a gas cylinder mechanism allowing infinite recline positioning, the Oasis recliner is a treat to relax in. Equally important to patient comfort is ease of operation for the care-giver. When developing the Oasis Recliner we gave extra consideration to this, which is evidenced in the features below.

- Ergonomically placed handles enable easier patient transfer
- Heavy duty design throughout for extra strength and product longevity
- Arm cap designed to assist patient ingress and egress
- Infinite recline positioning
- Side lever permits both patient and caregiver to operate leg extension
- Geometrically balanced mechanism substantially reduces effort required to operate leg extensions



Specifications:	
Product Code	5813
Seat Width	584 mm
Recline Depth	1880 mm
Seat Height	520 mm
Arm Height	692 mm
Overall Width	902 mm
Overall Height	1219 mm
SWL	225 kg

Regency Tilt & Recline

Packed full of features and with a 200kg weight capacity, this sleek recliner boasts a modern design appealing to wards and clinics. The fixed lounge style armrest has been uniquely designed to fit right up to most dining tables to increase versatility.

- 'Tilt in space' seat and leg rest action
- Gas assisted chair
- Euro style reclining backrest
- Elevating standard leg-rest
- Fixed lounge-style armrest
- Including push handle
- 4 x 100mm swivel braking castors
- Incorporating pressure management Multicush gel seating and air cushion backrest system
- Footrest option
- Australian made

* Also available in electric version



Specifications:	
Product Code	5024
Seat Width	550 mm
Seat Depth	520 mm
Seat Height	550 mm
SWL	200 kg



NDIS



NDIS

Regency Result3 Electric Recliner

Packed full of features and with a 200kg weight capacity, this sleek recliner boasts a modern design appealing to wards, aged care facilities and clinics. The fixed lounge-style armrest has been uniquely designed to fit right up to most dining tables to increase social interaction.

- Electric 'tilt-in-space' seat and leg-rest action
- Electric Euro-style reclining backrest
- Elevating standard leg-rest
- Incorporating pressure management Multicush gel seating and air cushion backrest system
- Fixed lounge-style armrest
- Auto lay-flat function, auto chair function, M3 customer set memory position
- 4 x 100 swivel braking castors
- 5 year warranty
- Australian made



Specifications:

Product Code	5025GH
Seat Width	550 mm
Seat Depth	520 mm
Seat Height	550 mm
SWL	200 kg

Regency Bahama Comfort Chair

Our Complex Care Chair is a popular choice for many healthcare environments. It's available in gas or electric formats, easy to maneuver and simple to clean for infection control. It's been carefully designed to keep your patients or residents social at meal times, fitting in closely to many dining tables.

- 'Tilt in space' seat and leg-rest action
- Euro style reclining backrest with built-in Bahama style head wings
- Elevating standard leg-rest
- Incorporates pressure management Multicush gel seating with air cushion backrest system
- Fixed lounge style armrest
- Includes push handle
- 180kg safe working load
- 5 year limited warranty
- Available in a huge range of Healthcare grade upholstery colours
- Australian made



Specifications:

Product Code	5022GH
Seat Width	550 mm
Seat Depth	570 mm
Seat Height	550 mm
SWL	180 kg

Also Available:

Regency Bahama Comfort Chair with Supports	
Product Code	5023GH

NDIS



NDIS



Cocoon G1 Chair

The lift function of the COCOON G1 armchair is very appreciated by people who have difficulty getting up, but who also find this function a real help when sitting down. Thanks to its coatings that meet fire-retardant standards and are easy to clean, the COCOON G1 can be used both at home and in the community. It consists of removable and washable zipped quilt modules at 30°. This electric lift chair is equipped with removable and / or retractable armrests to facilitate seated – standing transfers (evolving chair in the event of pathological development).

- Optimal hygiene: replaceable, removable and machine washable duvet modules at 30°
- Soft and ultra-comfortable duvets filled with siliconized hollow fiber
- 4 possible positions: Seated, Relaxed (TV), Rest and Lifter for standing transfer
- Different biocompatible coverings available: PVC/PU with M2 fire retardant treatment - mottled fabric - water-repellent velvet
- Different colours variable available depending on the coating
- Guarantees: chair: 2 years - cylinder and mechanics: 5 years
- Safe with superior manoeuvrability
- Adjustable footrests and handle heights
- One or Two Motors Available
- French manufacturing and design.



Specifications:

Product Code	5572
Seat Width	440 mm
Seat Depth	450 mm
Overall Width	930 mm
Seat Height	520 mm
Overall Height	1100 mm
Chair Weight	59 kg
SWL	160 kg

Cocoon G2 Chair

The COCOON G2 chair is an ultra-comfortable electric lift chair, removable and/or retractable multi-position to facilitate seated – standing transfers (evolving chair in the event of pathological development). The lift function of the COCOON G2 armchair is very appreciated by people who find it difficult to get up, but who also find this function a real help when sitting down. Thanks to its PVC coatings to fire-retardant standards, the COCOON G2 can be used both at home and in the community. Using the lifter is very simple. With a simple press of the remote control, the chair returns to a seated position and the chair begins its gentle ascent. You can release the pressure on the remote control at any time, your chair stops.

- It consists of removable and washable zipped quilt modules at 30°
- 4 possible positions: Seated, Relaxed (TV), Rest and Lifter for standing transfer
- Soft and ultra-comfortable duvets filled with siliconized hollow fiber
- Different biocompatible coverings available: PVC/PU with M2 fire retardant treatment - mottled fabric - water-repellent velvet
- 8 possible colors (variable depending on the coating)
- Different seat widths available
- Options available include: mobility kit - patient lift access kit - side shelf - magazine holder - High Density foam - head protection - battery kit
- Guarantees: chair: 2 years - cylinder and mechanics: 5 years
- French manufacturing and design.



Specifications:

Product Code	5575
Seat Width	440 mm
Seat Depth	450 mm
Overall Width	730 mm
Seat Height	520 mm
Overall Height	1100 mm
Chair Weight	59 kg
SWL	160 kg



Configura Chair

The Configura Chair has been designed as a highly adjustable seat that ensures a range of clients are kept comfortable and supported. The chair is adjustable in many ways - change your seat height, seat width and back rest to maximise comfort. Upholstered legrest with calf support

- Castors for easy manoeuvrability
- Pressure reducing fabric on backrest and armrest
- 'Reduced Shear' recline to maximise pelvic stability
- Vapour permeable fabric as standard on seat base, legrest, backrest and armrest
- Adjustable pillow-backrest
- Right hand control with rise setting

Optional Features:

- Different levels of pressure care seating cushions
- Different lateral backrest supports for different postural requirements
- Different control options
- Different headrest options



Specifications:


Product Code	5550M
Size	Medium
Seat Width	510 mm
Seat Height	460 mm
Armrest Height	180 mm
Overall Height	1140 mm
Overall Width	760 mm
SWL	160 kg

Also Available:

Size	Small	Large
Seat Width	460 mm	560 mm
Overall Width	710 mm	810 mm
SWL	130 kg	160 kg



Adjustable Lateral Support Backrest



Comfort Lateral Support Backrest



Multi-adjustable Pillow Backrest



Comfort Profiled Headrest

05.

Regency Chairs

Regency Care Chair	34
Regency R2700 Flotation Chair	35
Regency R2900 Flotation Chair - Low Seat	36

NDIS



Regency Care Chair

The 5050G Regency Care Chair is designed to give the client comfort by combining the original tilt system with an independent multi-positional back rest. Its many features take into account the needs of both carer and user of the chair. The Regency 5050G derives its status from a proven and successful track record across a broad spectrum of the industry.

- **Other sizes available.** Contact the team today to find out more.
- **Drop down armrests and side wings:** Regency armrests/wings fold right away to allow chair to be positioned close to bed for supine transfers. Also making sling insertion much easier.

- **Manoeuvrability:** Regency chairs have two front directional lock castors, two locking rear castors and a horizontal push handle reducing the effort to manoeuvre.
- **Gel cushion:** The guesswork has been removed ensuring the seat has the correct pressure to support the resident with the highest level of comfort.
- **Leg and footrest foldaway feature:** The uniquely designed non removable footplates prevents persons from sliding forward lowering the number of repositions required.
- **Manoeuvrable castor base:** 175mm castors reduces effort to move the chair.

Colours available:



Tangerine



Purple



Green



Mink



Lagoon



Mustard Pattern



Koala Pattern



Arcadia Marina



Bouquet Navy



Specifications:

Product Code	5050G
Seat Width	460 mm
Seat Depth	520 mm
Seat Height	600 mm
Backrest Length	700 mm
Armrest Height	175 mm
Legrest Length	480 - 650 mm
Overall Width	660 mm
Overall Depth	1040 mm
Overall Height	1310 mm
SWL	180 kg



Regency R2700 Flotation Chair

This Regency R2700 Series Care Chair is designed to be carer and client friendly by incorporating the original tilt seat system into a multi-positioning, multi-functional and manoeuvrable care chair. This particular model has a narrow seat width and shorter seat depth for the small adult or paediatric client.

- Superior range of tilt and recline in comparison with many products on the market
- Low seat height
- Pressure management seating – gel cushion in seat, air in backrest

- Forward seat tilt for ease to stand, legrest angles under seat
- Elevating legrest complete with fold out, extendible footplate
- Low skin shear 2-way stretch fabric on seating surfaces
- Drop down armrests and swingaway headwings for access
- Adjustable head rest cushion
- 125mm castors, braking on rear, directional lock on front
- Quality Australian made construction

Colours available:



Tangerine

Purple

Green



Mink

Lagoon

Mustard Pattern



Koala Pattern

Arcadia Marina

Bouquet Navy



Specifications:

Product Code	5030
Seat Width	410 mm
Seat Depth	420 mm
Seat Height	520 - 550 mm
Backrest Height	700 mm
Armrest Height	175 mm
Legrest Length	470 - 550 mm
Overall Width	660 mm
Overall Depth	1050 mm
Overall Height	1280 mm
SWL	180 kg



Regency R2900 Flotation Chair - Low Seat

The R2900 Series Regency Care Chair is designed to give the patient comfort by combining the original tilt system with an independent multi-positional back rest. Its many features take into account the needs of both carer and user of the chair.

The Regency R2900 derives its status from a proven and successful track record across a broad spectrum of the industry.

- Reclining backrest – full recline to supine position
- Tilting seat
- Pressure management seating – gel cushion in seat, self-air adjusting back, profiled head cushion
- 2-position armrests
- Forward seat tilt for ease to stand, legrest angles under seat
- Elevating contoured legrest complete with fold out, extendible footplate
- Low skin shear 2-way stretch fabric on seating surfaces
- Drop down armrests and swingaway headwings – side transfer friendly
- Covers replaceable
- Patient lifter friendly castor base
- 5 year limited warranty
- 175mm castors, braking on rear, directional lock on front
- Quality Australian made construction

Colours available:



Tangerine



Purple



Green



Mink



Lagoon



Mustard Pattern



Koala Pattern



Arcadia Marina



Bouquet Navy



Specifications:

Product Code	5040
Seat Width	460 mm
Seat Depth	520 mm
Seat Height	600 mm
Backrest Height	680 mm
Armrest Height	175 mm
Legrest Length	410 mm
Overall Height	1260 mm
SWL	180 kg

06.







Sertain Chairs

Sertain S4325 Treatment Chair	39
Sertain Hi-Lo Care Chair with Arms & Euro Back	39
Sertain Care Chair Electric Hi-Lo	40
Sertain Bariatric Hi-Lo Chair	41
Sertain S4540 Electric Hi-Lo Chair	41
Sertain I.V. Pole	42
Sertain Oxygen Bottle Carrier	42
Sertain Layflat Bar Option	42
Meal Tray	42

Choosing the right Care Chair /Hi-Lo Chair

We're passionate about ensuring that our products are specially tailored to meet the need of each individual and carer. That's why we've built a chair with parts that can be selected to give you the most out of your seating. Build your Sertain Chair below!

Build your Chair!

 Choose Model	Light Hi-Lo	Hi-Lo Sertain	
	Regency Results	Bahama Recliner	
	Treatment Chair	Bariatric Hi-Lo	
 Choose Size	Narrow	Standard	
	Large	Extra Large	
 Armrest Style	No Meal Tray Arm	Standard Removable Arm	
	Drop Arm	Procedure Arm	
	Procedure tray Arm		
 Backrest Style	40 Standard Series with Wings	40 Standard Series without Wings	
	45 Euro Series	50 Bahama Series	
 Legrest Style	Standard Legrest - No footplate	 Base Style	Standard Base
	Deluxe Legrest - Fold away footplate		Central Lock Base
 Accessories	Sertain IV Pole	Sertain Oxygen Bottle Carrier	
	Sertain Layflat Bar	Sertain Meal Tray	

Choose your Colour!



Could not find the information you were after?

Contact our team for a consultation for more Sertain Chair Options!

NDIS



Sertain S4325 Treatment Chair

The Sertain HILO Care Chair with procedure tray arms has been specifically designed to assist with cannulation. The HILO function of the chair allows for adjustment to a comfortable position for both the nurse and patient, while the armrest simplifies cannulation and blood collection.

- Electric 4 function (backrest, seat-tilt & leg-rest)
- Auto lay flat and auto chair function
- No head wings
- Procedure tray arm option
- Standard leg-rest
- Multicush pressure management seating system
- Hi-Low range: 550-850mmH
- Australian made

*Available in other colours and styles

NDIS



Sertain Hi-Lo Care Chair with Arms & Euro Back

Boasting a sleek Euro style backrest to add an extra touch of class to your ward, this Sertain HILO Care Chair is a popular choice for many hospitals wards. Available in a huge range of colours to compliment your ward, it promises an aesthetic touch without compromising on functionality.

- Electric 4 function (seat-tilt, backrest, Hi-Lo seat & legrest)
- Auto lay flat and auto chair function
- No head-wings – Euro-style backrest
- 160kg SWL (300kg SWL option also available)
- Standard removable armrest with no meal tray provision
- Standard push handle
- Hi-Low range: 550-850mmH
- Australian made

*Available in other colours and styles



Specifications:

Product Code	51281
Seat Width	585 mm
Seat Depth	570 mm
Seat Height	540 - 840 mm
SWL	160 kg



Specifications:

Product Code	5026GH
Seat Width	520 mm
Seat Depth	520 mm
Seat Height	550 - 850 mm
SWL	160 kg



Sertain Care Chair Electric Hi-Lo

The Sertain® Care Chair Series is there when you need it, to reposition the occupant, ease the transfer from bed to chair, increased access for wound dressing and procedures, manoeuvrable in small spaces.

- Electric hi-lo ; low enough to stand from with ease, but high enough to achieve a comfortable working height for nursing staff
- Electric reclining backrest – full recline for supine position
- Electric tilting seat and elevating legrest
- Handset operated with preset auto chair and auto layflat functions
- Pressure management seating – gel cushion in seat, air in backrestin seat, air in backrest
- Fold out, extendible footplate
- Low skin shear 2-way stretch fabric on seating surfaces
- Removable two-position armrests
- Swingaway three position headwings
- Adjustable head rest cushion
- 150mm castors, 3 x braking, 1 x directional lock
- Battery back up system
- Quality Australian made product



Specifications:

Product Code	5142
Seat Width	490 mm
Seat Depth	520 mm
Seat Height (Tilt/Normal)	540 - 940 mm
Backrest Height	700 mm
Armrest Height	150 mm
Legrest Length	440 mm
Overall Width	740 mm
Overall Depth	1080 mm
Overall Height	1200 - 1520 mm
SWL	160 kg





Sertain Bariatric Hilo Care Chair

For the bariatric client, this model will support up to 450kg, with functionality, comfort and the needs of both user and carer in mind. Sertain Bariatric Hilo Care Chair - Electric hi-low, backrest recline, tilting seat, auto chair & auto layflat functions, standard legrest - no foot plate, 2 position swing away head wings Rio Style, removable arms, pressure care seating (gel seat/air back), polysoft 2-way stretch vinyl, With NO Meal Tray Provision

- 150mm Dual Twin Wheel Castors
- Push Handle
- 450kg Weight Capacity
- Seat Size: 920W x 600D x 550Hmm
- 4 Columns
- Hi-Low Range: 550-850Hmm
- Australian Made

Sertain S4540 Electric Hi-Lo

For the bariatric client, this model will support up to 300kg, with functionality, comfort and the needs of both user and carer in mind. The electric operation of hi-lo adjustments mean transfers and patient care are able to be carried out at the correct working height by staff.

- Electric hi-lo
- Electric tilting seat and elevating legrest
- Handset operated with preset auto chair and auto layflat functions
- Pressure management seating – gel cushion in seat, air in backrest
- Fold out, extendible footplate
- Low skin shear 2-way stretch fabric on seating surfaces
- Removable two-position armrests
- Swingaway three position headwings, Rio Style
- Adjustable head rest cushion
- 150mm castors, 3 x braking, 1 x directional lock
- Battery back up system
- Quality Australian made product



Specifications:

Product Code	5150
Seat Width	920 mm
Seat Depth	600 mm
Seat Height	550 mm
Overall Height	550 - 850 mm
SWL	450 kg

Specifications:

Product Code	5147
Seat Width	580 mm
Seat Depth	520 mm
Seat Height (Tilt / Normal)	550 - 950 mm
Backrest Height	700 mm
Armrest Height	160 mm
Legrest Length	440 mm
Overall Width	790 mm
Overall Depth	1080 mm
Overall Height	1335 mm
SWL	350 kg



Sertain I.V. Pole

The IV Pole option allows the intravenous drip dependent client to stay seated in the comfort of their chair.

Specifications:

Product Code 5109



Sertain Layflat Bar Option

Supports chair at lie-flat position, during CPR emergency procedures. At time of manufacture.

Colour Options:



Specifications:

Product Code 5103



Sertain Oxygen Bottle Carrier

Oxygen Bottle Carrier, including bracket for mounting

Specifications:

Product Code 5108



Meal Tray

A popular option for the Sertain chair is the standard tray. Ideal for meals and activities Simple locking device secures tray in various positions.

- T-Mould Edge
- Cut-out means tray can be positioned close to user

Specifications:

Meal Tray to Suit 5120 5007S

Meal Tray To Suit 5131 5007S4

Meal Tray to Suit 5135 5007LG

07.

Bariatric

Embark Deluxe Extracare Day Chair	44
Utilite Extra Care Day Chair	44
Staebel Bariatric Sitting Chair	45
Deluxe Bariatric Bedside Commode	45
Eclipse Bariatric Manual Wheelchair	46
Staebel Gusto Heavy Duty Wheelchair	46

NDIS



Embark Deluxe Extracare Day Chair

The Extra Care Embark Deluxe has an impressive 250kg weight capacity and extra wide 60cm width, making it ideal for the heavier client. It still provides all the benefits of the Embark Deluxe, e.g. adjustable seat depth and height, and comes standard with tilt and move wheels.

- High backrest incorporating lumbar and lateral support
- Padded armrests
- Comfortable vinyl upholstered seat and back
- Adjustable seat height and depth
- Powder-coated steel frame with reinforcement for increased load



Specifications:

Product Code	5931BN3
Seat Width	600 mm
Seat Depth	450 - 530 mm
Seat Height	470 - 620 mm
Backrest Height	765 mm
Arm Height	220 mm
Overall Width	800 mm
Overall Depth	740 - 820 mm
Overall Height	1180 - 1330 mm
SWL	250 kg

NDIS



Utilite Extra Care Day Chair

The Extra Care Utilite is a heavy duty, multi purpose adjustable seat height and angle chair. Ideal for use in an area where adjustable seating is required, for example, lounge, dining room, bedroom, activity room, etc. This low back chair provides comfort, stability and practicality at an affordable price.

- Extra wide frame – 550mm seat width
- Shallow curved low backrest
- Padded armrests
- Comfortable vinyl upholstered seat and back
- Adjustable seat height to accommodate any user
- Powder-coated steel frame



Specifications:

Product Code	5911
Seat Width	550 mm
Seat Depth	430 mm
Seat Height	450 - 630 mm
Backrest Height	450 mm
Arm Height	230 mm
Overall Width	670 mm
Overall Depth	420 mm
Overall Height	850 - 1030 mm
SWL	200 kg



Staebel Bariatric Sitting Chair

Specially developed for hospital use, with a strong frame and high weight capacity, vinyl upholstery and tilt and move wheels to cater for the Bariatric patient.

- Adjustable height
- High Contoured backrest incorporating lumbar support
- Contoured seat
- Padded drop down, removable arms for side transfer
- Adjustable seat depth
- Standard with tilt and move wheels for ease of manoeuvrability

Deluxe Bariatric Bedside Commode

This bedside commode offers a comfortable and stable support, and is aesthetically pleasing with various colour options available. A skirt around the frame below seat disguises the seat, bowl and lid.

- 650mm seat width
- Padded armrests provide user comfort and support while seated, and when sitting and standing
- Commode bowl is removed easily and lid is provided to maintain a hygienic environment
- Powder coated frame for durability and infection control
- Fold away seat provides a comfortable day chair when commode function not in use
- Vinyl upholstery - avoid infection, andeasily cleaned



Specifications:

Product Code	5966
Seat Width	780 mm
Seat Depth	600 mm
Seat Height	460 - 540 mm
Backrest Height	675 mm
Arm Height	170 mm
Overall Width	900 mm
Overall Height	1000 - 1085 mm
Width Between Arms	790 mm
SWL	425 kg



Specifications:

Product Code	7963CH
Seat Width	650 mm
Seat Height	420 - 570 mm
Seat Depth	440 mm
Overall Width	805 mm
SWL	210 kg

NDIS

NDIS



Eclipse Bariatric Manual Wheelchair

The Eclipse is a rigid chair designed by 7 Canadian engineers for real Bariatric use, and ensures complete safety for the care staff attending the occupant.

Ultra-light in design and extended castor positioning for perfect weight balance, results in manual handling excellence for care staff when transporting the patient or resident through the facility.

- Removable armrest and leg rest to assist in preferred transfer.
- Rigid frame for strength and reducing frame movement and fatigue.
- Overall width is only 19cm extra on the seat width, keeping it streamlined for going through doorways
- Solid seat pan to overcome slung seat sagging issues.
- Anti-tip wheels for complete user safety.
- Colour options for easy identification of Bariatric equipment.



Specifications:

Product Code	4865
Seat Width	500 - 1200 mm
Seat Depth	400 - 550 mm
Seat Height	350 - 510 mm
Backrest angle	85 - 100°
Net/Gross Weight	20 kg
SWL	450 kg

Staebel Gusto Heavy Duty Wheelchair

The Staebel Gusto heavy duty wheelchair caters for the bariatric client. Featuring padded vinyl upholstery, heavy duty folding frame, 24 inch rear wheels and inset foot plates for client comfort, this wheelchair is ideal in the aged care and hospital environment. It comes standard with urethane tyres, eliminating the need to maintain tyre pressure. The swing away leg rests are removable and also adjust in length. The armrests will remove to enable side transfers if necessary. Footplates are alloy for longevity and ease of cleaning.

- Removable, swingaway, adjustable length legrests
- Removable armrests
- Padded vinyl upholstery
- 24" rear mag type wheels, urethane tyres



Specifications:







Product Code	46108	46608	47108
Seat Width	610 mm	660 mm	710 mm
Seat Depth	460 mm	520 mm	550 mm
Seat Height	520 mm	520 mm	520 mm
Backrest Height	430 mm	480 mm	520 mm
Armrest Height	220 mm	235 mm	220 mm
Overall Width	850 mm	880 mm	930 mm
Overall Depth	1100 mm	1150 mm	1170 mm
Overall Height	930 mm	1000 mm	990 mm
Fold Width	350 mm	350 mm	350 mm
Net/Gross Weight	26 / 30 kg	36 / 43 kg	36 / 46 kg
SWL	240 kg	270 kg	290 kg

08.

Cushion Matrix

The cushion matrix shows the relationship between our cushions and different variables (for example; weight limit, skin integrity or sizes) in a data set in grid format. Essentially, the cushion matrix is a quick guide which allows the reader to compare each cushion to help find the appropriate cushion for either their needs or patient needs in a time efficient matter.

Cushion Matrix

	Weight Capacity	Cushion Profile Height	Size Range	
			Width - inches	Width - mm
3210 Budget 	100 kg	50 mm		
SAC 	160 kg	25 mm		
Action Pilot 	Unlimited	30 mm		
Action Professional 	Unlimited	45 mm		
Relax Easy Low Profile 	130 kg	50 mm	14 x 16" to 22 x 18"	360 x 400mm to 550 x 450mm
Relax Easy High Profile 	190 kg	85 mm	16 x 16" to 20 x 20"	400 x 400mm to 500 x 500mm
Duogel 	130 kg	85 mm	14 x 14" to 24 x 18"	360 x 360mm to 600 x 450mm
Duo Profile 	130 kg	85 mm	16 x 16" to 18 x 18"	400 x 400mm to 450 x 450mm
Gel Cell 	130 kg	40 mm	15 x 14" to 20 x 18"	380 x 360mm to 500 x 450mm

PI Rating	Thigh Support	Contoured seat	Skin Integrity			Boney Prominence	Suitable for Lounge chairs
			Skin Intact	Skin Broken			
Prevention	✗	✗	✓	✗	✗	✗	
Stage 2	✗	✗	✓	✗	✗	✓	
Stage 3	✓	✗	✓	✗	✗	✓	
Stage 4	✓	✗	✓	✓	✓	✗	
Stage 2	✗	✗	✓	✗	✓	✓	
Stage 3	✗	✗	✓	✓	✓	✗	
Stage 3	✓	✓	✓	✓	✓	✗	
Stage 2	✓	✓	✓	✗	✓	✗	
Stage 4	✗	✓	✓	✓	✓	✓	

We are
Here to Help



Let's have a chat!

If you see something you like, you can book in your consultancy with one of our friendly team members. See below information for ways you can get in contact...

1800 267 267 | www.activemedicalsupplies.com.au | hello@activemedicalsupplies.com.au

activeMEDICAL
SUPPLIES > SOLUTIONS > SUPPORT

P. 1800 267 267 | **F.** 07 3899 8825
E. hello@activemedicalsupplies.com.au

www.activemedicalsupplies.com.au