

Standard Transfer Bench



Designed for individuals with restricted mobility, the Standard Transfer Bench provides safety and independence while transferring in and out of the bath tub. This durable aluminium transfer bench is adjustable in height.

Product Benefits:

- 3 piece plastic non-slip seat with drain holes
- Hand rail for added support while transferring
- Reversible backrest to suit every situation
- Welded frame

Specifications:

Product Code	1018J
DVA	B5910
Seat Width	670 mm
Seat Depth	410 mm
Seat Height	445 - 570 mm
Bath Rim Clear	370 - 495 mm
Seat-Back Height	410 mm
Overall Width	790 mm
Overall Depth	540 mm
SWL	140 kg