

Shower Stool with Arms & Back



The shower stool with arms and backrest, provides security when entering or exiting the shower and assist the user on and off the seat.

Product Benefits:

- Compact size for the smaller shower recess
- Padded seat and backrest for user comfort
- Adjustable seat height for the individual's requirements

Specifications:

Product Code	1569
Seat Width	370 mm
Seat Depth	290 mm
Seat Height	520 - 670 mm
Width Between Arms	640 - 840 mm
Armrest Height	640 - 840 mm
Overall Width	570 mm
Overall Depth	450 mm
SWL	125 kg