

# Staebel P5500 Homecare Bed



# Staebel P5500 Homecare Bed

The Staebel P5500 Homecare Bed is designed for easy transportation and storage. Ideal for short term cases, the P5500 bed can be assembled and dismantled with ease through the incorporation of a simple lock system which is secured with a pin, and folds compactly for storage. Consequently, the P5500 is an ideal home care bed, as well as a short-term care alternative.

The bed is supreme in its design targeted at easy transportation & storage. This is the ideal bed for transfer into homes compared to the conventional bed. Unique fold-up design with pin slot system that allows bed to be transported simply from one location to another. Can be assembled and dismantled with ease and stored compactly.

## Product Features:


- Electric Hi/Low
- Electric Kneebend
- Electric Backrest
- 100mm Breaking Castors
- Ironstone Powdercoat with Oak Head & Footboards
- Trendelenberg
- Manual Vascular Support



**Electric**  
Backrest



**100mm**  
Breaking Castors



**Ironstone**  
Powdercoat



**Oak Head &**  
Footboards



**Transportable Design**  
for easy transport & Storage

## Optional Features – Added at time of manufacture

- Battery Backup

## Accessories – Added after manufacture

- Self help pole & Triangle
- Horizontal Rails
- Clamp on Support Rail
- StandEzy – Straight Handle
- StandEzy – Hook Handle
- Padded Bed Rail Protector
- IV Pole
- Reading Light
- Clip on Extension



## Specifications:

Product Code	7020	7021
Size	Single	King Single
Overall Dimensions	930 x 2190 mm	1070 x 2190 mm
Sleeping Surface	900 x 2000 mm	1070 x 2000 mm
Height Range	295 – 700 mm	295 – 700 mm
Back / Knee Angle	66° / 35°	66° / 35°
TR / ATR	11° / -11°	11° / -11°
SWL	200 kg	200 kg