


Secura 1 Bed



Secura 1 Bed

The Secura 1 is a robust bed for demanding environments, providing a sleeping surface which can reduce harm to both the patient and carers. There are no electrics included, hence creating a safe solution for your facility, with a SWL of 200kg and durable framework. Suitable for your medium-high care patients.

Product Features:

- No electrics
- Fixed height solid frame
- Secure bolt on head and footboards
- Strong robust construction
- Enclosed mattress supports
- Anti-tamper screws
- Natural Oak 25mm laminate head and footboards
- Ironstone Powdercoat




Enclosed
Mattress Supports



Anti Tapper
Screws



Natural Oak
Laminated Head
and Footboards



Ironstone
Powdercoat

Optional Features - Added at time of manufacture

- Floor mounting plates
- Custom powdercoat colours (APO Grey shown)
- Extra length
- Head and Footboard Designs

Accessories - Added after manufacture

- Mental Health Mattress

Specifications:

Size	Single	King Single
Overall Dimensions	920 x 2070 mm	1090 x 2070 mm
Sleeping Surface Dimensions	900 x 2000 mm	1070 x 2000 mm
Mattress Base Height	450 mm	450 mm
SWL	200 kg	200 kg