

Horzion 3004 Change Table



The Horizon 3004 Change Table is ideal for children and small adults up to 90kg. The foot control operation provides flexible and easy to control height adjustment, while the single tower lifting column conceals electrics within the column for safe and smooth hi-low functioning.

Product Features:

- Electric Hi/Low
- Foot Control Operated
- 125mm Breaking and Directional Castors
- Charcoal Vinyl Top
- APO Grey and Ironstone Powder Coated Frame

Optional Features (added at time of Manufacturer):

- Hand Control (as shown)
- Custom Colour Padded Top
- Custom Colour Powder Coated Frame



Accessories Available:

- 7831 Padded Side Rails (as shown)
- Lap Belt



Specifications:

Product Code	7830
Width	700 mm
Length	1400 mm
Min Height	480 mm
Max Height	880 mm
SWL	90 kg