

Horizon 3250 Change Table



The Horizon 3250 Change Table is ideal for multiple environments. Designed and manufactured in Australia, our 3250 Change Table would be a pivotal asset to any school, sports facility, medical rooms or assisted accommodation. The lift function is a vertical lift, and for added flexibility, the backrest is still available as an option on this change table.

Product Features:

- 125mm braking and directional castors
- Battery backup
- Operated by a handset
- Charcoal padded vinyl top
- Ironstone base frame powdercoat
- APO Grey top frame powdercoat

Optional Features (added at time of Manufacturer):

- Backrest (pictured on the left)
- Central Locking Castors
- Face Hole
- Foot Control
- Custom Colour Padded Top
- Custom Colour Powder Coated Frame



Accessories Available:

- 7831 Padded Side Rails (as shown)
- Lap Belt



Specifications:

Product Code	7805	7840
Overall Dimensions	700 x 1880 mm	700 x 1880 mm
Height Range	450 - 900 mm	450 - 900 mm
Back Angle	-	75°
SWL	250 kg	250 kg