


Alrick EN9 Floorline Bed



Alrick EN9 Floorline Bed


The Alrick EN9 Floorline Bed is the first ever universal bed that suits low care, high care and bariatric care, this innovative three-in-one bed has a low height of 100mm and no towers! Suitable for care of elderly or dementia patients, through to acute and bariatric settings, the Australian made EN9 Endless Bed has you covered. One bed, endless possibilities!

Product Features:

- Electric Hi Low
- Electric Auto Regression Backrest
- Electric Kneebend
- Electric Forward & Reverse Trendelenburg
- Extremely Low Minimum Position
- Vertical Lift
- Auto-Regression Backrest
- Ironstone base and the APO Grey top as standard
- Removeable Slender Arch Natural Oak Head & Footboards
- Australian Owned and Made Product
- Safe Working load from 200-250kg



Chair Position
push of a button



Bluetooth
handset must be ordered separately



75mm
Dual Castors



Vertical Lift
Range

Optional Features- Added at time of manufacturer

- 300kg Handset (Single & King Single only)
- Custom Powdercoat Colours (as shown)
- Bluetooth Handset
- Dual Breaking
- Head & Footboard Options
- Underbed Light
- Battery Back Up
- Attendant Control
- CPR Backrest



Accessories- Added after manufacturer

- Mattresses
- Reading Light
- Emergency Power Supply
- Horizontal Rails
- Combo Grab Handle
- Self Help Pole & Triangle
- I.V Pole
- T-Bar Wall Protection
- StandEzy Bed Pole
- Bed Extension



Optional Headboards - Custom Options Available



Laminate
Curved Style




Upholstered
Jersey Headboard &
Plain Footboard




Timber
Eltham Style


Standard Colours



Natural Oak
Laminate




Ocean
Upholstered




Ash
Upholstered



Burnt Walnut
Timber



Beech
Timber



Saddle
Timber

Specifications:

Size	Single	King Single	Double	Queen	King
Overall Dimensions	920 x 2115 mm	1090 x 2115 mm	1390 x 2115 mm	1520 x 2115 mm	1820 x 2115 mm
Sleeping Surface Dimensions	900 x 2000 mm	1070 x 2000 mm	1370 x 2000 mm	1500 x 2000 mm	1800 x 2000 mm
Hi - Lo Position	100 - 725 mm	100 - 725 mm	100 - 725 mm	100 - 725 mm	100 - 725 mm
TR / ATR Angle	15° - 15°	15° - 15°	15° - 15°	15° - 15°	15° - 15°
Back / Knee Angle	62° / 35°	62° / 35°	62° / 35°	62° / 35°	62° / 35°
SWL	200 kg	200 kg	250 kg	250 kg	250 kg