

Alrick Nino Bed



Alrick Nino Bed

Our Alrick Nino Bed suits low-care and high-care. This innovative three-in-one bed has a low height of 100mm and a high height of 725mm. Available in a wide range of powdercoat colours and with the option of higher doors/sides, or doors on both sides.


The NINO allows the user to perceive the world around them while at the same time providing protection and security. The transparent Perspex merges with the room allowing occupants to have an unrestricted view of their immediate surroundings yet still protects against falls and injuries. One bed, endless possibilities!

Product Features:

- 2 Choices of powdercoat colour
- Chair Position at the push of a button
- Auto-regression Back Rest
- Trendelenburg and Anti Trendelenburg
- 75mm Castors
- 750mm High Doors (one side standard)
- Vertical Lift
- Battery Back Up System
- Australian Owned and Made
- Safe Working Load of 200kg



Transparent
Perspex Inserts



Bluetooth Connectivity
handset must be ordered separately



75mm
Twin Castors



Low Position
of 100mm

Powdercoat Colours Available:

Use this guide when choosing your 2 powdercoat colours. Custom options available



Colour 2
Doors / Sides

Colour 1
Bed base and End panels



Optional Features - Added at time of manufacture

- Custom Powdercoat colours
- Padded Sides
- Head & Footboard Options
- Bluetooth Handset
- Meshtex Sides
- Vertical Uprights on Doors (as shown)
- Underbed Light
- Optional 1100mm Height Sides



Specifications:

Size	Single	King Single
Overall Dimensions	2125 x 1025 mm	2125 x 1195 mm
Sleeping Surface Dimensions	2000 x 900 mm	2000 x 1070 mm
Low Position	100 mm	100 mm
High Position	725 mm	725 mm
Underbed Clearance	100 mm (at lowest point)	100 mm (at lowest point)
Back / Knee Angles (°)	62° / 35°	62° / 35°
Trendelenburg Angles (°)	- 15° / 15°	- 15° / 15°
SWL	200 kg	200 kg