

Staebel ComfortCare Bed Base Sheet



The ComfortCare is an in-bed sliding system for the care of the long term or 24 hour patient. The amazing ComfortCare offers patients 100% comfort and carers the ability to manoeuvre and reposition them in all directions with ease and then secure that position when needed.

Product Benefits:

- Moves a patient in all directions
 - The locking panel will hold the patient steady
 - Offers patients comfort and dignity
 - Patients can remain on this indefinitely, either directly or under a cotton bed sheet
 - Carers can re-position and manoeuvre patients without the need to disturb or handle directly
 - Reduces risk of injury to carers
 - Enhances patient care and clinical outcomes
 - Suitable for use with hoists
- Polycotton system, with low-friction inner linings
 - Two panel system enables up/down manoeuvres, turning, swivelling, sitting up safely, transferring and turning into prone position
 - Integral locking panel means safe position can be secured, even in sitting
 - Suitable for single carers
 - Washable at 40°C to 80°C and low temperature tumble dry
 - Spare panels available

Specifications:

Product Code	210296
Width	900 mm
Length	2000 mm
SWL	225 kg