

# Staebel Perfect Positioner



Originally designed to help reposition a tetraplegic gentleman in a moulded wheelchair, this unique handling device has been redesigned and now enables carers to easily and comfortably reposition the user (once hoisted into the wheelchair) accurately at the back of the chair and ensure their pelvic tilt is correct.

## Product Benefits:

- Especially useful in moulded and specialised wheelchairs
- Correct sitting position in wheelchairs
- Accurate pelvic tilt
- One carer or two

## Specifications:

<b>Product Code</b>	210340
<b>Width</b>	650 mm
<b>Length</b>	2600 mm