

# EasyGlide Transfer Board



The flexible three-part sliding board is easy to place in front of the wheels of a wheelchair to protect against soiling during a sliding transfer. It is especially suitable for levelling out height differences or for bridging short distances.

## Specifications:

<b>Product Code</b>	215020
<b>Width</b>	330 mm
<b>Length</b>	600 mm
<b>Thickness</b>	5 mm
<b>SWL</b>	135 kg

## Product Benefits:

- Anti-slip protection on the underside for safety and security.
- The long sides can be folded to protect clothing from the wheelchair's wheels.
- Easy to transport
- Made of strong polyethylene and is available in different sizes.

## Sizes Available:

<b>Product Code</b>	215050
<b>Width</b>	330 mm
<b>Length</b>	750 mm
<b>Thickness</b>	6 mm