



# Allegro Salsa Sling



The Salsa sling is for use with the Salsa stand up hoist. It features optimum back support and under arm padding for a safe and comfortable lift and also allows for the easy removal of clothing for toilet transfers.

It features a torso band for security and stability with Buckle (b) or Velcro (v) fastening.

	<b>S</b>	<b>M</b>	<b>L</b>	<b>XL</b>
<b>A</b>	730	880	1060	1220
<b>B</b>	260	290	370	410
<b>C (b)</b>	900	960	1020	1220
<b>C (v)</b>	1140	1180	1320	1440



<b>Buckle Fastening</b>		<b>Velcro Fastening</b>	
<b>S</b>	227701	<b>S</b>	27720
<b>M</b>	27702	<b>M</b>	27705
<b>L</b>	27703	<b>L</b>	27706
<b>XL</b>	27704	<b>XL</b>	27707