

Salsa 200 Standing Lifter



**Compact Active.
Standing Aid.**

Salsa 200 Standing Lifter

Designed as an active standing aid allowing a single carer to raise a patient into a comfortable secure standing position. The Salsa 200 is ideal for transferring from bed to chair, or chair to commode. Allegro Concepts standing slings are manufactured to suit the Salsa, allowing superior access for dressing and personal hygiene.

With a small footprint and low overall weight, the Salsa 200 is highly maneuverable. A high capacity safe working load of 200kg means the Salsa is equally at home in an institutional or home environment.

Designed and tested to comply with:
AS ISO 10535-2002 Hoists for the transfer of disabled persons – Requirements and test methods.

Product Benefits:

- Australian made quality product
- Electric leg adjustment
- Patient friendly lift arms enable multiple grip positions during lift
- Removable footplate for ease of transfers
- Height and angle adjustable contoured knee pad with removable calf strap
- Comes complete with battery pack and charger
- T.G.A. 136 065
- Safe Working Load of 200 kg



Specifications:

| | | | |
|-----------------------------------|----------|---------------------------------------|-------------------|
| Product Code | 27700 | Under Base Height | 55 mm |
| Base Length | 1010 mm | Lifter Weight | 51 kg |
| External Base Width Closed | 665 mm | Minimum Storage Height | 1190 mm |
| External Base Width Open | 1030 mm | Footplate | Removable |
| Internal Base Width Closed | 565 mm | Kneepad | Height Adjustable |
| Internal Base Width Open | 930 mm | SWL | 200 kg |
| Base Width Adjustment | Electric | Also Available: | |
| Under Bed Clearance | 130 mm | Salsa Handle Kit, After Market | |
| Minimum Height Lift Arm | 880 mm | Product Code | 27THK01 |
| Maximum Height Lift Arm | 1685 mm | | |