

Sonata 150 Lifter



**Superb value.
Even better quality.**

Sonata 150 Lifter

Product Features:

A lightweight, compact mobile patient lifting hoist with innovative and ergonomic design the Sonata 150 is an easy to use everyday lifter. With a safe working load of 150 kg the Sonata offers an outstanding lift capacity in a small patient lifter. It utilises the SPACELINK™ concept, a breakthrough for optimised patient space, that enhances floor lift capability and increases maximum lifting height.

With maintenance free bushes in all the major moving parts this hoist is an Australian Made example of quality and strength that is hard to find in other patient lifters.

The Allegro range of patient lifting equipment including lifters and accessories has been designed and tested to comply with: AS ISO 10535-2002 Hoists for the transfer of disabled persons – Requirements and test methods.

Product Benefits:

- Manual leg adjustment is a standard feature.
- Supplied complete with battery, in-built battery charger and instruction manual.
- Sling attachment sold separately.
- Safe Working load of 150kg



Specifications:

Product Code	2880	Hoisting Range*	870 mm
Base Length	1005 mm	Minimum Storage Height	1260 mm
External Base Width Closed	605 mm	Under Base Height	55 mm
External Base Width Open	1070 mm	Maximum Hoist Reach	875 mm
Internal Base Width Closed	510 mm	Minimum Hoist Reach	555 mm
Internal Base Width Open	975 mm	Hoist Reach at Maximum Height	555 mm
Base Width Adjustment	3 Position Manual	Mass of Heaviest Part (Mast Assembly)	22 kg
Under Bed Clearance	130 mm	Lifter Weight	41 kg
Height of Yoke	420 - 1735 mm	SWL	150 kg