

# Staebel Deluxe Mobile Overbed Table



The Viva Overbed Table offers the latest in smart and modern design, with an easy clean, seamless PE top to promote ease of use and satisfy infection control requirements. Featuring smooth lines, no joins or hard to reach contours, this innovative design makes previously hard to maintain items easier and quicker to clean.

## Product Benefits:

- Easy gas lift height adjustment, non-tilt
- PE easy clean top, 34mm deep, max 20kg load.
- Base height 98mm
- Four castors, including two with brakes
- Supplied fully assembled

## Specifications:

<b>Product Code</b>	4210
<b>Top Width</b>	900 mm
<b>Top Depth</b>	430 mm
<b>Minimum Height</b>	725 mm
<b>Maximum Height</b>	1050 mm
<b>SWL</b>	20 kg

# Real Life Benefits of the Staebel Deluxe Mobile Overbed Table



1

**Ultra-low base** to fit under as many beds as possible.

2

**High quality castors** to remove any jerky motion when the meal/patient belongings are on the table.

3

**Gas assist raise and lower** making the operation one-handed and extra easy for staff.

4

One-piece **injection moulded top** for ultimate lifespan and infection control/cleaning.

5

**Extra long top**, allowing extra space for patient requirements, whilst reaching further across the bed.

6

Ability to **vinyl wrap** with your own unique message.