

Staebel Sovereign Wheelchair



**The most universal
wheelchair - Adjustable
to suit any user.**

Staebel Sovereign Wheelchair

Product Features:

The Sovereign wheelchair by can be adjusted to suit almost any user. With much consideration in the design from the perspective of OTs, this chair is suited to the vast range of needs required by our clients.

Rather than numerous script forms for various users, this chair meets the needs of almost any patient. All wheelchair parts suit any adjustments made - no additional parts required.



Product Benefits:

- Adjustable & removable legrests
- Flip back & removable armrests
- Adjustable armrest height & length
- Seat depth, height & angle adjustment
- Backrest height, tension & angle adjustment
- Crash Tested
- Safe Working Load ranges from 125kg to 170kg



Specifications:

| Product Code | 4340C-37 | 4340C-40 | 4340C-43 | 4340C-46 | 4340C-49 | 4352C-52 | 4352C-55 | 4352C-58 |
|--|--------------|--------------|--------------|--------------|--------------|--------------|--------------|--------------|
| Seat Width | 430 mm | 460 mm | 430 mm | 460 mm | 490 mm | 520 mm | 550 mm | 580 mm |
| Anterior Seat Height | 420 mm | 445 mm | 420 mm | 445 mm | 450 mm | 470 mm | 475 mm | 500 mm |
| Posterior Seat Height | 410 mm | 430 mm | 410 mm | 430 mm | 450 mm | 470 mm | 490 mm | 510 mm |
| Seat Depth | 400 mm | 425 mm | 400 mm | 425 mm | 450 mm | 475 mm | 475 mm | 500 mm |
| Backrest Height (2cm increments) | 400 - 460 mm | 400 - 460 mm | 400 - 460 mm | 400 - 460 mm | 400 - 460 mm | 400 - 460 mm | 400 - 460 mm | 400 - 460 mm |
| Backrest Angle (5° increments) | 69° - 94° | 69° - 94° | 69° - 94° | 69° - 94° | 69° - 94° | 69° - 94° | 69° - 94° | 69° - 94° |
| Lower Leg Length (2cm increments) | 400 - 480 mm | 400 - 480 mm | 400 - 480 mm | 400 - 480 mm | 400 - 480 mm | 400 - 480 mm | 400 - 480 mm | 400 - 480 mm |
| Armrest Height (2cm increments) | 180 - 260 mm | 180 - 260 mm | 180 - 260 mm | 180 - 260 mm | 180 - 260 mm | 180 - 260 mm | 180 - 260 mm | 180 - 260 mm |
| Armrest Length (4cm increments) | 330 - 450 mm | 330 - 450 mm | 330 - 450 mm | 330 - 450 mm | 330 - 450 mm | 330 - 450 mm | 330 - 450 mm | 330 - 450 mm |
| Footplate Width | 105 - 155 mm | 105 - 155 mm | 105 - 155 mm | 105 - 155 mm | 105 - 155 mm | 105 - 155 mm | 105 - 155 mm | 105 - 155 mm |
| Gross Weight | 17 kg | 17 kg | 17 kg | 17 kg | 17 kg | 21 kg | 21 kg | |
| SWL | 125 kg | 125 kg | 125 kg | 125 kg | 125 kg | 170 kg | 170 kg | 170 kg |