

Deluxe Alloy Wheelchair



Versatility for every Patient.

Deluxe Alloy Wheelchair

The Deluxe Alloy Wheelchair, is a wheelchair that has been designed to accurately meet the needs of a hospital. A light-weight wheelchair with the benefit of many additional features helps us assist patients requiring transfers.

The improvement compared to previous designs (ranging from 25 – 30kg gross weight) is considerable. Much input from patients and nurses has assisted in the development of these superior quality wheelchairs.

Product Benefits:

- Flip back armrests
- Swingaway, detachable legrests
- Fold down backrest
- 175mm castors reduce effort required to move the chair.
- Quick release rear wheels
- 24" rear wheels with urethane tyres
- Aluminium frame
- Padded nylon upholstery
- Complies with AS3695:2008



Flip back armrests




Fold down backrest




Height Adjustable
Armrests




Oxygen Carrier



Double Bracing



Urethane Tyres



No Upholstery
Screws



IV Pole

Specifications:

Product Code	4629
Seat Width	460 mm
Seat Depth	400 mm
Seat Height	480 mm
Backrest Height	430 mm
Armrest Height	230 mm
Overall Width	650 mm
Overall Depth	1060 mm
Overall Height	960 mm
Fold Width	320 mm
Net/Gross Weight	11 / 14 kg
SWL	125 kg

Sizes Available:

Product Code	45129	45629
Seat Width	510 mm	560 mm
Seat Depth	400 mm	400 mm
Seat Height	490 mm	490 mm
Backrest Height	430 mm	430 mm
Armrest Height	210 – 280 mm	210 – 280 mm
Overall Width	700 mm	750 mm
Overall Depth	1000 mm	1000 mm
Overall Height	960 mm	960 mm
Fold Width	300 mm	300 mm
Net/Gross Weight	14.5 / 17 kg	15 / 17.5 kg
SWL	150 kg	150 kg