

# Embark Bariatric Chair



The Embark Bariatric Chair was specially developed for Acute Care/Aged Care residence use, with a strong frame and high weight capacity, vinyl upholstery and tilt and move wheels to cater for the bariatric patient.

## Product Benefits:

- Adjustable height
- High Contoured backrest incorporating lumbar support
- Contoured seat
- Padded drop down, removable arms for side transfer
- Adjustable seat depth
- Standard with tilt and move wheels for ease of manouvability

## Specifications:

<b>Product Code</b>	5966
<b>Seat Width</b>	780 mm
<b>Seat Depth</b>	600 mm
<b>Seat Height</b>	460 - 540 mm
<b>Backrest Height</b>	675 mm
<b>Arm Height</b>	170 mm
<b>Overall Width</b>	900 mm
<b>Overall Height</b>	1000 - 1085 mm
<b>Width Between Arms</b>	790 mm
<b>SWL</b>	425 kg