

# Trust Care Outdoor Walker



It is light and elegant, both folded and unfolded. It has an X frame cross folding system which makes it easy to fold, takes up little space and stands alone when folded. Big front and rear wheels make it easy to get over curbs and other obstacles. Even more important is the spring system in the wheels that gives a smooth movement compensating for stiff joints.

## Product Benefits:

- Soft, comfortable rubber handles that are height adjustable
- Easy to manage brake which can be locked
- Comfortable seat in imitation leather for strength.
- Elegant design in every aspect.
- Available in Black/Silver (60120) or Mother-of-Pearl/Silver - pictured (60150)

## Specifications:

<b>Product Code</b>	60150
<b>Handle Height</b>	780 - 960 mm
<b>Width</b>	580 mm
<b>Walker Weight</b>	6.4 kg
<b>SWL</b>	100 kg