

Bariatric Height Adjustable Crutches



The Bariatric Crutches offer push button height adjustment, providing support, strength and comfort for consumers up to 250kg.

Specifications:

Product Code	6020
Crutch Height	940 - 1150 mm
Patient Height	1219 - 1372 mm
SWL	250 kg

Sizes Available:

Product Code	6021	6022
Crutch Height	1130 - 1340 mm	1330 - 1540 mm
Patient Height	1372 - 1575 mm	1575 - 1778 mm
SWL	250 kg	250 kg