

Coopers Double Adjustable Crutches



The Coopers Elbow Crutches offer push button height adjustment, manufactured in a lightweight anodised aluminium, double adjustment for individual fitting of hand grips and cuffs. The triangular shape and angle of the handle reflect the natural shape of the hand and result in easier use.

Specifications:

Product Code	6120
Cuff to Hand Grip	180 – 255 mm
Hand Grip Height	555 – 810 mm
SWL	180 kg

Product Benefits:

- Tip size 22mm
- Supplied as a pair
- Fitted with anti-rattle bushes.

Sizes Available:

Product Code	6121	6122
Cuff to Hand Grip	200 – 275 mm	200 – 275 mm
Hand Grip Height	685 – 935 mm	840 – 1095 mm
SWL	180 kg	180 kg