

# Coopers Comfort Crutches



The Coopers Elbow Crutches with comfy handles feature moulded handles with a unique design giving pressure distribution characteristics greatly superior to any traditional crutch. Adjustment is made easy with push button operation, and double adjustment for individual fitting of hand grips and cuffs is incorporated.

## Product Benefits:

- Fitted with anti-rattle bushes.
- Supplied as a pair

## Specifications:

<b>Product Code</b>	6141
<b>Cuff to Hand Grip</b>	210 - 290 mm
<b>Hand to Grip Height</b>	710 - 960 mm
<b>SWL</b>	125 kg