

# Walking Stick Swan Neck



This adjustable height aluminium walking stick assists the user's mobility. The "swan neck" shaped handle helps to transmit the weight of the user vertically through the body of the stick. It provides exceptional support and comfort, eliminating slipping and encouraging correct posture. It also includes an anti rattle silencer and adjusts in 25mm increments.

## Product Benefits:

- Padded foam handle with carrying loop
- 100kg max user weight

## Specifications:

<b>Product Code</b>	6220
<b>Minimum Height</b>	730 mm
<b>Maximum Height</b>	990 mm