

Quad Walking Stick



These offset base Quad sticks are available in two sizes of low and profile base. Designed with a swan neck which helps to transmit the weight of the user vertically through the body of the stick. Rotate handle 180°; for left or right hand use. Excellent for users requiring a higher level of support. Includes anti rattle silencer.

Product Benefits:

- Adjustable height to suit user requirements
- Swan neck handle allows for even weight distribution
- Soft foam hand grip for added comfort
- Low profile offset base does not obstruct user's step
- Locking collar prevents rattle
- Lightweight aluminium and steel construction

Specifications:

Product Code	6250
Base Size	Small
Height	720 - 970 mm
SWL	100 kg

Sizes Available:

Product Code	6012J
Base Size	Large
Height	720 - 970 mm
SWL	100 kg