

Cruiser Folding Walker



The Cruiser is a popular seat walker manufactured from durable lightweight aluminium. 6" castors make this an ideal outdoor walker. But is also compact enough for use indoors.

Features include a padded seat, bag under seat, lockable brakes and padded, folding backrest for comfort while seated. The adjustable height handles are fitted with orthopaedic hand grips, therefore providing the user with excellent support. The walker also folds flat for transport.

Product Benefits:

- Height adjustable ergonomic shaped hand grips
- Lightweight aluminium construction
- Lockable handbrakes
- 152mm castors

Specifications:

Product Code	6730
Seat Width	440 mm
Seat Height	560 mm
Handle Height	820 - 930 mm
Overall Width	570 mm
Walker Weight	6 kg
SWL	125 kg