

Galaxy Seat Walker



Constructed from heavy duty aluminium, this bariatric seat walker can cater for those up to 220kg. It features a removable padded backrest and ergonomically shaped handgrips for user comfort. The handbrakes are lockable for stability. A wire basket is included also.

Product Benefits:

- Extra wide seat width and reinforced frame
- Folds for easy storage and transportation
- Height adjustable handles
- Handbrakes can be locked for added stability
- Supplied with removable backrest, basket and padded seat

Specifications:

Product Code	6736
Overall Width	710 mm
Width Between Handles	610 mm
Handle Height	770 - 920 mm
Walker Weight	11.1 kg
SWL	220 kg