

# Stimulator Forearm Support Walker



**We've thought of the functions  
and practicality to efficiently get  
your walking and rehab underway.**

# Stimulator Forearm Support Walker

## Product Benefits:

"The best rehab program and the best rollator will guarantee you the best results – is this what you're looking for? Our Stimulator Rehab Forearm Walker has been specifically designed for the needs of the hospital environment.

- Very lightweight, easy for the frail, weak patient to use
- High quality castors with super smooth rolling and minimal friction even on carpet, taking away any shuddering or wobbling
- Easy, simple and very gentle lever height adjustment
- Soft PU multi adjustable arm pads for any patient with any specific position requirement
- Multi adjustable hand grips, for the arthritic, ABI or stroke patients.
- Overall width designed to fit through nearly any doorway
- Four swivelling castors for ultimate manoeuvrability
- Individual brakes on each castor, as well as hand brakes
- Modern streamlined design gives an up-to-date image to your rehab and physio areas
- Limitless adjustment in the arm and hand pads allows the grip to suit any patient.
- 4 pivoting castors, yet still has hand brakes. Safety and control at its best! (And..... all 4 castors still have Individual brakes as well!)
- Gas assist height adjustment – and two quick release levers. No more tri-bolts to frantically adjust and Lose. (This also safe guards against unexpected failure of which could result in a patient fall.)
- Slimline base to access nearly every doorway and bathroom.
- Highest quality castor for smoothest travel.
- Flat urethane arm pads – no more shaped/curved metal plates with vinyl pads that wear out.
- Modern design – looks smart, not mundane.



## Specifications:

<b>Product Code</b>	6755
<b>Overall Width</b>	660 - 790 mm
<b>Overall Height</b>	900 - 1300 mm
<b>Length</b>	810 mm
<b>Distance Between Arm Supports</b>	310 - 470 mm
<b>SWL</b>	200 kg