

Petite Seat Walker



The Petite seat walker is manufactured from durable lightweight aluminium. 6" castors make this an ideal outdoor walker. But is also compact enough for use indoors.

Product Benefits:

- Ideal for smaller people requiring mobility assistance
- Durable lightweight aluminium frame
- Folding for convenient transport and storage
- Height adjustable handles
- Ergonomically designed rubber hand grips
- Fitted with rear braking and front swivel castors
- Locking brake mechanism with reflectors
- Padded seat
- Supplied with reinforced vinyl bag
- Removable backrest

Specifications:

| | |
|-----------------------------|--------------------|
| Product Code | 6771 |
| Seat Width | 300 mm |
| Seat Depth | 300 mm |
| Overall Width | 500 mm |
| Overall Depth | 595 mm |
| Overall Folded WxDxH | 490 x 275 x 740 mm |
| Handle Height | 740 - 845 mm |
| Castor Size | 150 x 25 mm |
| Walker Weight | 6 kg |
| SWL | 175 kg |