

Navigator Forearm Walker



The navigator is foldable, staying firmly closed making it is easy to transport. The frame and seat are tested and approved on users with a maximum weight of 150 kg. Light and reliable brakes act on the TPE covered unbreakable wheels. The optimal and rigid design of the frame and the fore arm support ensures the best user comfort (lightness, turning circle and steadiness).

Product Benefits:

- The Navigator is available in size large.
- The Navigator is easily height adjustable in the range from 900 – 1150 mm

Specifications:

| | |
|----------------------|---------------|
| Product Code | 6782 |
| Overall Width | 550 mm |
| Overall Depth | 610 mm |
| Height | 900 – 1150 mm |
| Seat Height | 620 mm |
| Handle Height | 900 – 1150 mm |
| Walker Weight | 10.8 kg |
| SWL | 150 kg |