

Server HD Rollator



Server HD folds easily and stays that way with a simple locking system that doubles up as an ergonomically shaped carrying handle. Frame and seat are tested and approved for maximum user weight of 200 kg. The brake mechanism is light yet positive. TPE covered, unbreakable wheels.

Product Benefits:

- Frame and seat are tested and approved for maximum user weight of 200 kg.
- Server HD push handles are height adjustable with 10 positions, from 74 cm to 102 cm.
- User comfort is ensured with an optimal frame design and ergonomic push handles.

Specifications:

| | |
|-----------------------------|---------------|
| Product Code | 6784 |
| Width | 690 mm |
| Depth | 680 mm |
| Height | 740 - 1020 mm |
| Seat Width | 550 mm |
| Depth (When Folded) | 680 mm |
| Height (When Folded) | 860 mm |
| Width (When Folded) | 225 mm |
| SWL | 200 kg |