

Grande Seat Walker



Product Benefits:

- Ideal for those requiring mobility assistance
- Durable lightweight aluminium frame
- Folding for convenient transport and storage
- Height adjustable handles
- Ergonomically designed rubber hand grips
- Fitted with rear braking wheels and front swivel castors
- Locking brake mechanism with reflectors
- Padded seat
- Supplied with reinforced vinyl bag
- Removable backrest

Specifications:

Product Code	6861
Seat Width	440 mm
Seat Depth	300 mm
Seat Height	550 mm
Overall Width	570 mm
Overall Depth	670 mm
Overall Folded WxDxH	570 x 350 x 815 mm
Handle Height	815 - 925 mm
Castor Size	194 x 40 mm
Walker Weight	7.4 kg
SWL	200 kg