

# Flexi-Link Arm Reacher



Efficient, safe and simple to operate. Our patented Flexi-Link connection system is today's most advanced device for connecting the portable lift to the track trolley. Its removable wand allows carers to move around freely as there are no handles in place and travel of the hoist is not impeded by lights, fans or doors. Available in either 60cm or 90cm lengths.

## Specifications:

<b>Product Code</b>	7422FL
<b>Length</b>	600 mm

## Also Available:

<b>Product Code</b>	7423FL
<b>Length</b>	900 mm