

Ave2 Birthing Bed



Brings expectant mothers maximum safety and comfort.



Ave2 Birthing Bed







The AVE 2 birthing bed brings expectant mothers maximum safety and comfort. Thanks to the AVE 2, the birth takes place naturally according to their expectations. For healthcare staff the bed is a reliable and practical tool which improves the quality and facilitates the care of the woman in labour at all stages of birth.

Product Benefits:

- · Trendelenburg and CPR
- Extra battery backup
- ABS easily disinfect
- Lockout Control
- Generous 150mm Castor size
- No scissor-lift column system
- Central brake easy to access
- Easy transport with directional castor and pushing handles
- Removable headboard for anesthesiologist (optional)
- Nurse Call (optional)
- Tilting holders for IV pole

Specifications:	
Product Code	77120
Total length of the lying area	2132 mm
Minimum length in the default setting	1570 mm
Full Width	980 mm
Minimum Height	600 mm
Maximum Height	1050 mm
Back Section	-10°/+70°
Seat Section	0°/+18.5°
Trendelenburg Position	0°/-10°
Leg Rests Vertical	0°/+135°
Leg Rests Horizontal	0°/+60°
Foot Section Adjustment Angle	0°/-22°
Maximum Patient Load	240 kg
Overall Product Weight	236 kg
Lifting Pole Bearing Capacity	75 kg







SWL

