

# Standard Beside Commode



The standard bedside commode offers padded armrests for comfort and features adjustable height, a removable padded seat cover and a plastic toilet seat with bowl and lid.

## Product Benefits:

- Value for money
- Stackable design
- 125kg Safe Working Load

## Also Available:

### Standard Beside Commode with Padded Arms

**Product Code** 7952PA

## Specifications:

<b>Product Code</b>	7952
<b>Overall Width</b>	575 mm
<b>Overall Height</b>	940 - 1090 mm
<b>Overall Depth</b>	600 mm
<b>Seat Width</b>	440 mm
<b>Seat Height</b>	500 - 650 mm
<b>Seat Depth</b>	440 mm
<b>Back Height</b>	475 mm
<b>Toilet Seat</b>	475 - 625 mm
<b>SWL</b>	125 kg