

Horizon 4000 Patient Transport Trolley



Horizon 4000 Patient Transport Trolley is designed for use as both a day surgery bed and a recovery trolley. It is very easy to manoeuvre because it is equipped with locking anti static castors. The Trendelenburg permits the carer to regulate the blood circulation of the patient while being transferred from operating room to ward.

Product Features:

- Electric Hi Low
- Electric Backrest
- Electric Forward & Reverse Trendelenburg
- 150mm Twin wheel Anti-static Central Locking Castors
- Push Handles
- Battery Backup
- Corner Buffers
- APO Grey top frame powdercoat
- Ironstone base frame powdercoat

Optional Features (added at time of Manufacturer):

- Kneebreak
- Attendant Control
- CPR Emergency Release
- Pull out Extension

Accessories Available:

- Horizontal Rails
- Self Help Pole & Triangle
- Vertical Rails (as shown)
- I.V.Pole
- Horizontal Rails
- TBR2000 T-Bar Wall Protection Buffer
- Clamp on Support Rail



Specifications:

Product Code	7720
Overall Dimensions	840 x 2230mm
Sleeping Surface Dimensions	720 x 2025 mm
Height Range	390 - 830 mm
Underbed Clearance	160 mm
Back Angle	75°
TR / ATR	16° / 16°
Bed Weight	120 kg
SWL	250 kg