

Staebel EN7 Floorline Bed



Staebel EN7 Floorline Bed



The Staebel EN7 Floorline Bed is one of the lowest beds available on the market and designed for intensive use. The bed lowers to floor level to prevent injuries and accidents caused by falls from beds. All actuators and mechanisms are concealed beneath the platform, allowing the bed to be constructed with no end pedestals! This allows visibility for users, and nursing staff access for procedures.




Chair Position
push of a button



Bluetooth
handset must be ordered separately



75mm Dual Castors



APO Grey Powdercoat

Product Features:

- Extremely Low Minimum Position
- Auto-Regression Backrest
- 75mm Castors
- APO Grey Powdercoated Frame
- Chair Position at the Push of a Button
- Removeable Slender Arch Head & Footboard in Natural Oak
- Vertical Lift

Optional Features- *Added at time of manufacturer*

- Underbed Light (As shown)
- Bluetooth Handset
- 300kg Handset
- Battery Back Up
- Dual Breaking
- Attendant Control
- CPR Backrest



Specifications:

Product Code	7753	7755	7758
Size	Single	Wide	King Single
Overall Dimensions	920 x 2115 mm	1020 x 2115 mm	1090 x 2115 mm
Sleeping Surface Dimensions	900 x 2000 mm	1000 x 2000 mm	1070 x 2000 mm
Hi - Lo Position	100 - 725 mm	100 - 725 mm	100 - 725 mm
TR/RTR Angle	15° - 15°	15° - 15°	15° - 15°
Back / Knee Angle	62° / 35°	62° / 35°	62° / 35°
Bed Weight	110 kg	110 kg	110 kg
SWL	200 kg	300 kg	300 kg