

Multicare ICU Bed



**A combination of
advanced innovations
& functional features.**

Multicare Bed

The Multicare is an accomplished member of the critical care team. It simplifies interventions in critical care situations improving patient care through treatment and rehabilitation whilst providing comfortable working conditions for nursing staff. A combination of advanced innovations with functional features ensures that the Multicare exceeds your expectations.

Product Benefits:

- Automatic Lateral tilting (both sides to 30 degrees)
- Multi Protect 4-part siderails
- Multiboard features an LCD touchscreen (in both siderails)
- Mattress platform tilt and extension adjustment
- Tilt indicators (backrest, lateral tilt, TR/ATR tilt)
- Programmed positions: CPR, orthopaedic chair, 30° HOB, mobilisation, Trendelenburg
- Activation GO button, safety STOP button
- Service data interface, user communication setup, quick reference guide on LCD touchscreen
- Battery back-up indicator (2 batteries)
- The Multicare also includes patient control panel in the side guards
- The Multicare can be customised to suit your specific needs: there are four different kits
- Mobi Lift handles – assists the patient to exit the bed safely.
- I-brake – an automatic brake that is activated after 10 seconds if the bed remains unbraked and is plugged in. It protects against uncontrolled movement, impacts and falls caused by leaning on an unbraked bed.
- I-Drive – fifth castor for safe transporting and manoeuvring of the bed, with automatic retraction of the castor to the chassis. If the bed is plugged in, it is automatically retracted. The fifth castor runs from it's own battery.



Specifications:

Product Code	74200
External Dimensions (Shortest Mattress Platform Position)	1050 x 2150 mm
External Dimensions (Longest Mattress Platform Position)	1050 x 2370 mm
Mattress Dimensions	2080 x 860 mm
Maximum Mattress Height	230 mm
Height Adjustment	440 - 820 mm
Lateral Tilt	±30°
TR / ATR Angle	13° / 16°
Bed Weight	210 kg
SWL	250 kg

