

Image 3 Bed



**Designed to focus on
patient's safety and
mobilisation.**

Image 3 Bed



Product Benefits:

- Image 3, electrically adjustable bed respects the growing demands of the patient and nursing staff. Its attractive design will find universal use in a whole spectrum of hospital care. Image 3 deals with the patient's safety and mobilisation in a unique way: combining the SafeFree® siderails and reduced low height.
- Image 3 was designed with focusing on 100% safe recovery of patient in the bed. Safety is the key characteristics of Image 3. Reduced height combined with telescopic SafeFree® siderails create absolutely safe environment for the patient, that can be adjusted according to their specific needs.
- Also included is a complex system of supportive automatic features called 3-STOP STRATEGY. It contains functions for safe sleeping, independent mobilisation for quicker recovery, convenient care and cleaning procedures for hospital staff. This bed contains extended features which prevent pressure injuries and contribute to their easier treatment. Mattress platform profiling called Ergoframe® provides needed comfort for the patient.
- Furthermore, this bed is designed with emphasis to maximally ease everyday activities such as nursing procedures, transfers and other routine tasks. Ergonomics, intuitive controls, smooth manipulation with bed and easy cleaning bring relief to hospital staff.
- Image 3 is elegant and its simple design with clear lines impresses with pleasant, home looking atmosphere. Though it is very functional, it brings ergonomics and user-friendly handling. This way, the Image 3 can be easily used in any hospital area of care.



Specifications:

Product Code	74801
External Dimensions	1040 × 2200 mm
Clearance Height	130 / 155 mm
Mattress Dimensions	2000 x 860 x 900 mm
Mattress Platform Height Adjustment	280 - 800 mm / 300 - 820 mm
Height of Siderails (Split)	400 mm
Height of Siderails (Single Collapsible)	390 mm
Maximum Mattress Height (Split Siderails)	180 mm
Maximum Mattress Height (Single Collapsible Siderails)	170 mm
Bed Weight	145 kg
Backrest Angle	0 - 70°
Thighrest Angle	0 - 34°
TR / ATR Angle	0 - 14°
SWL	265 kg