

Eleganza 2 Bed with Junior Kit



**A system designed for
maximum safety for
smaller patients.**


Eleganza 2 Bed with Junior Kit



Product Benefits:

Maximum safety for small patients. Functionality and simple controls. This unique hospital bed, thanks to the side rail safety concept, is suitable for adult and paediatric patients.

- Patients from 90 to 220cm in height and can safely use the bed.
- Integrated Mobilift innovation with support handles
- Highest possible level safety level to eliminate falls or injuries and to indicate safe bed positioning
- Attractive design with innovative column construction
- 'Smart' patient monitoring features: safety monitor, bed monitor and vital monitor
- Ease of transport with fifth castor
- Increased patient comfort and ease of cleaning
- Patient control of some functions
- Practical and ergonomic accessories available



Specifications:

Product Code	76618
External Dimensions	995 x 2230 mm
Bed Extension	310 mm
Bed Height (Hi-Lo Adjustment)	395 – 775 mm
Internal Dimensions / (Compatible with Mattress Width)	2000 x 900 mm / (900 and 860 mm)
Max. Mattress Height	180 mm
Height of Split Plastic Side Rails	405 mm
Weight of Bed	140 kg

Clearance of Undercarriage	150 mm with 150 mm Tente castors
Ergoframe®	110 mm
Maximum Backrest angle	70°
Maximum Thigh Rest angle	35°
Calf Rest Angle	21°
TR/RTR Angle	16° / 16°
Lifting Pole Bearing Capacity	75 kg
SWL	250 kg