

Tom 2 Cot



**A safe & unique cot
designed for easy access.**

Tom 2 Cot



A state-of-the-art children's cot must be safe and be able to provide excellent access for nursing staff. The electric positioning of TOM 2 children's cot responds to these and further demands in its original concept of telescopic siderails, column construction mechanism and above standard safety parameters.

The Tom 2 bed fully complies with the new standard for hospital children's beds. DIN 32623:2009-11 and its construction guarantees perfect protection against children fall. Also complies with standards: EN 60601-1, EMC 60601-1-2, EMC 60601-1-6

Product Benefits:

- 3 height positions for side rails, 0cm, 40cm, 80cm
- Backrest retracts when it raises, meaning zero gap
- Integrated collapsible lockable side rails
- 4 x Tente Motion 125mm castors
- Central operated castors (one castor with steer system)
- 4 x holders for IV in corners of bed
- Undercarriage cover
- Removable head/foot ends
- Accessory holders: Pair of accessory rails
- Battery Backup
- Horizontal Wall/corner bumpers
- Quality European product

Specifications:

Product Code	76758
External Dimensions	1554 x 83 mm
Recommended Mattress Size	1420 x 700 / 100 mm
Mattress Platform Height (Electrically Adjusted)	580 - 830 mm
TA / ATR Tilt	+15° / -15°
Siderail Height Above Mattress Platform	50 mm / 400 mm / 800 mm
Autoregression of Backrest	120 mm
Maximum Gap Between Bckrest and Head End	25 mm
Clearance	130 mm
Maximum Backrest Angle	40°
Bed Weight	80 kg
SWL	75 kg

Also Available:

Tom 2 Cot Mattress

Product Code	76751
---------------------	-------