

Staebel Floorline Bed



**Introducing the first ever
universal bed that suits low
care, high care & bariatric care.**

Staebel Floorline Bed



From care of frail elderly or dementia patients, through to acute and bariatric settings, the patent pending Australian made Staebel Floorline Bed has you covered.

Product Benefits:

- Unique individual safe working load of 200kg or 300kg (simply swap the handsets)
- Low height of 100mm and high height of 725mm
- Electric kneebend
- Electric forward & reverse Trendelenburg
- Auto regression backrest
- 75mm twin wheel castors
- Customisable head and footboards (fini blond laminate)
- Universal solutions for low, high and bariatric care
- Chair position at the push of a button



Specifications:

Product Code	7753
Size	Single
Overall Width	920 mm
Overall Length	2080 mm
Hi - Lo Position	100 - 720 mm
TR/RTR Angle	15° - 15°
Back / Knee Angle	62° / 35°
Bed Weight	110 kg
SWL	200 kg

Size Available:

Product Code	7758
Size	King Single
Overall Width	1090 mm
Overall Length	2080 mm
Hi - Lo Position	120 - 720 mm
TR/RTR Angle	15° - 15°
Back / Knee Angle	62° / 35°
Bed Weight	110 kg
SWL	300 kg