

Alrick 2400 Acute Bed



Alrick 2400 Acute Bed

The Alrick 2400 Acute Bed is engineered with a robust design to ensure safety and stability. It boasts an impressive safe working load capacity of 400kg. This weight capacity ensures that the bed can safely accommodate individuals of varying weights and provides peace of mind to both users and caregivers.

The bed's sturdy construction, reinforced frame, and high-quality materials contribute to its exceptional weight capacity. Whether you require additional support due to medical conditions or simply prefer the reassurance of a higher weight limit, the Alrick 2400 Acute Bed delivers strength and durability from King Single width up to King Double width.

Product Features:


- APO Grey & Ironstone Powder Coated Frame
- Steel Frame in Ironstone with Classic White Laminate Infills
- Electric Hi-Low
- Electric Auto Regression Backrest
- Electric Kneebreak
- 150mm Twin Wheel Breaking Castors
- Trendelenburg and Anti-Trendelenburg
- Battery Back Up
- Manual Vascular Support
- Backlit Handset
- CPR Release
- Bluetooth connectivity suit the HB200 - handset to be purchased separately

Optional Features:

- Bluetooth handset
- Head & footboard options
- Pull out extension
- Attendant control
- USB Charging Ports
- Custom Powdercoat colours

Accessories Available:

- Horizontal Rails
- Clamp on Support Rail
- Self Help Pole & Triangle
- I.V. Pole
- T-Bar Wall Protection Buffer



Specifications:

Size	King Single	Extra Wide	Double	Queen	King Double
External Dimensions	1155 x 2195 mm	1285 x 2195 mm	1455 x 2195 mm	1585 x 2195 mm	1885 x 2195 mm
Sleeping Surface	1070 x 2000 mm	1200 x 2000 mm	1370 x 2000 mm	1500 x 2000 mm	1800 x 2000 mm
High Position	830 mm	830 mm	830 mm	830 mm	830 mm
Low Position	390 mm	390 mm	390 mm	390 mm	390 mm
Backrest / Knee Angle	77° / 40°	77° / 40°	77° / 40°	77° / 40°	77° / 40°
TR / ATR Angle	-16° / 16°	-16° / 16°	-16° / 16°	-16° / 16°	-16° / 16°
SWL	400 kg	400 kg	350 kg	350 kg	350 kg